

INSTALLATION QUALIFICATION

INSTRUMENT ISSUE MODEL NO	:	Turbochem 100
SERIAL NO	:	4611-1086
SOFTWARE VERSION	:	1.2.0.437
MANUFACTURER	:	AWARENESS TECHNOLOGIES INC, USA
DISTRIBUTOR	:	CPC DIAGNOSTICS PVT. LTD.,
		No. 324, KHADI GRAMMODYOG COMPLEX,
		AVVAI SHANMUGAM SALAI, GOPALAPURAM
		CHENNAI-86
CUSTOMER NAME	:	ALOKNANDA DIAGNOSTIC,NEW DELHI.
DATE OF INSTALLATION	:	21/01/2013 NOSTIC
ENGINEER (FOR CPC DIAGNOSTIC	CS PVT.	LTD.)
		2 Simple
SIGNATURE	:	
NAME	:	Mr. R.SURAJ KUMAR,
DESIGNATION	:	HEAD-CUSTOMER SUPPORT
COMPANY		CPC DIAGNOSTICS PVT. LTD.
CUSTOMER (ALOKNANDA DIAGNO	OSTIC,N	EW DELHI.)
SIGNATURE	:	

NAME

DESIGNATION

AM,



INSTALLATION REQUIREMENTS

Installation Environment Conditions

Place should not be subjected to direct sunlight Place should be flat and minimum of dust.

Power Supply

RAM

Voltage 220-240V AC 1.0 KVA Online UPS.

Temperature And Humidity Conditions

The temperature of the analyzer room should be between $18 - 30^{\circ}C$ The temperature fluctuation during the analysis should be $\pm 2^{\circ}C$

Dimensions Of System Main Units (Unit: cm)

TURBOCHEM100 Width 53cm x Length 50cm x Height 40cm Weight 15Kg DPR Width700 x Depth 700 x Height 1430

ANALYSER PC QUALIFICATION

POWER SUPPLY : 230/50 Hz AS PER REQUIREMENTS

2 GB

UPS : AVAILABLE

CURRENT : AS PER REQUIREMENTS

COM PORTS AVAILABLE : YES
ONLINE POSSIBLE : YES

ONLINE POSSIBLE : YES

PROCESSOR SPEED : 1.46GHz

STORAGE CAPACITY : 500 Gb

MONITOR : 15"

PROCESSOR : INTEL

OPERATING SYSTEM : WINDOWS 8

SOFTWARE CD : AVAILABLE

PRINTER PORT : AVAILABLE

VERIFIED BY: Mr. R.SURAJ KUMAR

REVIEWED BY:

(FOR ALOKNANDA DIAGNOSTIC, NEW DELHI.)

SIGNATURE:

New No.324, Old No.237/1, Khadi Gramodyog Complex, Avvai Shanmugam Salai, Gopalapuram, Chennai - 600 086. Phone: +91 44 2346 0161(4 lines), Fax: +91 44 2346 0165

CIN No.: U52599TN1991PTC021095

SIGNATURE:

Corporate Office: #9, V Floor, Gokul Tower,

No.9&10 C.P. Ramaswamy Road, Alwarpet, Chennai - 600 018.

Phone: +91 44 24993989, Fax: +91 44 24993357

E-Mail: sales@cpcdiagnostics.in Website: www.cpcdiagnostics.in



PERFORMANCE QUALIFICATION

FILTER VOLTAGES – All wave length voltage should be within 2 to 10 volts

Filters(nm)	Volts(v)
340	5.97
405	6.41
505	6.89
545	6.95
580	7.08
630	7.19

PRECISION CHECK

This module is highlighting the reproducibility and pip petting accuracy of the system. Measuring criteria

- 1. One end point chemistry
- 2.One kinetic chemistry
- 3.Large sample chemistry
- 4. Small sample chemistry
- 5.R1 only
- 6.R1 + R2

BY DOING PRECISION CHECK WITH GLUCOSE, ALT, CREATININE ABOVE REQUIREMENT WILL BE COMPLETED.

DETAILS OF REAGENTS, CALIBRATORS AND CONTROLS

REAGENTS

: Jeev reagents

CALIBRATOR

: Identi Chemistry Callibrator (LOT NO 27189)

CONTROL

: Randox Hum LEVEL 2-844un

HERE BY I'm VERIFYING THAT THE CONTROL RESULT OF THE EXECUTED PARAMETERS ARE WITHIN THE SD RANGE AND CV RANGE.

VERIFIED BY: Mr. R.SURAJ KUMAR

REVIEWED BY:

(FOR ALOKNANDA DIAGNOSTIC, NEW DELHI.)

SIGNATURE:

SIGNATURE:

New No.324, Old No.237/1, Khadi Gramodyog Complex, Awai Shanmugam Salai, Gopalapuram, Chennai - 600 086. Phone: +91 44 2346 0161(4 lines), Fax: +91 44 2346 0165

CIN No.: U52599TN1991PTC021095

Corporate Office: #9, V Floor, Gokul Tower,

No.9&10 C.P. Ramaswamy Road, Alwarpet, Chennai - 600 018.

Phone: +91 44 24993989, Fax: +91 44 24993357

E-Mail: sales@cpcdiagnostics.in Website: www.cpcdiagnostics.in



OPERATIONAL QUALIFICATION

ANALYSER ALIGNMENT CHECK	REMARK		
 □ PROBE Vs RACK MOVER □ PROBE Vs CUVETTES □ PROBE Vs WASHCUP □ PHOTOMETER VS CUVETTE CARRIER 	OK OK OK OK		
RACK MOVEMENT CHECK			
□ RACK MOVEMENT TIME	OK		
CUVETTE CARRIER / READ SHUTTLE MOVEMENT CHECK			
□ CUVETTE CARRIER Y MOVEMENT TIME □ READ SHUTTLE X MOVEMENT TIME	OK OK		
CUVETTE CARRIER / LAMP			
□ LAMP CONDITION □ CUVETTE TEMPERATURE	GOOD OK		
PROBE MECHANISM			
 □ PROBE X MOVEMENT TIME □ PROBE Z MOVEMENT TIME □ PROBE WASHING □ PROBE LIQUID LEVEL DETECTION 	OK OK OK OK		
SYRINGE CONDITION			
□ SYRINGE MOVEMENT □ SYRINGE LEAK	OK NO		

VERIFIED BY: Mr. R.SURAJ KUMAR

SIGNATURE:

REVIEWED BY:

(FOR ALOKNANDA DIAGNOSTIC, NEWDELHI.)

SIGNATURE:

New No.324, Old No.237/1, Khadi Gramodyog Complex, Avvai Shanmugam Salai, Gopalapuram, Chennai - 600 086. Phone: +91 44 2346 0161(4 lines), Fax: +91 44 2346 0165 CIN No.: U52599TN1991PTC021095

Corporate Office: #9, V Floor, Gokul Tower,

No.9&10 C.P. Ramaswamy Road, Alwarpet, Chennai - 600 018.

Phone: +91 44 24993989, Fax: +91 44 24993357

E-Mail: sales@cpcdiagnostics.in Website: www.cpcdiagnostics.in