


cobas[®]
Life needs answers

cobas Integra 400 plus Analyzer

Qualification Service
Installation Qualification / Operation Qualification (v.1.0)





cobas Integra 400 *plus*



General Information

Country: **India**

Customer Name: **Care Diagnostics**

Customer Address: **3rd Cross, oppo. ICICI Bank, Neeladri Nagar, Electronic city, Phase-1 Bengaluru-100**

Person Responsible for Quality Assurance: **Dr. Meenu.EV**

System Information

cobas Integra 400+: **S/N 420027**

Host provider: **NA**

Software Version: **3.6**

Installation Information

Installation Start Date: **18.02.2020**

First Installation: **yes**

Reconfiguration: From: **specify** To: **specify**

Relocation: From: **specify** To: **specify**

Roche Responsible Representative : **Mr. Sumith**



Installation Qualification:

This document forms the basis of the Qualification Services Certificate. It certifies that the instrument is installed according to the manufacturer's specifications. The report presents and documents the test procedures, the documentation, reference and acceptance criteria used to verify that the system is installed according specifications. The report demonstrated that all installation qualification criteria have been met satisfactorily.

Notice: The following tests are to be carried out by trained Roche personnel only.

Purpose: The purpose of this test is to confirm that the instrument was delivered undamaged and installed correctly.

Test #	Test	Pass Fail	Signature Date
IQ.1.1	Operator's Manual available	Pass	20.02.2020
IQ 1.2	Environmental parameters met	Pass	20.02.2020
IQ 1.3	Instrument delivered undamaged and complete	Pass	20.02.2020 20.02.2020
IQ 1.4	Transport locking successfully removed	Pass	20.02.2020
IQ 1.5	All connections correctly installed	Pass	20.02.2020
IQ 1.6	Instrument positioned according to Installation Manual	Pass	20.02.2020 20.02.2020
IQ 1.7	Instrument boot process successfully	Pass	20.02.2020
IQ 1.8	Checksum according to specification	Pass	20.02.2020
IQ 1.9	Mechanical adjustments complete	Pass	20.02.2020
IQ 1.10	Auxiliary components positioned	Pass	20.02.2020
IQ 1.11	Instrument installation check	Pass	20.02.2020
IQ 1.12	Host communication settings checked	yes	20.02.2020

Test #	Test	Pass Fail	Signature Date
IQ.2	Installation Qualification for cobas Integra 400+	yes	20.02.2020



Deviation Report: Any discrepancies found during the installation must be documented in the space below. Roche personnel will then investigate the deviation and decide upon the most appropriate action to be taken.

Deviation #1

Investigation

Action taken

Deviation resolved satisfactorily? specify

Deviation #2

Investigation

Action taken

Deviation resolved satisfactorily? specify

Deviation #3

Investigation

Action taken

Deviation resolved satisfactorily? specify



Operational Qualification:

This document is the basis of the Qualification Service Certificate. It certifies that the instrument is operating according to the manufacture's specifications. This report presents and documents the test procedures, documentation, references and acceptance criteria used to verify that the specified system is operating according to the specifications. The report demonstrates that all operational qualification criteria have been met satisfactorily.

Notice: The following tests are to be carried out by trained Roche personnel only.

Purpose: The purpose of this test is to check that the modules are operating in accordance with the

Test #	Test	Pass Fail	Signature Date
OQ.1	Calibration successfully	Pass	Signature 21/02/2020
OQ.2	Quality Control successfully	Pass	
OQ.3	Accuracy check successfully	Pass	

Deviation Report: Any discrepancies found during the installation must be documented in the space below. Roche personnel will then investigate the deviation and decide upon the most appropriate action to be taken.

Deviation #1

Investigation

Action taken

Deviation resolved satisfactorily? specify

Deviation #2

Investigation

Action taken

Deviation resolved satisfactorily? specify



Conclusion

All test results are acceptable.

yes

Any deviation or non-conformances observed have been recorded as a deviation and the relevant forms completed.

yes

All acceptance criteria have been met. This equipment is deemed acceptable and the unit is approved for its intended use.

yes

Comments

Instrument working satisfactorily

Completed by Roche Representative

Date

21/02/2020

Print Name Loy Mascarenhas
Application Specialist

Signature

Reviewed by Customer Contact

Date

22/02/2020

Print Name Dh. Meenu E. V

Signature

Reviewed by Customer Quality Assurance

Date

Print Name _____

Signature _____



Installation Qualification for cobas[®] Integra 400 plus

Description

IQ.1.1	Operator's Manual available	
	Check that a copy of the latest version of the Operator's Manual is available.	Pass
IQ 1.2	Environmental parameters	
	Ambient temperature in the lab is between 15° and 32 °C	Pass
	Relative Humidity maximum of 50% at 32 °C and non-condensing	Pass
	Bacteria free, deionized water < 10 cfu/ml	Pass
	Water conductivity 1.0 µS/cm or less	Pass
	Dust and Vibration free	Pass
	Instrument is not exposed to direct sunlight	Pass
IQ 1.3	Instrument delivered undamaged and	
	All covers are undamaged	Pass
	All accessory boxes are delivered	Pass
	Instrument does not show any external damage	Pass
IQ 1.4	Transport locking successfully removed	
	Unpacking of the different modules and accessories without damage to units	Pass
IQ 1.5	All connections correctly installed	
	Power supply voltage at the customer facility:	230V ±2V
	UPS system available:	yes
	Voltage fluctuation less than 230 ±5V	Pass
	Grounding less than 1.0 V	Pass



IQ 1.6	Instrument positioned according to Installation Manual	
	System layout is according to the description in the	Pass
IQ 1.7	Instrument boot process successful	
	IP address configuration successful	Pass
	System Configuration successful	Pass
	First system boot-up	Pass
IQ 1.8	Checksum according to specification	
	Version of installed cobas Integra 400+ software	3.5.2
	Installation of country language successful	yes
IQ 1.9	Mechanical adjustments complete	
	All mechanical adjustments are carried out	Pass
IQ 1.10	Auxiliary components positioned	
		Pass
IQ 1.11	Instrument installation check	
	Print function	yes
	Rack/Sample barcode read check	Pass



IQ 1.12	Host communication settings checked	
	Check Host settings according to Host manual	yes
	Check Host communication	yes

Deviation Report: Any discrepancies found during the installation must be documented in the space below. Roche personnel will then investigate the deviation and decide upon the most appropriate action to be taken.

Deviation #1

Investigation

Action taken

Deviation resolved satisfactorily? specify

Deviation #2

Investigation

Action taken

Deviation resolved satisfactorily? specify

Deviation #3

Investigation

Action taken

Deviation resolved satisfactorily? specify



Installation Qualification for cobas® Integra 400 plus:

Description

IQ.2.1	Function check of Integra400+ module according to specifications	
	System layout is according to the description in the manual	Pass
	Integra 400+ is installed according to the installation manual and using official tools	Pass
IQ.2.2	Mechanical adjustments complete	
	All mechanical adjustments for the different Integra 400+ mechanical parts are carried out	Pass
IQ.2.3	Auxiliary components positioned	
	Wash solutions are installed at the Integra 400+	Pass
	ISE electrodes are installed	not applicable
	ISE solutions are installed	not applicable
	Probe B & Probe C (Reagent & Sample) pipettors installed	Pass
IQ.2.4	Instrument installation check	
	Air water Calibration	Pass
	Prime Fluid System	Pass
	Analyzer Rotor (Reaction) temperature 37°C ± 0.5°C	Pass



Carry out Instrument Check according to Method Sheet of the cobas c pack INSTC (Art. No. 04851013 190) (attached printout)

Pass

ISE Check 20 times (attached printout)

not applicable

IQ 2.5 Assay installation

Download of applications from TAS (attached list of applications)

Pass

Load corresponding reagent c-packs

Pass

Deviation Report: Any discrepancies found during the installation must be documented in the space below. Roche personnel will then investigate the deviation and decide upon the most appropriate action to be taken.

Deviation #1

Investigation

Action taken

Deviation resolved satisfactorily? specify

Deviation #2

Investigation

Action taken

Deviation resolved satisfactorily? specify



Operational Qualification:

Notice: The steps described in OQ.1 have to be carried out after a new system installation and after any repair action which requires additional calibration.

If the service action does not affect the measurement performance, only apply steps OQ.2 and OQ.3 of the Operation Qualification.

Description

OQ 1 Calibration

Calibration of all photometric parameters successful (attached printout)

yes

Calibration of all ISE parameters successful (attached printout)

not applicable

OQ.2 Quality Control

Specify the type of control used:

Preci control (Roche)

QC of all photometric parameters within acceptable range (attached printout)

yes

QC of ISE parameters within acceptable range (attached printout)

not applicable

OQ.3.1 Accuracy check for ISE

Perform test with analytical reagents

	Number of det.
Na	21
K	21
Cl	21

Sample solution: PNU (code 300 / Cat. No. 10171735, 10171743, 10651257).
Fill 21 Hitachi cups with PNU (code: 300) and perform Na, K and Cl tests. Calculate the CV.

Accuracy check for ISE was within acceptable range not applicable



OQ.3.2 Accuracy check for Photometric Assays

Perform test with analytical reagents

	Number of det.
2-point/end-point Assay	21
Rate A Assay	21

Sample solution: Precicontrol (Roche Control)
Fill cobas cups with Precicontrol & perform 21 determinations of each parameter.

Accuracy check for Photometric Assays was within acceptable range **yes**

OQ.3.2 Precision check for Photometric Assays

Perform test with analytical reagents

	Number of det.
2-point/end-point Assay	21
Rate A Assay	21

Sample solution: Precicontrol (Roche Control)
Fill cobas cups with Precicontrol & perform 21 determinations of each parameter.

Precision check for Photometric Assays was within acceptable range **yes**

Deviation Report: Any discrepancies found during the installation must be documented in the space below. Roche personnel will then investigate the deviation and decide upon the most appropriate action to be taken.

Deviation #1

Investigation

Action taken

Deviation resolved satisfactorily? **specify**



Attachments

1. Precision - Intra & Inter assay
2. Accuracy Check



(ISO 9001:2015 Certified)

Instrument Commissioning Report

Name of Institution: Care Diagnostics
 Complete Address: 3rd Cross, oppo. ICICI Bank, Neeladri Nagar, Electronic city, Phase-1.
 City: Bengaluru Pin Code: 560100
 Person Incharge : Dr.Meenu EV REXIS ORD : ORD-0014381412
 Instrument Model: cobas Integra 400Plus Serial No.: 4220027
 SW Version : 3.6 Training Completed On : 21-02-2020
 Call Received Date: 19-02-2020 Call Attended Date: 20-02-2020
 Travel Hours : 9+2+2+9 Work Hours: 9+9


System Configuration & Preliminary application work:

(Please check in the Box)

- | | | | |
|---|-------------------------------------|---------------------------------|-------------------------------------|
| A. System configuration & programming | <input checked="" type="checkbox"/> | E. Assay Calibration | <input checked="" type="checkbox"/> |
| B. Hardware Overview | <input checked="" type="checkbox"/> | F. QC Run | <input checked="" type="checkbox"/> |
| C. Software Overview | <input checked="" type="checkbox"/> | G. User Maintenance | <input checked="" type="checkbox"/> |
| D. Sample & QC Processing | <input checked="" type="checkbox"/> | H. Precision Check & Evaluation | <input checked="" type="checkbox"/> |
| I. Basics of Telephone Trouble Shooting | <input checked="" type="checkbox"/> | | |

Training Feedback:

(Please rate as per the Scale)



- | | |
|--------------------------------------|---------------------------------|
| Structure | <input type="text" value="10"/> |
| Objective Reached | <input type="text" value="10"/> |
| Ratio of Lecture to Practical time | <input type="text" value="10"/> |
| Quality of Lecture | <input type="text" value="10"/> |
| Completeness of Information | <input type="text" value="10"/> |
| Knowledge | <input type="text" value="10"/> |
| Openness to Questions | <input type="text" value="10"/> |
| Overall Training Satisfaction | <input type="text" value="10"/> |

Comments

Signature of Trainer [Signature] Signature & Seal of Customer [Signature] Date 22/2/2020



Name of Trainer Loy Mascarenhas Trainer's Comment: Instrument Working Satisfactory

"In case of Complaints/Inquiries please reach our Customer Support Centre at 1800-123-7599/044-30413900"



Roche Professional Services
Engineering Installation Report

Name of Institution: CARE DIAGNOSTICS
 Complete Address: 3RD CROSS, NEELADRI NAGAR, ELECTRONIC CITY PHASE 1
 City: BENGALURU Pin Code: 560100 Telephone No: 5610702227
 Person In-Charge: SUDHARSHAN.M Email ID: Care2lab@gmail.com
 Instrument Model: COBAS INTEGRAL 400 PLUS Serial No: 420027
 Software Ver: VGR: 3.6.2 Cobas Link Serial No: NA
 Department (kindly \checkmark): Pathology Blood Bank Other (mention) _____
 Installation Commencement Date: _____ Date of Installation: _____
 Travel Hours _____ Work Hours _____ REXIS ORD No ORD-0014381410

◆ Transportation Damage & Discrepancy Report: (please update after inspection)


A. List of Missing items: - NIL -
 B. Notes on any damages: - NIL -
 C. General Remark: - NIL -

◆ Installation checks & Summary: (please indicate status)

A. Power Supply:	OK <u>P-N: 230.2V, P-E: 230.6V, N-E: 0.4V</u>
B. Software Installation & Boot-up:	OK
C. Mechanical Adjustments & Setup:	OK
D. Fluidic Adjustment & Setup:	OK
E. System Checks & Status Print Outs:	OK

Installation Feedback:

(Please rate your Satisfaction Level as per the Scale)



- Pre-Installation Information
- Service Personnel's response in general
- Product installation & its timely Completion
- Overall Installation Satisfaction

Comments



Signature of Engineer [Signature] Signature & Seal of Customer [Signature] Date 20/2/2020

Name of Engineer SUMITH.P Engineer's Comment INSTALLATION COMPLETED SUCCESSFULLY

"In case of Complaints/Inquiries please reach our Customer Support Centre at 1800-123-7599/044-30413900"
 Version No.: 05 RPS/IN STL/SOP03/R/EIRBI/01 Page 1 of 1

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 11:06

MA2.1 Check Cuvette Transport Performance

Description

PURPOSE

Allows moving cuvettes from the Cuvette Buffer to the Waste Box during a selected cycle time and duration.
Allows checking system capacity in an ongoing test.
Allows manually checking the Cuvette Transport System LBs and sensor

Cycle Time [3..15 Sec] 7

Wait Cycles [0..65535 Cyc] 0

Duration [1..65535 Sec] 120

Status Sensors and LB's:

Cuvette Reservoir LB Level not okay

Cuvette Bowl Sensor no Cuvette detected

Cuvette Buffer max LB Level not okay

Cuvette Buffer min LB Level okay

Test Status okay

Time elapsed [s] 120

Theoretical Nr. of Cuvettes to be transported 17

Actual Nr. of Cuvettes transported 18

Start

Logfile... Print Cancel

OK

INTEGRA 400 plus

System Version: 3.6.2.1904
System Serial Number: 420027

Cui Version: 3.6.2

Cui Build: 1904

Wednesday, February 19, 2020 / 11:27

M24.1 Get Abs Air/Water Calibration Values

Description

PURPOSE

The Abs Photometer Calibration determines the Air/Water Correction values. These values are needed for result calculations. Displays the number of cuvettes outside the range for every measurement procedure.

Abs Air Water Calibration

Number of Cuvettes 50

Number of Measurements 20

Calibrate

Nr. of Cuvettes outside Range

Air Measurements	0
H2O Measurements	0
H2O - Air Measurements	0
Total	0

Accepted Nr. of Cuvettes outside Range 10 (not more then 20%)

Store Calibration Data

Abs Air/Water Calibration: okay

View Calibration Data...

Logfile... Print Cancel

ADK

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 11:28

M24.1 Get Abs Air/Water Calibration Values - Data

Description

PURPOSE

The Abs Photometer Calibration determines the Air/Water Correction values. These values are needed for result calculations. Displays the number of cuvettes outside the range for every measurement procedure.

Wavelength	A_CUV (H2O)	A (H2O)
1. 340 nm	0.0969	0.0137
2. 378 nm	0.0959	0.0154
3. 409 nm	0.0966	0.0174
4. 480 nm	0.0997	0.0214
5. 512 nm	0.1009	0.0228
6. 520 nm	0.1012	0.0231
7. 800 nm	0.0135	-0.0386
8. 552 nm	0.1022	0.0241
9. 583 nm	0.1033	0.0251
10. 629 nm	0.1046	0.0264
11. 652 nm	0.1046	0.0267
12. 659 nm	0.1042	0.0264

Data Source
New Values

Creation Date:
19-FEB-2020 11:28:17

Into Database
 Into a File

Store Values

Logfile... Print Cancel

OK

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 11:44

M13.1 Check Workstation Accuracy

Description

PURPOSE

Allows verifying the adjustment of all four workstations in an ongoing test. The check results are displayed as okay or failed.

Number of Cycles (min. 20) Start

Accuracy Check WS In/Out	okay
Accuracy Check Workstation B	okay
Accuracy Check Workstation C	okay
Accuracy Check Workstation FP	okay
Accuracy Check all Workstations	okay

View Data...

Logfile... Print Cancel

OK

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 11:44

M13.1 Check Workstation Accuracy - Data

Description

PURPOSE

Displays the measurements of all 12 wavelengths after checking workstation accuracy.
The measurements are in the <ABSWSAccuracy.log> file.

Wavelength	StdDev WS InOut <0.0003 Abs	StdDev WSB <0.0003 Abs	StdDev WSC <0.0003 Abs	StdDev WSFP <0.0003 Abs
1. 340 nm	0.0002	0.0002	0.0000	0.0003
2. 378 nm	0.0000	0.0000	0.0000	0.0002
3. 409 nm	0.0000	0.0000	0.0000	0.0002
4. 480 nm	0.0000	0.0000	0.0000	0.0002
5. 512 nm	0.0000	0.0000	0.0000	0.0002
6. 520 nm	0.0000	0.0000	0.0000	0.0002
7. 800 nm	0.0000	0.0000	0.0000	0.0000
8. 552 nm	0.0000	0.0000	0.0000	0.0002
9. 583 nm	0.0000	0.0000	0.0000	0.0002
10. 629 nm	0.0000	0.0000	0.0000	0.0002
11. 652 nm	0.0000	0.0000	0.0000	0.0002
12. 659 nm	0.0000	0.0000	0.0000	0.0002

Logfile... Print Cancel

AK

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 11:51

M13.2 Check Analyzer Cuvette Handling

Description

PURPOSE

Allows verifying Analyzer handling in an ongoing mechanical test. The WS In/Out inserts a new Cuvette into the rotor every machine cycle. Up to 65530 Cuvettes can be cycled.
This procedure runs without any Abs measurement. The results are

Enable WSB
 Enable WSC
 Enable WS FP

Number of Cuvettes used: 10

Status Cuvette Handling Check: okay

Start

Print Cancel

OK

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 11:53

M12.2 Check Analyzer Rotor Init Position

Description

PURPOSE

Allows verifying the Analyzer Rotor adjustment and checking the Analyzer Rotor Init Position.
Duration 1 min, 20 sec.

Check No.	No. of steps
1. Check	3
2. Check	3
3. Check	4
4. Check	4
5. Check	4
6. Check	4
7. Check	4
8. Check	4
9. Check	4
10. Check	4

Tolerance:
No. of steps
allowed range 3 to 5
optimum is 4

Rotor Init Position Check okay

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 12:01

M72.3 Check Init Position Accuracy

Description

PURPOSE

Allows carrying out an Init Position accuracy check for the Main Init and probes. Displays the difference between the first and subsequent checks. At least 10 sequences should be run. Duration: 5 min, 20 sec. (10 sequences)

Initialize Robotic Transfer

No. of Transfer sequences: 10 Start Transfer Check

Axis	Main Init difference [steps]	Probe difference [steps]
X	2	2
Y	1	2
Z	4	4

Check failed after sequence:

Status Transfer Accuracy Check: okay

Logfile... Print Cancel

Atk

INTEGRA 400 plus

System Version: 3.6.2.1904
System Serial Number: 420027

Cui Version: 3.6.2

Cui Build: 1904

Wednesday, February 19, 2020 / 13:33

Results Sample -Statistics

Statistics for all calculated Tests on Patient:

Patient ID CHECK B
Name CHECK B

# Tests	16		
CHKBR	Mean	SD	CV%
[Abs]	0.40605	0.001	0.33
Range	Low	High	
Extremes	0.40385	0.40847	
68% Confidence	0.40470	0.40741	
95% Confidence	0.40335	0.40876	

Test:

Handwritten signature

INTEGRA 400 plus

System Version: 3.6.2.1904
System Serial Number: 420027

Cui Version: 3.6.2

Cui Build: 1904

Wednesday, February 19, 2020 / 13:33

Results Sample - Statistics

Statistics for all calculated Tests on Patient:

Patient ID CHECK B
Name CHECK B

# Tests	16		
CHKBS	Mean	SD	CV%
[Abs]	1.22513	0.004	0.34
Range	Low	High	
Extremes	1.21853	1.23084	
68% Confidence	1.22093	1.22933	
95% Confidence	1.21672	1.23354	

Test:

CHK

INTEGRA 400 plus

System Version: 3.6.2.1904
System Serial Number: 420027

Cui Version: 3.6.2

Cui Build: 1904

Wednesday, February 19, 2020 / 13:33

Results Sample - Statistics [X]

Statistics for all calculated Tests on Patient:

Patient ID CHECK B
Name CHECK B

# Tests	16		
CHKBSR	Mean	SD	CV%
[Abs]	0.38728	0.001	0.30
Range	Low	High	
Extremes	0.38552	0.39021	
68% Confidence	0.38613	0.38843	
95% Confidence	0.38499	0.38957	

Test: [v]

OK

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 13:32

Results Sample - Statistics

Statistics for all calculated Tests on Patient:

Patient ID CHECK C
Name CHECK C

# Tests	16		
CHKCR	Mean	SD	CV%
[Abs]	0.40528	0.001	0.23
Range	Low	High	
Extremes	0.40371	0.40660	
68% Confidence	0.40435	0.40621	
95% Confidence	0.40342	0.40714	

Test:

Att

INTEGRA 400 plus

System Version: 3.6.2.1904

Cui Version: 3.6.2

Cui Build: 1904

System Serial Number: 420027

Wednesday, February 19, 2020 / 13:32

Results Sample - Statistics

Statistics for all calculated Tests on Patient:

Patient ID CHECK C
Name CHECK C

# Tests	16		
CHKCS	Mean	SD	CV%
[Abs]	1.24370	0.008	0.66
Range	Low	High	
Extremes	1.22921	1.25679	
68% Confidence	1.23548	1.25192	
95% Confidence	1.22726	1.26014	

Test:

AKH

INTEGRA 400 plus

System Version: 3.6.2.1904
System Serial Number: 420027

Cui Version: 3.6.2

Cui Build: 1904

Wednesday, February 19, 2020 / 13:32

Results Sample - Statistics

Statistics for all calculated Tests on Patient:

Patient ID CHECK C
Name CHECK C

# Tests	16		
CHKCSR	Mean	SD	CV%
[Abs]	0.38791	0.001	0.36
Range	Low	High	
Extremes	0.38483	0.38987	
68% Confidence	0.38652	0.38931	
95% Confidence	0.38512	0.39071	

Test:

OKL



CALIBRATION CERTIFICATE

HOSPITAL/LAB NAME CARE DIAGNOSTICS
ADDRESS, BANGALORE.
 KARNATAKA
ANALYSER COBAS INTEGRA 400 PLUS
SERIAL NO 420027
INSTALLATION PLACE LABORATORY.
DATE 20, FEBRUARY, 2020

This is to certify that the above mentioned analyzer was calibrated by the following procedure as per manufacturer's recommendations

1. Calibration of FP High Voltage and k/l factors
2. Workstation Accuracy check
3. Analyzer Rotor Init position check
4. Check Abs sensitivity 100% and 0%
5. Cuvette transport performance check
6. CHECK INIT POSITION ACCURACY

All the checks and calibration were done and system is being used under recommended conditions.
Next calibration is due on 19th February 2021.

For Roche Diagnostics India Pvt. Ltd.

Sumith P
Sr.Technical Service Specialist
Bangalore



CALIBRATION CERTIFICATE

Customer Name : Care Diagnostics
3rd cross, Neeladri Nagar,
Electronic City Phase I, opposite to ICICI Bank,
Bengaluru, Karnataka 560100

Instrument Model : DIESTRO ELECTROLYTE INSTRUMENT

Instrument Serial No. : 2157

Calibration Done on : 06/03/2020

Calibration Due on : 05/03/2021

THE FOLLOWING CHECKS AND SETUPS ARE CARRIED OUT FOR DIESTRO ELECTROLYTE INSTRUMENT FOR PERFORMANCE QUALIFICATION.

TRI LEVEL QC Exp Date 2021/10 Level -1

PARAMETER	MEASURED (MMOL/L)	RANGE (MMOL/L)
Na+	121.1	110 – 122 MEAN (116)
K+	2.01	1.85- 2.16 MEAN(2.00)
Cl-	83.0	71 – 85 MEAN (77)

TRI LEVEL QC Exp Date 2021/10 Level -2

PARAMETER	MEASURED (MMOL/L)	RANGE (MMOL/L)
Na+	140.6	132 – 147 MEAN (140)
K+	3.98	4.13 – 4.51 MEAN(4.47)
Cl-	100.2	89 – 104 MEAN(99)

TRILEVEL QC

Exp Date 2021/10

Level -3

PARAMETER	MEASURED (MMOL/L)	RANGE (MMOL/L)
Na+	157.3	153- 160 MEAN (161)
K+	8.03	6.61 - 7.70 MEAN (7.15)
Cl-	120.1	107 - 125 MEAN(116)

TWO POINT CALIBRATION FOR ELECTRODES

PARAMETER	CAL VALUES (MMOL/L)	SLOPE (MMOL/L)
Na+	57.61	50 - 64
	57.30	
K+	61.04	50 - 64
	60.63	
Cl-	50.86	40 - 64
	49.19	

REPEATABILITY TEST FOR SAMPLE


PARAMETER	CAL VALUES (MMOL/L)	RANGES (MMOL/L)
Na+	134.5	135 - 148
	134.4	
K+	4.87	3.50 - 5.30
	4.80	
Cl-	99.8	99 - 107
	98.8	

REMARKS: ALL MEASURED VALUE IS WITHIN THE TARGET RANGE FOR ALL THREE LEVELS OF CONTROLS AND CALIBRATION.

Protocol performed by: JK Biomed

Name : KALAIARASAN TR
Signature : 
Designation: SERVICE MANAGER
Date : 06/03/2020

Customer Authorization:

Name : SUDHARSAAN M
Designation: LAB MANAGER
Signature : 
Date : 09/03/2020