



UF Commissioning Report

Customer:	CHRISTIAN INSTITUTE OF HEALTH SCIENCES & RESEARCH	Date:	12 September, 2029
Model:	UF-4000	Engineer:	JAKIR HUSSAIN MONDAL
Serial No :	11929	A.Specialist:	Jayant Kumar Samal
Date Install:	12-Sep-20	BL Date:	

1 BACKGROUND VERIFICATION

1.1 Before performing adjustments to the Urinalysis analyzer, ensure there is no background error appeared to the analyzer.

Table 1 - result of background check

Parameters	Results	Acceptable Range	Status
RBC	0.0	$\leq 1.0 / \mu\text{L}$	Pass
NL RBC	0.0	$\leq 1.0 / \mu\text{L}$	Pass
WBC	0.0	$\leq 1.0 / \mu\text{L}$	Pass
WBC Clumps	0.00	$\leq 1.0 / \mu\text{L}$	Pass
EC	0.0	$\leq 1.0 / \mu\text{L}$	Pass
Squa.EC	0	$\leq 1.0 / \mu\text{L}$	Pass
Non SEC	0	$\leq 1.0 / \mu\text{L}$	Pass
Tran EC	0	$\leq 1.0 / \mu\text{L}$	Pass
RTEC	0	$\leq 1.0 / \mu\text{L}$	Pass
CAST	0	$\leq 1.0 / \mu\text{L}$	Pass
Hy. CAST	0	$\leq 1.0 / \mu\text{L}$	Pass
Path.CAST	0	$\leq 1.0 / \mu\text{L}$	Pass
BACT	0	$\leq 5.0 / \mu\text{L}$	Pass
X'TAL	0	$\leq 1.0 / \mu\text{L}$	Pass
YLC	0	$\leq 1.0 / \mu\text{L}$	Pass
SPERM	0	$\leq 1.0 / \mu\text{L}$	Pass
MUCUS	0	$\leq 1.0 / \mu\text{L}$	Pass
SF_TC	123	$\leq 3000 \text{ Count}$	Pass
CW_TC	43	$\leq 300 \text{ Count}$	Pass
CB_TC	103	$\leq 3000 \text{ Count}$	Pass

Paste the Result of Background result screen

3 OPTICAL AXIS VERIFICATION

Material Used : Latex DUKE 4207A (7um) (P/N: CX472678)
 Latex DUKE 4010A (1um) (P/N: 52123572)
 Latex Fluorescent C47410(Marketing name: A-16500) (P/N: BP783854)

3.1 Optical Axis adjustment result must be done by executing FINE [Optical Axis Adjustment]. The results are shown here:

Table 3 - Laser Latex run results

Parameter	Latex Count Result				CV Result			
	Target	Results	Limit	Status	Results	Limit	Status	
7µm SF_FSC_P	1st	143	± 50	Pass	1.7	< 2.8 %	Pass	
	2nd	143		Pass			Pass	
	3rd	143		Pass			Pass	
1µm CB_FSC_P	1st				Yes	Check for single peak	Pass	
	2nd				Yes		Pass	
	3rd				Yes		Pass	
2.5µm SF_FLH_P	1st	108	+ 60	Pass	3.1	< 7.5 %	Pass	
	2nd	108		Pass			Pass	
	3rd	107		Pass			Pass	
2.5µm SF_SSH_P	1st	102	± 60	Pass	4.8	< 7.5 %	Pass	
	2nd	102		Pass			4.5	Pass
	3rd	102		Pass			4.5	Pass
2.5µm SF_DSS_P	1st	128	± 60	Pass	Yes	Check for single peak	Pass	
	2nd	127		Pass			Yes	Pass
	3rd	128		Pass			Yes	Pass

SENSITIVITY ADJUSTMENT

4.1 Material Used : UF-Calibrator (P/N: CN383000)
 Lot No. : UA0077
 Expiry date : 16 October, 2020

4.2 UF Calibrator Assay Sheet UF Calibrator Assay Sheet

Table for UF Calibrator Assay value

Parameter	Target	Lower Range	Upper Range
Calibration Particle /uL	843.6	801.4	885.8
Conductivity ms/cm	34.9	32.8	37.0
Parameter	Lower Range	Upper Range	Target
SF_FSC_P	150.3	153.3	151.8
SF_FSC_W	43.2	44.2	43.7
CW_FSC_P	102.7	105.7	104.2
CB_FSC_P	163.2	166.2	164.7
SF_FLL_P	109.2	115.2	112.2
SF_FLH_P	109.2	115.2	112.2
SF_DSS_P	169.7	179.7	174.7
CW_SSH_P	108.1	114.1	111.1
CW_FLH_P	61.0	65.0	63
CW_DSS_P	21.2	23.2	22.2
CB_SSH_P	201.4	217.4	209.4
CB_FLL_P	105.5	125.5	115.5
SF_SSH_P	189.7	199.7	194.7
SF_SSL_P	189.7	199.7	194.7
CW_SSL_P	21.4	23.4	22.4
CW_FLL_P	61.0	65.0	63.0
CB_FLH_P	105.5	125.5	115.5

4.3 Sensitivity adjustment is done by executing [Sensitivity adjustment]. The results are shown here:

4.4 Go to sensitivity adjustment

Table 6 - result of sensitivity

Parameter SED [Ch]	U	Target	U/L	Results	Acceptable Range(±)	Status
SF_FSC_P		151.8		151.9	1.5	Pass
SF_FSC_W		43.7		44.0	0.5	Pass
CW_FSC_P		104.2		104.5	1.5	Pass
CB_FSC_P		164.7		164.8	1.5	Pass
SF_FLL_P		112.2		112.2	5.0	Pass
SF_FLH_P		112.2		111.9	3.0	Pass
SF_DSS_P		174.7		172.8	5.0	Pass
CW_SSH_P		111.1		110.4	3.0	Pass
CW_FLH_P		63.0		61.5	2.0	Pass
CW_DSS_P		22.2		21.8	1.0	Pass
CB_SSH_P		209.4		210.6	3.0	Pass
CB_FLL_P		115.5		114.0	10.0	Pass
SF_SSH_P		194.7		190.8	5.0	Pass
SF_SSL_P		194.7		192.0	3.0	Pass
CW_SSL_P		22.4		22.3	1.0	Pass
CW_FLL_P		63.0		61.8	2.0	Pass

* Due to average of the UF Calibrator target value was been round-up to one decimal point, adjustment final result will be base on the status result.

Analyzer UF 4000
 File No 1
 Lot No UK0083

QC Data
 Material UF-CONTROL-L
 Exp Date 2020/10/12

SYSTEM

Date Time	BACT	Cond
2020/09/12 20 18 00	218.5	6.4
2020/09/12 20 21 15	223.5	7.6
2020/09/12 20 23 28	230.5	7.9

MEAN	224.1	7.3
SD	6.02	0.79
CV (%)	2.6	10.8
UL	260.1	11.8
TARGET	204.7	8.4
LL	143.3	5.0
Unit	/uL	mS/cm

A3636

00-14_Bulldozer

2020/09/12 21 44 20

1/1

Table 11 - Result for Control QC - H

UF CONTROL		UF Control -H Assay Value				Status
Parameter	Result	Target	Limit	Min	Max	
RBC uL	188.0	206.2	41.2	165.0	247.4	Pass
WBC uL	835.3	842.2	168.4	673.8	1010.6	Pass
EC uL	77.5	76.6	38.3	38.3	114.9	Pass
CAST uL	15.58	18.79	9.40	9.39	28.19	Pass
BACT uL	801.8	808.9	202.2	606.7	1011.1	Pass
COND mS/cm	35.0	34.8	3.5	31.3	38.3	Pass
SF_FSC_P CH	152.3			136.5	166.9	Pass
SF_FSC_W CH	44.0			41.8	46.2	Pass
SF_FLH_P CH	111.7			80.1	148.7	Pass
SF_SSL_P CH	189.1			80.1	235.8	Pass
SF_DSS_P CH	170.0			141.5	212.3	Pass
CW_FSC_P CH	104.6			93.7	114.6	Pass
CW_FLH_P CH	93.9			64.6	124.2	Pass
CW_SSH_P CH	108.4			89.2	131.8	Pass
CW_SSL_P CH	23.7			19.5	29.3	Pass
CW_DSS_P CH	27.2			22.5	31.8	Pass
CB_FSC_P CH	164.5			147.5	180.3	Pass
CB_FLL_P CH	174.7			116.2	223.5	Pass
CB_SSH_P CH	205.5			178.2	241.1	Pass

Background Check

Sample No

BACKGROUNDCHECK0000001

Parameters	Blank	Limit	Unit
RBC	0.0	1.0	/ul
WBC	0.0	1.0	/ul
WBC Clumps	0.0	1.0	/ul
EC	0.0	1.0	/ul
Squar.TC	0.0	1.0	/ul
Roun.SEC	0.0	1.0	/ul
Tran.EC	0.0	1.0	/ul
RTEC	0.0	1.0	/ul
CAST	0.00	1.00	/ul
Hy.CAST	0.00	1.00	/ul
Path.CAST	0.00	1.00	/ul
BACT	0.0	5.0	/ul
X.TAL	0.0	1.0	/ul
YLC	0.1	1.0	/ul

Parameters	Blank	Limit	Unit
SPERM	0.0	1.0	/ul
TRICUS	0.00	1.00	/ul

Parameters	Blank	Limit	Unit
SP.TC	123	3000	COUNT
CM.TC	43	300	COUNT
CB.TC	103	3000	COUNT

Close

Processing



URINE

2 MIXING MOTOR VERIFICATION



2.1 Mixing Motor speed needs to be confirmed before executing the adjustment. Execute [Mixing Motor] adjustment and result is shown here:

Table 2 - Mixing motor speed adjustment result

Mixing motor type	Reading	Acceptable range	Status
SFch Reactor Mixing Motor	799	750-850	Pass
CRch Reactor Mixing Motor	800	750-850	Pass



- All laser adjustments result must be attached here:



Optical Axis Adjustment

Fine Adjustment



SF ch
TotalCount

1938

CR(BACT) ch
TotalCount


381

	SF_FSC_P	CB_FSC_P
MFV[ch]	143.0	33.0
CV[%]	17	31.9

HIST_LATEX_SF_FSC_P



HIST_LATEX_CB_FSC_P



F100-4000000

Optical Axis Adjustment

Fine Adjustment



SF ch
TotalCount


258

CR(BACT) ch
TotalCount

1794

	SF_FSC_P	CB_FSC_P
MFV[ch]	26.0	142.0
CV[%]	35.8	2.2

HIST_LATEX_SF_FSC_P



HIST_LATEX_CB_FSC_P



SF_FLH_P/
SF_SSH_P

Modify
Adj. Value

Analysis

Cancel

F100-4000000

Sensitivity Adjustment

SFch	CR(WBC)ch	CR(BACT)ch			
SF FSC P	Target	Gain			
151.9	151.8	1.375			
SF SSH P	Target	Gain	2.800	high	▼
190.8	194.7	0.800			
SF FLH P	Target	Impressed Voltage			
111.9	112.2	2675	high	▼	
SF DSS P	Target	Impressed Voltage			
172.8	174.7	2840			
SF FSC W	Target	Sheath Pressure			
44.0	43.7	0.107			

Analysis

Scattergram
Histogram

Confirm

Cancel

Sensitivity Adjustment



SFch	CR(WBC)ch	CR(BACT)ch			
SF FSC P	Target	Gain			
151.9	151.8	1.375			
SF SSH P	Target	Impressed Voltage			
192.0	194.7	1795	low	▼	
SF FLH P	Target	Gain	0.387	low	▼
112.2	112.2	4.800			
SF DSS P	Target	Impressed Voltage			
172.8	174.7	2840			
SF FSC W	Target	Sheath Pressure			
44.0	43.7	0.107			



Analysis

Scattergram
Histogram

Confirm

Cancel






CERTIFICATION

9.1 We certify that the UF Automated Urine Particle Analyzer S/N : **11929** has been successfully calibrated in accordance with the manufacturer's recommendations.

Report and Calibration Performed By :


Report Reviewed and Accepted By :



Signature (Engineer 1)

Name: Jakir Hussain Mondal

Date: 12 September, 2020



Signature (Customer)

Name: _____

Date: 14/9/2020

Dr. Nito Yepthomi
MBBS Pathology
Reg. No. NMCP-18/01006
Christian Institute of Health Sciences & Research
Dimapur, Nagaland

Signature (Engineer 1)

Name: _____

Date: _____

Company Stamp : (Vendor)

Company Stamp : (Customer)