

TULIP DIAGNOSTICS (P). LTD. (INSTRUMENT DIVISION)

BAMBOLIM, ST-INEZ & VERNA - GOA

Ref SOP No: SOP/ID/012
CC No: CC16/ID/003

Installation Report

| | |
|--|--|
| Name & Address of Customer: GI pathology lab Services. | Instrument Model: rounrel penta.2.0 |
| Contact Person: Mr. Gudappa Khodake. | Instrument Sr. No.: 0841230200217. |
| Contact No.: 9673521051. | Installation Date: 01/02/2021. |
| Power supply details: Instrument connected to: Online UPS Yes/No <input checked="" type="checkbox"/> Offline UPS Yes/No <input type="checkbox"/> | Invoice No.: |
| Supply input to machine from UPS: LN: 241 NE: 3V LE: 244. | Invoice Date: |
| Supply input to UPS from mains: LN: - NE: - LE: - | AMC Due Date: |
| Accessories Provided: All Accessories provided. | |
| Person Trained: Mr. Gudappa, Khodake | |
| Demonstrated: - yes - | |
| Remarks If any: - NO - | |
| Note: The above machine is satisfactory installed at the site and staff adequately trained. * Received all accessories as per packing list | |
| Name of Department/Head: Signature & Stamp: Date: | Company's Sales & Service Engineer: Signature: Sachin Dongre Date: 01/03/2021. |
| Gitaniali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim P.O. Goa. Pin: 403 202. Tel: 0832-2458546/51, Fax: 0832-2458544, E-mail: sales@tulipgroup.com, Website: www.tulipgroup.com | |

T TULIP DIAGNOSTICS (P) LTD.

Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim P.O. Goa. Pin: 403 202. Tel: 0832-2458546/51, Fax: 0832-2458544.
Email ID: ccgintd@tulipgroup.org. Website: www.tulipgroup.com

Warranty Certificate

| | | | |
|---------------------|--------------------------------------|--------------------|----------------|
| Customer Name: | G.I. Pathology Lab Services. | | |
| Address: | Nyan Apartment, Rajyog Society, PUNE | | |
| Contact Person: | 9099029851 | Contact No: | — |
| Instrument Model: | COUNCEL pent9.2.0 | Instrument Sr. No: | 0841230200217. |
| Invoice No. & Date: | | Installation Date: | 01/03/2021. |

Warranty Condition

- To activate the warranty, the customer must fill the installation report & warranty certificate and return the slip to Tulip Diagnostics (P) Ltd.
- Tulip Diagnostics (P) Ltd., guarantees that all its instruments are free from manufacturing defects or faults.
- Tulip undertakes repair or substitute free of charge replacement of spare part which may be found to have manufacturing defects.
- The warranty does not cover to defects of parts which are subject to wear & tear.
- Repair & interventions carried out during the period of the warranty do not extend or renew the period of warranty.
- The repairs of the instrument will be carried out by Tulip's authorized engineer only.
- Tulip reserves the right to recall the instrument for repair at the head office if major/frequent problem has been observed in the instrument.

Termination Of Warranty

The warranty shall be terminated at the end of the warranty period & also in the following cases:

- Where attempts to make repairs or alterations have been made by unauthorized person &/or with spare parts which are not originals.
- Alteration have been made to the serial number of the product on the certificate or on the instrument.
- The instrument is transferred to a new location, without prior written approval from Tulip Diagnostics (P) Ltd.

Validity & Duration

- This warranty shall be considered valid only on the condition that this certificate is accompanied by installation report and other purchase documents.
- The warranty is valid for a period of ___ months from the date of installation or ___ months from the date of invoice whichever is earlier.

Elements Not Covered Under This Warranty

The following damage & faults are not covered under this warranty:

- Damage deriving &/or originating from an insufficient or inadequate electric circuit or from the area where the instrument is set up & used.
- Breakdowns caused by careless handling, imprudence, lack of expertise & in any case caused by lack of skill or any degree of negligence on the part of the operator.
- Damage, defects & faults deriving from unexpected events, accidents during transport by the purchaser, due to FORCE MAJORE & in any case, due to situation which can in no way be attributed to manufacturing &/or material defects.
- Replacement of the defective instrument is never foreseen.
- Tulip shall accept no responsibility whatsoever for damage either directly or indirectly to persons or materials from the use of the instrument.
- Consumable items such as fuses, lamps, all tubing's, etc. are not covered under this warranty condition.

Governing Law

In case of any disputes the Courts of Panaji, Goa will have the jurisdictions in all matters.

Customer Signature & stamp

Name:

Date:



Company's Representative Signature & Stamp

Name:

Date:

A handwritten signature in black ink, appearing to read 'Sachin Dange', is written above the company's name and date.

Sachin Dange

01/03/2021.



QUALITY CONTROL SHEET

| | |
|--------------|--------------------|
| Description | COUNCELL PENTA 2-0 |
| Serial No. | 0841230200217 |
| QC Pass Date | 22/02/2021 |

| | |
|---------------|------------------------------|
| PART 1 | TEST APPARATUS |
| REAGENTS | DILUENT = 20200106 |
| | Differential Lyse = 20190909 |
| | LH Lyse = 20190927 |
| QC MATERIAL | CBC DH2101 |

| | | |
|--|---|---|
| PART 2 | QUALITATIVE TEST | |
| Circle which appropriate | | |
| 1. Chassis/ Housing – verify physical integrity cleanliness | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 14. HgB Lamp – Verify condition |
| 2. AC Plug - Verify integrity | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 15. Fittings / Connectors – examine all cables and connectors |
| 3. Mount/ Fastener – verify physical integrity | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 16. Audible Signals – Verify operation and level |
| 4. Strain Relief – verify physical integrity at both ends of line cord | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 17. Indicators / Display – Verify proper illumination and operation |
| 5. Line cord - verify proper insulation and integrity | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 18. Optics / Measuring Area- Inspect |
| 6. Circuit Breakers/ Fuses – Verify correct rating | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 19. Alarms – Check all warning alarms |
| 7. Cables – Verify physical integrity | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 20. Initialization Process – Verify |
| 8. Control / Switches – Verify operation | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 21. User Settings – Verify |
| 9. Hydraulics and Pneumatic System – Verify | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 22. Labeling – Verify caution and warning label |
| 10. Tube / Hoses – Verify condition | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 23. User Calibration – Verify |
| 11. Electrodes and Transducers – Verify condition and operation | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 24. Accessories - Verify as appropriate |
| 12. Probes - Verify Integrity | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | 25. Others (Specify) |
| 13. Fan / Motor / Compressor – Verify condition and operation | <input checked="" type="radio"/> PASS <input type="radio"/> FAIL <input type="radio"/> NA | |



| PART 3 | | QUANTITATIVE TASKS | |
|-----------|-------------------------|--------------------|------------------|
| A. | Gain Factor (%) | | Remark |
| LAS | 95 | | <u>PASS</u> FAIL |
| MAS | 155 | | <u>PASS</u> FAIL |
| WAS | 120 | | <u>PASS</u> FAIL |
| Width | 120 | | <u>PASS</u> FAIL |
| B. | Set Value | Blank Voltage | Remark |
| RBC | 39 | - | <u>PASS</u> FAIL |
| HgB | 85 | 4-52 | <u>PASS</u> FAIL |
| C. | Calibration Coefficient | | |
| Parameter | Coefficient | | Remark |
| WBC | 100 | | <u>PASS</u> FAIL |
| RBC | 100 | | <u>PASS</u> FAIL |
| HgB | 100 | | <u>PASS</u> FAIL |
| MCV | 100 | | <u>PASS</u> FAIL |
| PLT | 100 | | <u>PASS</u> FAIL |



| D. | Measured Values | | |
|----------------|------------------|---------------------------|-------------|
| NORMAL CONTROL | | | |
| Parameter | Obtained Value | Range | Remark |
| WBC | 8.27 $10^9/L$ | (7.21 - 9.21) $10^9/L$ | (PASS) FAIL |
| Neu | 4.64 $10^9/L$ | (3.87 - 5.67) $10^9/L$ | (PASS) FAIL |
| Lym | 2.36 $10^9/L$ | (1.56 - 2.96) $10^9/L$ | (PASS) FAIL |
| Mono | 0.46 $10^9/L$ | (0.00 - 0.90) $10^9/L$ | (PASS) FAIL |
| Eoso | 0.91 $10^9/L$ | (0.00 - 1.66) $10^9/L$ | (PASS) FAIL |
| Baso | 0.11 $10^9/L$ | (0.00 - 0.20) $10^9/L$ | (PASS) FAIL |
| Neu % | 0.561 % | (0.449 - 0.689) % | (PASS) FAIL |
| Lym % | 0.285 % | (0.196 - 0.356) % | (PASS) FAIL |
| Mono % | 0.056 % | (0.000 - 0.110) % | (PASS) FAIL |
| Eoso % | 0.098 % | (0.000 - 0.202) % | (PASS) FAIL |
| Baso % | 0.013 % | (0.000 - 0.024) % | (PASS) FAIL |
| RBC | 4.56 $10^{12}/L$ | (4.34 - 4.82) $10^{12}/L$ | (PASS) FAIL |
| HgB | 13.2 g/dL | (12.7 - 13.9) g/dL | (PASS) FAIL |
| MCV | 93.7 fL | (89.4 - 99.4) fL | (PASS) FAIL |
| PLT | 236 $10^9/L$ | (209 - 289) $10^9/L$ | (PASS) FAIL |



| D. | | Measured Values | |
|-------------|------------------|---------------------------|-------------|
| LOW CONTROL | | | |
| Parameter | Obtained Value | Range | Remark |
| WBC | 3.43 $10^9/L$ | (2.93 - 3.93) $10^9/L$ | (PASS) FAIL |
| Neu | 1.71 $10^9/L$ | (1.30 - 2.20) $10^9/L$ | (PASS) FAIL |
| Lym | 1.29 $10^9/L$ | (0.91 - 1.61) $10^9/L$ | (PASS) FAIL |
| Mono | 0.23 $10^9/L$ | (0.00 - 0.42) $10^9/L$ | (PASS) FAIL |
| Eoso | 0.20 $10^9/L$ | (0.00 - 0.42) $10^9/L$ | (PASS) FAIL |
| Baso | 0.06 $10^9/L$ | (0.00 - 0.14) $10^9/L$ | (PASS) FAIL |
| Neu % | 0.498 % | (0.388 - 0.628) % | (PASS) FAIL |
| Lym % | 0.375 % | (0.278 - 0.458) % | (PASS) FAIL |
| Mono % | 0.068 % | (0.000 - 0.126) % | (PASS) FAIL |
| Eoso % | 0.059 % | (0.000 - 0.122) % | (PASS) FAIL |
| Baso % | 0.017 % | (0.000 - 0.040) % | (PASS) FAIL |
| RBC | 2.45 $10^{12}/L$ | (2.26 - 2.62) $10^{12}/L$ | (PASS) FAIL |
| HgB | 5.8 g/dL | (5.4 - 6.2) g/dL | (PASS) FAIL |
| MCV | 77.5 fL | (72.4 - 82.4) fL | (PASS) FAIL |
| PLT | 49 $10^9/L$ | (30 - 70) $10^9/L$ | (PASS) FAIL |




| D. | Measured Values | | |
|--------------|------------------|---------------------------|-------------|
| HIGH CONTROL | | | |
| Parameter | Obtained Value | Range | Remark |
| WBC | 18.31 $10^9/L$ | (16.17 - 21.17) $10^9/L$ | (PASS) FAIL |
| Neu | 11.82 $10^9/L$ | (9.89 - 14.39) $10^9/L$ | (PASS) FAIL |
| Lym | 3.53 $10^9/L$ | (2.18 - 4.88) $10^9/L$ | (PASS) FAIL |
| Mono | 1.04 $10^9/L$ | (0.00 - 1.98) $10^9/L$ | (PASS) FAIL |
| Eoso | 1.92 $10^9/L$ | (0.00 - 4.00) $10^9/L$ | (PASS) FAIL |
| Baso | 0.18 $10^9/L$ | (0.00 - 0.40) $10^9/L$ | (PASS) FAIL |
| Neu % | 0.645 % | (0.530 - 0.770) % | (PASS) FAIL |
| Lym % | 0.193 % | (0.119 - 0.259) % | (PASS) FAIL |
| Mono % | 0.057 % | (0.000 - 0.106) % | (PASS) FAIL |
| Eoso % | 0.105 % | (0.000 - 0.214) % | (PASS) FAIL |
| Baso % | 0.010 % | (0.000 - 0.022) % | (PASS) FAIL |
| RBC | 5.18 $10^{12}/L$ | (4.88 - 5.48) $10^{12}/L$ | (PASS) FAIL |
| HgB | 16.4 g/dL | (15.8 - 17.4) g/dL | (PASS) FAIL |
| MCV | 102.2 fL | (97.8 - 107.8) fL | (PASS) FAIL |
| PLT | 469 $10^9/L$ | (409 - 529) $10^9/L$ | (PASS) FAIL |

Checked By:

Approved By:

File No. **1** Lot No. **DH2101** Level **Normal** Exp.Date **05/03/2021**
 Mode **Whole Blood** QC Type **CBC 5D** QC Sample ID **1**

| Para | Result | Unit | Para | Result | Unit |
|------|--------|--------------------|--------|--------|---------------------|
| WBC | 8.27 | 10 ⁹ /L | RBC | 4.56 | 10 ¹² /L |
| Neu# | 4.64 | 10 ⁹ /L | HGB | 13.2 | g/dL |
| Lym# | 2.36 | 10 ⁹ /L | HCT | 0.427 | fl |
| Mon# | 0.46 | 10 ⁹ /L | MCV | 93.7 | pg |
| Eos# | 0.81 | 10 ⁹ /L | MCH | 29.0 | g/dL |
| Bas# | 0.11 | 10 ⁹ /L | MCHC | 30.9 | fl |
| Neu% | 0.561 | | RDW-CV | 0.154 | |
| Lym% | 0.285 | | RDW-SD | 53.0 | |
| Mon% | 0.056 | | PLT | 236 | 10 ⁹ /L |
| Eos% | 0.098 | | MPV | 9.4 | fl |
| Bas% | 0.013 | | PDW | 16.0 | mL/L |
| | | | PCT | 2.20 | 10 ⁹ /L |
| | | | P-LCC | 62 | |
| | | | P-LCR | 0.264 | |



Setup
QC Table
QC Graph
Edit Result



Sample Analysis



Review



QC



Reagent Setup



Diluent




Print

File No. 2 Lot No. DH2101 Level Low Exp.Date 05/03/2021

Mode Whole Blood QC Type CBC 5D QC Sample ID 2

| Para | Result | Unit | Para | Result | Unit |
|------|--------|--------------------|--------|--------|---------------------|
| WBC | 3.43 | 10 ⁹ /L | RBC | 2.45 | 10 ¹² /L |
| Neu# | 1.71 | 10 ⁹ /L | HGB | 5.8 | g/dL |
| Lym# | 1.29 | 10 ⁹ /L | HCT | 0.190 | |
| Mon# | 0.23 | 10 ⁹ /L | MCV | 77.5 | fL |
| Eos# | 0.20 | 10 ⁹ /L | MCH | 23.6 | pg |
| Bas# | 0.06 | 10 ⁹ /L | MCHC | 30.4 | g/dL |
| Neu% | 0.498 | | RDW-CV | 0.173 | |
| Lym% | 0.375 | | RDW-SD | 49.0 | fL |
| Mon% | 0.068 | | PLT | 49 | 10 ⁹ /L |
| Eos% | 0.059 | | MPV | 10.0 | fL |
| Bas% | 0.017 | | PDW | 16.2 | mL/L |
| | | | PCT | 0.49 | mL/L |
| | | | P-LCC | 15 | 10 ⁹ /L |
| | | | P-LCR | 0.304 | |



Setup

QC Table

QC Graph


Edit Result



Mode:Whole Blood-CBC+DIFF

Service

11:54 22/02/2021




Sample Analysis



Review



QC



Reagent Setup



Diluent




Print

File No. 3 Lot No. DH2101 Level High Exp.Date 05/03/2021

Mode Whole Blood QC Type CBC 5D QC Sample ID 3


| Para | Result | Unit | Para | Result | Unit |
|------|--------|--------------------|--------|--------|---------------------|
| WBC | 18.31 | 10 ⁹ /L | RBC | 5.18 | 10 ¹² /L |
| Neu# | 11.82 | 10 ⁹ /L | HGB | 16.4 | g/dL |
| Lym# | 3.53 | 10 ⁹ /L | HCT | 0.529 | |
| Mon# | 1.04 | 10 ⁹ /L | MCV | 102.2 | fL |
| Eos# | 1.92 | 10 ⁹ /L | MCH | 31.6 | pg |
| Bas# | 0.18 | 10 ⁹ /L | MCHC | 30.9 | g/dL |
| Neu% | 0.645 | | RDW-CV | 0.141 | |
| Lym% | 0.193 | | RDW-SD | 53.4 | fL |
| Mon% | 0.057 | | PLT | 469 | 10 ⁹ /L |
| Eos% | 0.105 | | MPV | 9.6 | fL |
| Bas% | 0.010 | | PDW | 15.9 | |
| | | | PCT | 4.52 | mL/L |
| | | | P-LCC | 132 | 10 ⁹ /L |
| | | | P-LCR | 0.281 | |



DIFF
LAS MAS



WBC
LAS MAS



RBC
LAS MAS




PLT
LAS MAS

QC Table

QC Graph

Setup

Edit Result



Mode:Whole Blood-CBC+DIFF

12:03
22/02/2021

Installation Qualification
For
Councill Penta 2.0
Auto hematology Analyzer



Contents

Instrument Identification

Environment Condition

Instrument Delivery and Documentation

Instrument Safety

Assembly and Installation

Summary Report

Instrument Identification

Auto hematology analyzer (Councill penta 2.0)

| Sr. No. | Specification | As per design qualification | Specification received | Match (Yes/No) |
|---------|---------------------------|--|--|----------------|
| 1 | General Feature | The measurement methods used in this analyzer are: the electrical Impedance method for determining the RBC and PLT data; the colorimetric method for determining the HGB; laser-based flow cytometry for determining the WBC data. | Electrical. Impedance method. for RBC and PLT. Colorimetric. method for HGB. laser. based flow. cytometry for. WBC data. | YES. |
| 2 | Sample aspiration volume. | Whole Blood mode 20 μ L, pri-diluted mode 20 μ L of whole blood/capillary blood sample and 480 μ L of diluent | whole blood 20 μ L. pri-diluted 20 μ L. capillary/whole blood = 20 μ L. Diluent - 480 μ L. | YES. |
| 3 | LED lamp | 530nm LED | 530nm LED | YES. |
| 4 | Communication port | LAN Port | LAN Port | YES. |
| 5 | Printer | External ink tank and leaser printer | External ink. tank & laser printer. | YES. |
| 6 | Power requirement | AC100V~240V, Input Power \leq 200VA, 50/60HZ. | AC 100V~240V. Input power \leq 200VA | YES |
| 7 | Grounding | Good grounding <5V | Good grounding. 2.5V. | YES. |

Comments:

Environment Condition

Has the Instrument been adequately acclimatized since transport or storage?

YES:

NO:

| Operating Parameters | Specified Range | Conditions Met |
|----------------------------|--------------------------|--------------------------------|
| Ambient Temperature | 10 to 30 degrees | 27°C |
| Relative Humidity | ≤85% | YES |
| Supply Voltage & Frequency | 100 to 240V, 50/60 Hz AC | 42 41/5 244, NE 3V. |

Comments:

Instrument is in good condition -

Instrument Delivery and Documentation

Unpack the Instrument carefully. Are any items missing against the Packing List?

YES:

NO: ✓

If yes, state the missing Items: _____

Is there any damage to the Instrument or Accessories?

YES:

NO: ✓

If yes: Description of Damage: _____

Corrective Action: _____

Manufacturing informed: YES: NO:

Is all Standard Documentation Included?

Operating Manual: Present: ✓ Missing:

Warranty Card: Present: ✓ Missing:

Is all additional documentation Included:

Present: ✓ Missing:

Installation Qualification:

Operational Qualification:

Performance Qualification:

Comments:

All documents are present with instrument.


Instrument Safety

Manufacturer Safety Recommendations:

| | | |
|--|--------|-----|
| Allow enough space for the Instrument | YES: ✓ | NO: |
| Position on level, non-combustible surface Indoors | YES: ✓ | NO: |
| Area free from aggressive/ explosive chemical mixtures | YES: ✓ | NO: |
| Ensure correct power supply and all power switches Easily accessible | YES: ✓ | NO: |
| Keep the analyzer away from electromagnetic waves | YES: ✓ | NO: |
| Read Operator Manual before use | YES: ✓ | NO: |

Comments:

followed all the safety instructions.



Assembly and Installation

| Assembled and installed By: | User: | Specialized Engineer: |
|--|-------|-----------------------|
| Installation Procedure: N/A | | OK |
| Unpack and retain Packaging | | |
| Verifying the site before connecting to UPS power supply | | |
| Checking the instrument before performing the Test | | |

Comments:

All checks completed.

Summary Report

Instrument Name: *Coucell Penta. 210.*

Sr. No.: *0841230200217.*

Manufacturer: *Tulip Diagnostics Pvt Ltd [Color clinical system]*

Assessment of complete Installation Qualification:

No Deviations:

Deviations:


| Deviation | Impact on Operation | Justification for Acceptance |
|-----------|---------------------|------------------------------|
| | | |
| | | |
| | | |

Successful completion of the preceding activities and checks indicates that this Instrument is installed successfully as per the set protocol. The instrument has passed the Installation Qualification procedure and may now be released for use.

| | | | |
|--------------------------------|--------------------------------|--------------------------------|--|
| IQ Completed By: | <i>Mr. Sachin. Danghe.</i> | Date: <i>01/03/2021</i> | |
| Deviations Approved By: | | Date: | |
| IQ Approved By: | | Date: | |


Comments:

Installation completed.



Tulip Diagnostics (P) Ltd.

Operation Qualification
For
Council penta 2.0
Auto Hematology Analyzer.



Contents

Safety Tests

Pre-Run Checks

Functional Tests

Training

Summary Report.

Safety Tests

This instrument confirms to the following standards and is supplied fully tested by the Manufacturer.

Electrical equipment for measurement, control and laboratory use.
Instrument requirement for laboratory for Whole Blood testing.

The following electrical safety tests should be carried out by a competent person after Installation.

| Safety Tests | Pass | Fail | Date | Signed |
|----------------------|-------|------|------------|---------------|
| Earthing voltage | Pass. | | 01/03/2021 | Sachin Wangle |
| Input Supply voltage | Pass | | 01/03/2021 | Sachin Wangle |

Comments:

safety test pass-

Pre-Run Checks



Auto Hematology Analyzer (Councell Penta 2.0)

| Pre run check | OK | N/A |
|--------------------------|----|-----|
| Mechanics initialization | OK | |
| Fluidic initialization | OK | |
| Ambient temperature | OK | |
| HGB AD value | OK | |
| Background test | OK | |

Comments:

Functional Tests

Auto Hematology Analyzer (Councell Penta 2.0)

| Functional Test | Acceptance Criteria | Pass | Fail |
|---------------------|--|-------|------|
| Instrument power ON | <p>Log on the Auto hematology analyzer</p>  <p>After self test the instrument will display the following screen.</p>  | Pass | |
| Self Test | The instrument will also prompt the error status on the screen if any. | Pass. | |
| Background Test | $WBC \leq 0.2 \times 10^9 / L$, $RBC \leq 0.02 \times 10^{12} / L$, $HGB \leq 1g/L$, $PLT \leq 10 \times 10^9 / L$, $HCT \leq 0.5\%$ | Pass. | |

Comments:

Training

| Trained Operators | Read Operator Manual | Read SOP | Practical Training | Authorizing Signature |
|-------------------|----------------------|----------|--------------------|-----------------------|
| Suraj Rathod | YES. | YES. | YES. | Sachin Dangle |

| | | | | |
|--------------------|-----|-----|-----|---------------|
| Shalundkamble | YES | YES | YES | Sachin Dangle |
| Gundappa Khadke | YES | YES | YES | Sachin Dangle |
| Monika Veer | YES | YES | YES | Sachin Dangle |
| Ms. Kajal | YES | YES | YES | Sachin Dangle |
| Mr. Omkar | YES | YES | YES | Sachin Dangle |
| Omkarharihar | YES | YES | YES | Sachin Dangle |
| Prerajal pharbhade | YES | YES | YES | Sachin Dangle |
| Akshay Rathod | YES | YES | YES | Sachin Dangle |
| | | | | |

Comments:

Training completed -

Summary Report

Instrument Name: Councell Penta, 2.0

Sr. No.: 0841230200217

Manufacturer: Tulip Diagnostic Pvt Ltd.

Assessment of complete Operational Qualification:

No Deviations: ✓

Deviations:

| Deviation | Impact on Operation | Justification for Acceptance |
|-----------|---------------------|------------------------------|
| | | |
| | | |
| | | |
| | | |

Successful completion of the preceding activities and checks indicates that this Instrument is operating satisfactorily following delivery and installation. The instrument has passed the Operational Qualification procedure and may now be released for use.

| | | | |
|--------------------------------|---------------|--------------|-------------|
| OQ Completed By: | Sachin Dangle | Date: | 01/03/2021. |
| Deviations Approved By: | | Date: | |
| OQ Approved By: | | Date: | |


Comments:

Operation Qualification completed.



**PERFORMANCE QUALIFICATION
FOR
COUNCELL PENTA 2.0**

AUTO HEMATOLOGY ANALYZER



Contents

- 1) Performance Check
- 2) Routine Maintenance
- 3) Safety Tests
- 4) Summary Report

1) Performance Check

Auto Hematology Analyzer (Councell penta 2.0)

Test 1

Test Name: Sample Processing

Purpose: Ability to Process Samples

Method:

Run the control samples five times consecutively

Acceptance Criteria: Each of the results obtained above should be within the range as specified in the control chart.

Parameters Values for Verification:

Level 1 (Low Control)

RBC Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|---------------------------|------------------|------|------|
| 1. | (2.26 - 2.62) $10^{12}/L$ | 2.45 $10^{12}/L$ | Pass | |
| 2. | | | | |
| 3. | | | | |

WBC Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|------------------------|------------------|-------|------|
| 1. | (2.93 - 3.93) $10^9/L$ | 3.43 $10^9/L$ | Pass. | |
| 2. | | | | |
| 3. | | | | |

Hemoglobin:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|------------------|------------------|------|------|
| 1. | (5.4 - 6.2) g/dL | 5.8 g/dL | Pass | |
| 2. | | | | |
| 3. | | | | |

MCV:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|------------------|------------------|------|------|
| 1. | (72.4 - 82.4) fL | 77.5 fL | Pass | |
| 2. | | | | |
| 3. | | | | |

Platelet Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|------------------------------|-----------------------|------|------|
| 1. | (30 - 70) 10 ⁹ /L | 49 10 ⁹ /L | Pass | |
| 2. | | | | |
| 3. | | | | |

Level 2 (Normal Control)

RBC Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|-----------------------------------|----------------------------|------|------|
| 1. | (4.34 - 4.82) 10 ¹² /L | 4.56 x 10 ¹² /L | Pass | |
| 2. | | | | |
| 3. | | | | |

WBC Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|----------------------------------|-------------------------|------|------|
| 1. | (7.21 - 9.21) 10 ⁹ /L | 8.27 10 ⁹ /L | Pass | |
| 2. | | | | |
| 3. | | | | |

Hemoglobin:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|--------------------|------------------|------|------|
| 1. | (12.7 - 13.9) g/dL | 13.2 g/dL | Pass | |
| 2. | | | | |
| 3. | | | | |

MCV:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|------------------|------------------|------|------|
| 1. | (89.4 - 99.4) fL | 93.7 fL | Pass | |
| 2. | | | | |
| 3. | | | | |

Platelet Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|--------------------------------|------------------------|------|------|
| 1. | (209 - 289) 10 ⁹ /L | 236 10 ⁹ /L | Pass | |
| 2. | | | | |
| 3. | | | | |

Level 3 (High Control)

RBC Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|---------------------------|------------------|------|------|
| 1. | (4.88 - 5.48) $10^{12}/L$ | 5.18 $10^{12}/L$ | Pass | |
| 2. | | | | |
| 3. | | | | |

WBC Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|--------------------------|------------------|------|------|
| 1. | (16.17 - 21.17) $10^9/L$ | 18.31 $10^9/L$ | Pass | |
| 2. | | | | |
| 3. | | | | |

Hemoglobin:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|--------------------|------------------|------|------|
| 1. | (15.8 - 17.4) g/dL | 16.4 g/dL | Pass | |
| 2. | | | | |
| 3. | | | | |

MCV:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|-------------------|------------------|------|------|
| 1. | (97.8 - 107.8) fL | 102.2 fL | Pass | |
| 2. | | | | |
| 3. | | | | |

Platelet Count:

| Test | Control Values | Results Obtained | Pass | Fail |
|------|----------------------|------------------|------|------|
| 1. | (409 - 529) $10^9/L$ | 469 $10^9/L$ | Pass | |
| 2. | | | | |
| 3. | | | | |


Validation Team: Instrument division,

Name: Sachin Dangle

Designation: Senior sales & service Engineer,

Signature Sachin Dangle,

Date 05/03/2021.



Test 2

Test Name:

Tests for checking the performance of the instruments during analysis
Tests for checking long term performance of the instrument

Purpose:

The purpose of the above checks is to ensure the reliability of the results being obtained.

Method:

1. During Sample analysis:

To run control samples each time the instrument is used for sample analysis and verification of the results of the controls to be within the reference range to be established by performance of the precision experiments.

2. Long term Performance

This is to be checked by Levy Jennings plots to be updated once in six months


Validation Team: *Instrument division,*

Name *Sachin Dongke*

Designation *Senior Sales & Service Engineer -*

Signature *Sadhu Deyke*

Date *05/03/2021*



Routine Maintenance

Using probe cleaner daily

The cleanser soak should be performed:

- 1 During the installation and shut down of the machine
- 2 When the background results' abnormal , bad differential of scattergram or clogging.
- 3 Analyzer has been running for more than 24 hours.
- 4 More than 1000 samples were run within 24 hours.

2)Swab Cleaning and Maintenance

Repeatedly scrub the bottom end of the swab and the inside of its lower aperture with a Q-tip dampened with the diluted probe cleanser

Washing.

- 1.High background of WBC or HGB, clean WBC chamber
2. High background of RBC or PLT, clean RBC chamber
3. High background of scatter gram, clean Flowcell.
4. Clean Sampling probe when it's dirty

.Maintenance for other parts

Clean the dirty area with soft cloth with alcohol

Validation Team : *Instrument division*


Name: *Ravi Annigeli*

Designation: *Sales Offices.*

Signature

R. Annigeli

Tulip Diagnostics (P) Ltd.



Date:

3) Safety Tests

This instrument confirms to the following standards and is supplied fully tested by the Manufacturer.
Electrical equipment for measurement, control and laboratory use.
Instrument requirement for laboratory for Blood testing.

The following electrical safety tests should be carried out by a competent person after Installation.

| Safety Tests | Pass | Fail | Date | Signed |
|--------------|------|------|------------|---------------|
| Earthing | Pass | | 01/03/2021 | Sachin Donghe |
| UPS Supply | Pass | | 01/03/2021 | Sachin Donghe |

Comments:

4) Summary Report

Instrument Name: *concell Penta 2+0.*

Sr. No.: *0841230200217*

Manufacturer: *Tulip Diagnostic Pvt Ltd.*

Assessment of complete Performance Qualification:

No Deviations:

Deviations:

| Deviation | Impact on Operation | Justification for Acceptance |
|--------------|---------------------|------------------------------|
| | | |
| | | |
| | | |

Successful completion of the preceding activities and checks indicates that this Instrument is operating satisfactorily following delivery and installation. The instrument has passed the Operational Qualification procedure and may now be released for use.

| | | | |
|-------------------------|----------------------|-------|-------------------|
| PQ Completed By: | <i>Sachin Darghe</i> | Date: | <i>01/03/2021</i> |
| Deviations Approved By: | | Date: | |
| PQ Approved By: | | Date: | |

Comments:

control & safety test pass.

