

Date: 06-08-2021  
Effective Date: 06-08-2021

### Certificate of Calibration

**Customer Name: "MAEER'S VishwaRaj Hospital Laboratory, Loni Kalbhor, Pune"**

**Model : Automated Hematology Analyzer Sysmex XS-800i**

**Serial No. : 67641**

**Calibration Done Date: 06-08-2021**


**Next Calibration Due Date On or Before: 06-08-2022**

**Lab In-charge: . Dr. Aniruddha Garud**

*This is to certify that the above-mentioned product has been verified of calibration for CBC 6 parameters (WBC, RBC, HGB, HCT, MCV and PLT).*

Calibration at site performed by  
Er Name : Er. Khurshid M. Kazi  
Designation : Technical Director  
Soha Biomedical Engineers, Pune

Signature



Encl:

1. Assay Sheet of Calibrator.
2. Printouts


Date: 06-08-2021  
Effective Date: 06-08-2021

## Certificate of Inspection

1. Model: Automated Hematology Analyzer Sysmex XS-800i
2. Serial No.: 67641
3. Calibration Date: 06-08-2021
4. Material used: SCS-1000 (Lot No. 11940525, Expiry date: 15-Aug-2021)

Calibration for CBC 6 parameters using the measurement standard material (SCS-1000) was completed. The calibration result of 5 runs is summarized in the following table. Please refer to the attached sheets for the details.

*Chunika*




Soha Biomedical Engineers, Pune

### 5. BACKGROUND CHECK

PARAMETER	RESULT	Range
WBC	0.0	$0.3 \times 10^3 / \mu\text{L}$ or less
RBC	0.00	$0.02 \times 10^6 / \mu\text{L}$ or less
HGB	0.0	0.1 g/dL or less
PLT	4	$10 \times 10^3 / \mu\text{L}$ or less

*Shreshth*



SOHA BIOMEDICAL ENGINEERS  
PUNE



**6. PRECISION STUDY PERFORMED ON THE ANALYZER USING A BLOOD SAMPLE  
(ORIGINALS ATTACHED)**

SMP NO	WBC	RBC	HGB	HCT	MCV	PLT
P01	5.55	5.20	13.00	39.8	76.5	292
P02	5.6	5.29	13.00	40.40	76.4	285
P03	5.66	5.19	12.80	39.7	76.5	293
P04	5.87	5.22	12.90	39.8	76.2	285
P05	5.56	5.19	12.80	39.6	76.3	287
P06	5.72	5.24	12.80	40.0	76.3	281
P07	5.58	5.26	12.80	40.0	76.0	295
P08	5.62	5.22	12.80	39.7	76.1	291
P09	5.65	5.30	13.00	40.2	75.8	290
P10	5.67	5.29	13.00	40.1	75.8	296
Mean	5.65	5.24	12.89	39.93	76.19	289.50
SD	0.094	0.043	0.099	0.254	0.260	4.859
CV%	<b>1.671</b>	<b>0.815</b>	<b>0.771</b>	<b>0.636</b>	<b>0.341</b>	<b>1.678</b>
Acceptable CV%	Within 3.0%	Within 1.5%	Within 1.5%	Within 1.5%	Within 1.5%	Within 4.0%
Result	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>





**7. CALIBRATION DATA**

SMP NO/TIME	WBC	RBC	HGB	HCT	MCV	PLT
C00	7.22	4.51	12.4	37.0	82.0	248
C01	7.15	4.53	12.4	36.2	79.9	245
C02	7.53	4.51	12.4	36.1	80.0	235
C03	7.14	4.54	12.5	36.3	80.0	240
C04	7.19	4.58	12.5	36.7	80.1	246
<b>MEAN</b>	<b>7.25</b>	<b>4.534</b>	<b>12.44</b>	<b>36.46</b>	<b>80.40</b>	<b>242.8</b>
Acceptable Limits	7.046 - 7.626	4.499 - 4.639	12.29 - 12.48	35.79 - 37.40	79.08 - 81.11	240.0 - 257.6
<b>Result</b>	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>	<b>PASS</b>

*Phunji*



Sysmex Calibrator System Assay Sheet

11940525  
15-Aug-2021

LOT

Parameter	XE-Series		XT-Series		XS-Series* ✓	
	Assay Target	Acceptable Limits	Assay Target	Acceptable Limits	Assay Target	Acceptable Limits
WBC K/uL	7.752	7.446 - 8.059	7.979	7.664 - 8.294	7.336	7.046 - 7.626
RBC M/uL	4.587	4.516 - 4.657	4.482	4.413 - 4.551	4.569	4.499 - 4.639
HGB g/dL	12.49	12.40 - 12.59	12.30	12.21 - 12.40	12.39	12.29 - 12.48
HCT %	35.94	35.14 - 36.73	34.47	33.71 - 35.23	36.60	35.79 - 37.40
MCV fL	78.35	77.35 - 79.34	76.90	75.92 - 77.88	80.09	79.08 - 81.11
PLT K/uL	248.9	240.1 - 257.7	247.7	239.0 - 256.5	248.8	240.0 - 257.6

Parameter	K-4500 / K-1000 / K-800		pocH-100i**		KX-21		XP-Series	
	Assay Target	Acceptable Limits	Assay Target	Acceptable Limits	Assay Target	Acceptable Limits	Assay Target	Acceptable Limits
WBC K/uL	7.76	7.42 - 8.09	7.46	7.14 - 7.78	7.72	7.39 - 8.05	7.12	6.82 - 7.43
RBC M/uL	4.530	4.440 - 4.621	4.565	4.474 - 4.657	4.524	4.433 - 4.614	4.485	4.395 - 4.575
HGB g/dL	12.39	12.26 - 12.51	12.10	11.98 - 12.22	12.50	12.38 - 12.63	12.07	11.95 - 12.19
HCT %	33.34	32.60 - 34.07	35.64	34.85 - 36.43	33.87	33.12 - 34.62	33.69	32.94 - 34.43
MCV fL	73.59	72.78 - 74.40	78.07	77.21 - 78.92	74.87	74.05 - 75.69	75.10	74.28 - 75.93
PLT K/uL	256.0	243.2 - 268.8	244.9	232.6 - 257.1	275.9	262.1 - 289.7	268.5	255.0 - 281.9

SCS-1000 ASSAY TERM DEFINED

Assay Target – This is the assigned value for calibration.

Acceptable Limits – These limits represent the interval around the Assay Target that can be attributed to the expanded uncertainty of the total traceability chain.  
A calibrator mean (n=5) that falls within these limits indicates an accurately calibrated instrument.

\* XS-1000i/XS-800i – Assay target for WBC only for operation in CBC+Diff mode

\*\* pocH-100i – Assay Target for WBC only for systems operating under software version 00-18 and following



## SYSMEX XS-800I AUTOMATED HEAMATOLOGY ANALYZER

### INSTALLATION QUALIFICATION

For

“MAEER’s Vishwaraj Hospital Laboratory”

Marketed by:  
Transasia Bio-Medicals Ltd.,  
Transasia House,  
Chandivali Studio road,  
Andheri (E),  
MUMBAI – 400 072

1 | Page





UNMATCHED SERVICE

SINCE 1979...  
Page No.

## Table of Contents

Sr. No.	Contents	Page No.
I	Approval of the IQ procedure	3
II	Instructions	4
III	Scope	5
IV	Ancillary Information	6
V	Installation Qualification	8
VI	Installation Procedure	11
VII	Comments	13
VIII	System Certification	14




## I. Approval of the IQ procedure:

MAEER's Vishwaraj Hospital Laboratory and Transasia are jointly responsible for the installation of the system SYSMEX – HEMATOLOGY Analyzer, Model: XS-800I, Serial No. 67641 in the clinical lab of MAEER's Vishwaraj Hospital Laboratory as per the attached protocol.

### Protocol Performed By: Transasia Representative

Name : Mr. Swapnil Kamarikar  
Title : INSTALLATIONQUALIFICATION  
Company : TRANSASIA BIO-MEDICALS LTD.

Signature: 

Date: 6/8/21




### Validation Team from \_\_\_\_\_:

Name :  
Designation :  
Department :

### Customer Authorizations:

Name : Dr. Avinash Patel  
Title : INSTALLATION QUALIFICATION  
Site :

Signature: 

Date:



## II. Instructions

1. This document is to be completed at the time the system is shifted to its current location (new) and set up for operation.
2. An authorized TRANSASIA representative will check the system and enter the specific data as outlined in the appropriate Installation Qualification. Each result will be noted and dated.
3. Employee of MAEER's Vishwaraj Hospital Laboratory will verify each result and sign in the last page. The members of the validation team will carry this out.
4. ALL deviations from normal specification to include any problems with installation will be noted under COMMENTS. All resolution to such problems will also be noted in the COMMENTS section. Additional space is provided at the end of this protocol for the same.
5. This document contains proprietary information and is in no way to be copied, photographed or duplicated in any way without expressed written authorization by the Transasia Bio-Medicals Ltd., Transasia House, Mumbai.

Validation Team:

Name *Swapnil Kameniker*

Designation *Service Eng.*

Signature *[Handwritten Signature]*

Date *6/8/21*







### III. Scope

This Installation Qualification protocol will be performed on the SYSMEX-Hematology Analyzer, Model XS-800I Serial No. 67641 located in MAEER's Vishwaraj Hospital Laboratory. This Protocol will define the documentation that will be used to evaluate the instruments installation in accordance with the manufacture's specifications and intended use. Successful completion of this protocol will verify that the instrument identified has been installed in accordance with the intended usage.

Installation checks will also be performed to verify that the instrument has been installed with proper connections and utilities.


Trained, knowledgeable personnel will perform qualification studies.

Any exceptional conditions encountered during the qualification studies will be identified for review. Exceptional conditions will be investigated and the appropriate course of action determined. All documents will be initialed and dated.

Validation Team:

Name: Swapnil Kamenikar

Designation: Service Eng

Signature: 

Date: 6/8/21





## IV. Ancillary Information.

### a. Certification of Purchase Order Compliance

I certify to the best of my knowledge, the instrument is purchased under Purchase order No. \_\_\_\_\_, Dt. \_\_\_\_\_ sent against Quotation number \_\_\_\_\_ dt. \_\_\_\_\_ is in compliance with the specifications of the Purchase order.

Verified By : \_\_\_\_\_ Date : \_\_\_\_\_

### b. Utilities

Sr.No.	Utility	Yes / No	Verified By	Date
1.	Environmental condition as per requirement: (Ambient range of temperature 15 – 35 °C, air conditioning facility, non exposure to direct sunlight, non-interference from high frequency radio waves)	Yes / No	✓	
2.	Adequate space for installation : (Minimum in mm. W 450 X D 660 X H 450)	Yes / No	✓	
3.	Cellpack DCL, SULFOLYSER,,Lysercell WDF, Fluorocell WDF and Cell Clean	Yes / No	✓	
4.	Power Source Requirements*	Yes / No	✓	

\* Encircle applicable source

Validation Team :

Name Sneha Kankar

Designation Service Eng-  
Kans

Signature

Date 6/8/21



6 | Page



**c. The instrument has been verified for the following**

Sr.No.	Verification	Yes / No	Verified By	Date
1.	Instrument is identified	Yes / No	✓	
2.	Manufacturer's specifications are included	Yes / No	✓	
3.	Accessories / Consumables are listed	Yes / No	✓	
4.	Manufacturer's certificate of Compliance attached	Yes / No	✓	

Validation Team:

Name *Suresh Kumar*

Designation *Service Eng*

Signature *[Signature]*

Date *6/8/21*







## V. Installation Qualification

### A. Equipment Description

This Sysmex XS-800I is a fully automated Hematology analyzer for in vitro diagnostic use in clinical laboratories. The XS-800I provides accurate and precise test results for ( ) parameters.

Instrument identification		Verified by	Date
Equipment Name	Automated Hematology	✓	
Model	XS800i	✓	
Manufacturer	Sysmex Corporation	✓	
Marketed By	Transasia	✓	
Equipment #		✓	
Serial Number	67641	✓	
Size (in mm)	W 450 X D 660 X H 450	✓	
Power	AC 220 V	✓	
Frequency	50 - 60 Hz	✓	
Power Consumption	Less Than 250 VA	✓	

Validation Team:

Name *Suresh Kumar*

Designation *Service Eng*

Signature *[Signature]*

Date *6/8/21*



8 | Page





## Consumables:

Consumables such as Cellpack SULFOLYSER, Stomatlyser 4DL, Stroomtolyser 4DS and Cell Clean were supplied along with instrument.

Currently a sufficient stock of the same is being maintained Yes  No

## C. List of Manuals, Certificates and Drawings

Transasia provides the following with the instrument.

1. Operator's Manual

## D. Change Control Procedure

The instrument will not be altered, enhanced, modified or substituted for another system until a formal Change Control Authorization is approved from Transasia Bio-Medicals Ltd. and MAEER's Vishwaraj Hospital Laboratory.

## E. Maintenance

The instrument listed within this document will be placed under the control of the purchasing institution with respect to proper maintenance procedures as detailed in the operations manual Chapter 13

A trained analyst using the manuals provided with the instrumentation can perform simple maintenance. Upon expiration of the warranty period Transasia offers several levels of Maintenance Agreements and Performance Testing services to assist you in maintaining **GLP/GMP** compliance. Contacting your local representative and requesting the additional Service Agreement can supply additional information.

Validation Team:

Name *Swagunil Komeriker*

Designation *Service Eng*

Signature *[Signature]*

Date *6/8/21*







## F. Spare Parts

Transasia strongly recommends the end user maintain a basic of consumable parts on site to minimize down time due to minor failures. They have provided a list of such consumable parts and the same is also available in the Operator's Manual no.

## C. Equipment Logs

Title	Location	Verified by	Date

Sample page of the logbook is attached to this document

Effective date:

Validation Team:

Name *Sugan' Hemantkar.*

Designation *Service Eng.*

Signature *[Signature]*

Date *08/21*





## VI. Installation Procedure

(These had been performed at the time of original installation at the initial location)

1. **Check Before Installation**  
Refer to Chapter 7 of Sysmex XS800i Service Manual
3. **Grounding**  
Refer to Chapter 1 of Sysmex XS800i Service Manual
4. **Installation Environment & Space**  
Refer to Chapter 1 of Sysmex XS800i Service Manual
5. **Connect Air & Reagent Tubes**  
Refer to Chapter 7 of Sysmex XS800i Service Manual
6. **Connect Connection Cord & Power Cord**  
Refer to Chapter 7 of Sysmex XS800i Service Manual
7. **Turn Power On**  
Refer to Chapter 7 of Sysmex XS800i Service Manual

Validation Team:

Name *Sudipil Kanchan*  
Designation *Service Eng*  
Signature *[Signature]*  
Date *6/8/21*





UNMATCHED SERVICE  
SINCE 1979...

## VII. COMMENTS:

Validation Team:

Name *Sugoni Kamanikar*

Designation *Kan*

Signature

Date *6/8/21*







UNMATCHED SERVICE

SINCE 1979

## VIII. System Certification


Study data has determined that the system described in this document either meets all criteria outlined in this independently Installation Qualification Protocol, or exceptional conditions have been identified and documentation included. Exceptional conditions, if any, have been addressed. The system is ready for specified usage.

Report Performed By : Transasia Representative

Name : Mr. Swapnil Kamarikar

Title : INSTALLATION QUALIFICATION Signature: 

Company: TRANSASIA BIO-MEDICALS LTD. Date : 6/8/21



Customer Authorizations:

Name : *Dr. Anindha Kameel*

Title : INSTALLATION QUALIFICATION Signature:

Site : Date :

Name : *Dr. Anindha Kameel*

Title : INSTALLATION QUALIFICATION Signature: 

Site : Date :



Date: 6/8/21

Reagent Check done OK

Printer checked ✓

Analyzer switched ON at ✓

SELF CHECK performed ✓

RINSE CYCLE completed ✓

Background limits within acceptable range ✓

Analysis start time

Analysis end time

No. of samples analyzed

Shut down procedure done ✓

Analyzer switched OFF at ✓

Recorded by: Mr. Kiran Pandey Checked by: Dr. Anirudh Gaur

Date: 6/8/21



## SYMEX XS800i AUTOMATED HEMATOLOGY ANALYZER

### OPERATIONAL QUALIFICATION

For

“MAEER’s Vishwaraj Hospital Laboratory”

Marketed by:  
Transasia Bio-Medicals Ltd.,  
(ISO 13485 CERTIFIED)  
Transasia House,  
Chandivali Studio road,  
Andheri (E),  
MUMBAI – 400 072

1 | Page





## Table of Contents

Sr. No.	Contents	Page No.
I	Approval of the OQ procedure	3
II	Instructions	4
III	Scope	5
IV	Operational Qualification	6
V	System Certification	13



## I. Approval of the OQ procedure:

MAEER's Vishwaraj Hospital Laboratory and Transasia are jointly responsible for operational check of the HEMATOLOGY Analyzer, Model: XS800i, serial# 67641 in the clinical lab of MAEER's Vishwaraj Hospital Laboratory as per protocol attached.

### Protocol Performed by: Transasia Representative

Name : Mr. Swapnil Kamarikar  
Title : OPERATIONAL QUALIFICATION  
Company : TRANSASIA BIO-MEDICALS LTD.

Signature: *Kes*

Date: 6/8/21



### Validation Team from \_\_\_\_\_

Name :  
Designation :  
Department :

### Customer Authorization:

Name : *Dr. Anindha Gamed*  
Title : OPERATIONAL QUALIFICATION  
Site :

Signature: *Anindha*

Date:



## II. Instructions

1. The TRANSASIA representative will check each module and enter the specific data as described in the Operational Qualification. Each result will be noted and dated.
2. Employee of MAEER's Vishwaraj Hospital Laboratory will verify each result and sign in the last page. The member/s of the validation team will be responsible for the same.
3. Any deviations from the acceptance criteria detailed in this document will be noted in the **COMMENTS** section of the OQ protocol. All resolution to such problems will also be noted in the **COMMENTS** section, and must be resolved prior to issuance of a SYSTEM CERTIFICATION. This will be an additional cost to the purchasing institution (CUSTOMER). However this additional cost will be waived when this test is conducted at time of initial performance check of new instruments.
4. Any test data, which does not meet the specified acceptance criteria, will be submitted to the appropriate laboratory personnel for solution. All steps taken subsequently will be documented.
5. This document contains proprietary information and is in no way to be copied, photographed or duplicated in any way without expressed written authorization by the Product Manager at Transasia Bio-Medicals Ltd., Transasia House, Mumbai.

Validation Team:

Name *Swepril Kamnaker*

Designation *Service Eng.*

Signature *[Signature]*

Date *6/8/21*







### III. Scope

This Operational Qualification protocol will be performed on the Hematology Analyzer, Model XS-800I, Serial No:67641 located in MAEER's Vishwaraj Hospital Laboratory. This Protocol will define the documentation that will be used to evaluate the instruments installation in accordance with the manufacturer's specifications and intended use. Successful completion of this protocol will verify that the instrument identified is performing in accordance with the intended usage.

Trained, knowledgeable personnel will perform qualification studies.

Any exceptional conditions encountered during the qualification studies will be identified for review. Exceptional conditions will be investigated and the appropriate course of action determined. All documents will be initialed and dated.

Validation Team:

Name *Swepnil Kamaniker*

Designation *Service Eng.*

Signature *[Handwritten Signature]*

Date *6/8/21*





## IV. Operational Qualification

### a. Instrument Identification

Verified Date

1. Model Name XS800i  
2. Serial Number 67641

\_\_\_\_\_  
\_\_\_\_\_

b. Following is a list of tests to be performed and verified:

<u>Test No.</u>	<u>Test Name</u>	<u>Test Purpose</u>	<u>Verified Date</u>
1.	Whole Blood (WB) X- aspiration motor operation	To the WB aspiration motor operation	_____
2	Sheath Motor Test.	To check Operation of Sheath Motor	_____
3.	Diagnostic Test for Auto Sampler and Bar code reader	To Check operation of Auto sampler and Barcode reader.	_____

Validation Team:

Name Swapnil Komenkar

Designation Service Eng.

Signature

Date 6/8/21





## c. Operational Testing

### Test 1

**Test Name** : Whole Blood Aspiration Motor Test  
**Purpose** : To test the Aspiration Motor movement  
**Method** : Please follow the steps described in chapter 2, page 2.9.1 of Sysmex XN – 350 manual

	<u>PARAMETER</u>	<u>PASS</u>	<u>FAIL</u>
Parameter values for verification :	Whold Blood Aspiration Motor Test	✓	

Validation Team:

Name *Swapnil Kameniker*

Designation *Service Eng*

Signature *[Signature]*

Date *6/8/21*







## Test 2

**Test Name** : Sheath Motor Test  
**Purpose** : To test the Sheath Motor Operation Test.  
**Method** : Please follow the steps described in chapter 2, page 2.9.1 of Sysmex XS800i's manual No.

	<u>PARAMETER</u>	<u>PASS</u>	<u>FAIL</u>
Parameter values for verification :	Sheath Motor Motor Test	✓	

Validation Team:

Name *Sudhail Kamaniker*

Designation *Service Eng.*

Signature *[Signature]*

Date *6/8/21*





### Test 3

**Test Name :** Diagnostics Test for Auto Sampler & Barcode Operation  
**Purpose :** To test the Operation of Auto Sampler & Barcode.  
**Method :** Please follow the steps described in chapter 2, page 2.9.1&4 of Sysmex XS800i manual No.

	<u>PARAMETER</u>	<u>PASS</u>	<u>FAIL</u>
Parameter values for verification :	Sampler & Barcode Test	✓	

Validation Team:

Name *Swepnil Kamaniker*

Designation *service Eng -*

Signature *[Handwritten Signature]*

Date *6/8/21*





## d. Operational Procedure

### a. Certificate of Training

#### 1. Technician Training

This certifies that the technicians listed below have received basic user training in the following categories for the system described in this Installation Qualification.

Mr. Vijay Shahu who is certified by Transasia Bio-Medicals Ltd has conducted the training.

Sr.No.	Training Program	Initials	Date
1.	Instrument Setup		
2.	System Operation		
3.	Basic Troubleshooting & Maintenance		

#### 2. Operator Training

The users responsible for the operation of this instrument will be trained in the proper usage of the system. Training will focus on the basic operation and maintenance of the system. The training of the operators will be documented and the training records will be filed as indicated below:

Sr.No.	Operators	Location	Initials	Date
1	Kiran Pandit	VRH	Kiran	
2	Sachin Bankar	VRH	Sachin	
3	Krishna Kendre	VRH	Kandre	

Validation Team:

Name Swapnil Kamanikar

Designation Service Engg

Signature [Signature]

Date 6/8/21





## b. Customer SOP

Title	Number	Revision #	Effective Date	Location	Verified By	Date
Operating Procedure		NA				

Validation Team:

Name *Surendra Komeniker*

Designation *Service Eng*

Signature *[Signature]*

Date *6/8/21*







## COMMENTS:

Validation Team:

Name *Swapanil Kameniker*

Designation *Service Eng.*

Signature *[Handwritten Signature]*

Date *6/8/21*





## V. System Certification

Study data has determined that the system described in this document either meets all criteria outlined in this Operational Qualification Protocol, or exceptional conditions have been identified and documentation included. Exceptional conditions, if any, have been addressed. The system is ready for specified usage.

Report Performed By: Transasia Representative

Name : Mr. Swapnil Kamarikar

Title : OPERATIONAL QUALIFICATION Signature: *[Signature]*

Company: TRANSASIA BIO-MEDICALS LTD. Date : 6/8/21



### Customer Authorizations:

Name : *Dr. Anindha Gamed*

Title : OPERATIONAL QUALIFICATION Signature:

Site : Date :

Name : *Dr. Anindha Gamed*

Title : OPERATIONAL QUALIFICATION Signature: *[Signature]*

Site : Date :



# SYSMEX XS800i AUTOMATED HEMATOLOGY ANALYZER

## PERFORMANCE QUALIFICATION

For

“MAEER’s Vishwaraj Hospital Laboratory”

Marketed by:  
Transasia Bio-Medicals Ltd.,  
(ISO 9002 CERTIFIED)  
Transasia House,  
Chandivali Studio road,  
Andheri (E),  
MUMBAI – 400 072



UNMATCHED SERVICE  
SINCE 1979...

## Table of Contents:

Sr. No.	Contents	Page No.
I	Approval of the PQ procedure	3
II	Instructions	4
III	Scope	5
IV	Performance Qualification	6
V	System Certification	10





## I. Approval of the PQ procedure

MAEER's Vishwaraj Hospital Laboratory and Transasia are jointly responsible for conducting the Performance Check of the Hematology Analyzer, SYSMEX Model : XS800i, Serial No.67641 in the clinical lab of MAEER's Vishwaraj Hospital Laboratory as per the attached protocol.

### Protocol Performed By: Transasia Representative

Name : Mr. Vijay Shahu  
Title : PERFORMANCE QUALIFICATION  
Company : TRANSASIA BIO-MEDICALS LTD.

Signature: *Vijay*  
Date: 6/8/21



### Validation Team from :

Name :  
Designation :  
Department :

Name :  
Designation :  
Department :

### Customer Authorizations:

Name : *Dr. Anindha Gaud*  
Title : PERFORMANCE QUALIFICATION  
Site : Pune

Signature: *Anindha*

Date: 6/8/21

Name : *Dr. Anindha Gaud*  
Title : PERFORMANCE QUALIFICATION  
Site : Pune

Signature: *Anindha*

Date: 6/8/21



## II. Instructions

1. An authorized TRANSASIA representative will check for the performance of the instrument and enter the specific data as outlined in the Performance Qualification. Each result will be noted and dated.
2. Performance checks on a regular basis described in the Further Performance Checks (vide-infra) will be responsibility of the customer's personnel.
3. Employee of MAEER's Vishwaraj Hospital Laboratory will verify each result and sign in the last page. The members of the validation team will carry this out.
4. ALL deviations from the acceptance criteria detailed in this document will be noted in the COMMENTS section at the end of each PQ protocol. All resolution to such problems will also be noted in the COMMENTS section, and must be resolved prior to issuance of a SYSTEM CERTIFICATION. These will be an additional cost to the purchasing institution (CUSTOMER). However this additional cost will be waived when this test is conducted at time of initial performance check of new instruments.
5. Any test data that does not meet the specified acceptance criteria will be submitted to the appropriate laboratory personnel for solution. All steps taken subsequently will be documented.
6. This document contains proprietary information and is in no way to be copied, photographed or duplicated in any way without expressed written authorization by the Production Manager at Transasia Bio-Medicals Ltd., Transasia House, Mumbai.

Validation Team:

Name Vijay Saher  
Designation Vij App Spl.  
Signature  
Date 6/8/21



4 | Page



### III. Scope

This Performance Qualification protocol will be performed on the Hematology Analyzer, Model XS-800I, Serial No 67641 located in MAEER's Vishwaraj Hospital Laboratory. This Protocol will define the documentation that will be used to evaluate the instruments installation in accordance with the manufacture's specifications and intended use. Successful completion of this protocol will verify that the instrument identified is performing in accordance with the intended usage.

Trained, knowledgeable personnel will perform qualification studies.

Any exceptional conditions encountered during the qualification studies will be identified for review. Exceptional conditions will be investigated and the appropriate course of action determined. All documents will be initialed and dated.

Validation Team:

Name *Vijay Sahu*  
Designation *App. Spl.*  
Signature *Vijay*  
Date *6/8/21*



5 | Page





## IV. Performance Qualification

### a. Instrument Identification

Verified Date

1. Model Name SYSMEX XS800i \_\_\_\_\_  
2. Serial Number 67641 \_\_\_\_\_

b. Following is a list of tests to be performed and verified:

<u>Test No.</u>	<u>Test Name</u>	<u>Test Purpose</u>	<u>Verified Date</u>
02	Sample Processing	Ability to process samples	_____
03	Further Performance Checks	Regular Maintenance	NA

Validation Team:

Name *Vijay Sahu*  
Designation *App. Spl.*  
Signature *V.S.*  
Date *6/8/21*



6 | Page





## c. Performance Testing

### Test 1

**Test Name:** Sample Processing  
**Purpose:** Ability to Process Samples  
**Method:**

#### 1. Run the control samples five times consecutively

**Acceptance Criteria:** Each of the results obtained above should be within the range as specified in the control chart.

#### Parameters Values for Verification:

#### RBC Count:

Test	Control Values	Results Obtained	Pass	Fail
1.	4.2	4.1	✓	
2.	4.2	4.2	✓	
3.	4.2	4.2	✓	
4.	4.2	4.2	✓	
5.	4.2	4.2	✓	

#### Validation Team:

Name *Vijay sahu*  
Designation *App Spl.*  
Signature *Vij*  
Date *6/8/21*





### WBC Count:

Test	Control Values	Results Obtained	Pass	Fail
1.	6.1	6.2	✓	
2.	6.1	5.9	✓	
3.	6.1	6.0	✓	
4.	6.1	6.1	✓	
5.	6.1	6.1	✓	

### Hemoglobin:

Test	Control Values	Results Obtained	Pass	Fail
1.	12.1	11.9	✓	
2.	12.1	12.0	✓	
3.	12.1	12.0	✓	
4.	12.1	12.0	✓	
5.	12.1	12.0	✓	

### Platelet Count:

Test	Control Values	Results Obtained	Pass	Fail
1.	350	340	✓	
2.	350	345	✓	
3.	350	343	✓	
4.	350	333	✓	
5.	350	352	✓	

### Validation Team:

Name vijay sahu

Designation Appl. spl.

Signature

Date 6/8/21





## Test 2

### Test Name:

1. Tests for checking the performance of the instruments during analysis
2. Tests for checking long term performance of the instrument

### Purpose:

The purpose of the above checks is to ensure the reliability of the results being obtained.

### Method:

#### 1. During Sample analysis:

To run control samples each time the instrument is used for sample analysis and verification of the results of the controls to be within the reference range to be established by performance of the precision experiments.

#### 2. Long term Performance

This is to be checked by Levy Jennings plots to be updated once in six months

### Validation Team:

Name *Vijay Sahu*  
Designation *Appl. Spl.*  
Signature *Vij*  
Date *6/8/21*





## V. System Certification

Study data has determined that the system described in this document either meets all criteria outlined in this Performance Qualification Protocol, or exceptional conditions have been identified and documentation included. Exceptional conditions, if any, have been addressed. The system is ready for specified usage.

Report Performed By: Transasia Representative

Name : Mr. Vijay Shahu

Title : PERFORMANCE QUALIFICATION Signature: *Vij*

Company: TRANSASIA BIO-MEDICALS LTD. Date : 6/8/21



Customer Authorizations:

Name : *Dr. Avinodhar Ganesh*

Title : PERFORMANCE QUALIFICATION Signature:

Site : Date : 6/8/21

Name : *Dr. Avinodhar Ganesh*

Title : PERFORMANCE QUALIFICATION Signature:

Site : Date : 6/8/21



Item

Sample No.

P00

Rack:

Tube:


2021/08/06 20:57:35

Ward:


Dr.:

WBC	5.75	[10 <sup>3</sup> /uL]		
RBC	5.17	[10 <sup>6</sup> /uL]		
HGB	12.9	[g/dL]		
HCT	39.7	[%]		
MCV	76.8	[fL]		
MCH	25.0	[pg]		
MCHC	32.5	[g/dL]		
PLT	290	[10 <sup>3</sup> /uL]		
RDW-SD	34.4	[fL]		
RDW-CV	12.6	[%]		
PDW	12.5	[fL]		
MPV	11.2	[fL]		
P-LCR	34.9	[%]		
PCT	0.33	[%]		
NEUT	3.60	[10 <sup>3</sup> /uL]	62.6	[%]
LYMPH	1.84	[10 <sup>3</sup> /uL]	32.0	[%]
MONO	0.21	[10 <sup>3</sup> /uL]	3.7	[%]
EO	0.06	[10 <sup>3</sup> /uL]	1.0	[%]
BASO	0.04	[10 <sup>3</sup> /uL]	0.7	[%]


WBC



RBC



PLT




DIFF

Sample No. P01 Rack: Tube: 2021/08/06 20:59:03  
 Ward: Dr.:

WBC	5.55	[10 <sup>3</sup> /uL]		
RBC	5.20	[10 <sup>6</sup> /uL]		
HGB	13.0	[g/dL]		
HCT	39.8	[%]		
MCV	76.5	[fL]		
MCH	25.0	[pg]		
MCHC	32.7	[g/dL]		
PLT	292	[10 <sup>3</sup> /uL]		
RDW-SD	34.2	[fL]		
RDW-CV	12.4	[%]		
PDW	13.3	[fL]		
MPV	11.3	[fL]		
P-LCR	35.2	[%]		
PCT	0.33	[%]		
NEUT	3.44	[10 <sup>3</sup> /uL]	62.0	[%]
LYMPH	1.81	[10 <sup>3</sup> /uL]	32.6	[%]
MONO	0.21	[10 <sup>3</sup> /uL]	3.8	[%]
EO	0.06	[10 <sup>3</sup> /uL]	1.1	[%]
BASO	0.03	[10 <sup>3</sup> /uL]	0.5	[%]


WBC



RBC



PLT




DIFF

Sample No. P02 Rack: Tube: 2021/08/06 21:00:10

Ward: Dr.:

WBC	5.60	[10 <sup>3</sup> /uL]		
RBC	5.29	[10 <sup>6</sup> /uL]		
HGB	13.0	[g/dL]		
HCT	40.4	[%]		
MCV	76.4	[fL]		
MCH	24.6	[pg]		
MCHC	32.2	[g/dL]		
PLT	285	[10 <sup>3</sup> /uL]		
RDW-SD	34.2	[fL]		
RDW-CV	12.4	[%]		
PDW	13.3	[fL]		
MPV	11.3	[fL]		
P-LCR	35.3	[%]		
PCT	0.32	[%]		
NEUT	3.46	[10 <sup>3</sup> /uL]	61.8	[%]
LYMPH	1.79	[10 <sup>3</sup> /uL]	32.0	[%]
MONO	0.22	[10 <sup>3</sup> /uL]	3.9	[%]
EO	0.08	[10 <sup>3</sup> /uL]	1.4	[%]
BASO	0.05	[10 <sup>3</sup> /uL]	0.9	[%]


WBC



RBC



PLT



DIFF

Sample No. P03 Rack: Tube: 2021/08/06 21:01:10  
 Ward: Dr.:

WBC	5.66	[10 <sup>3</sup> /uL]		
RBC	5.19	[10 <sup>6</sup> /uL]		
HGB	12.8	[g/dL]		
HCT	39.7	[%]		
MCV	76.5	[fL]		
MCH	24.7	[pg]		
MCHC	32.2	[g/dL]		
PLT	293	[10 <sup>3</sup> /uL]		
RDW-SD	34.1	[fL]		
RDW-CV	12.4	[%]		
PDW	13.8	[fL]		
MPV	11.4	[fL]		
P-LCR	36.1	[%]		
PCT	0.33	[%]		
NEUT	3.54	[10 <sup>3</sup> /uL]	62.6	[%]
LYMPH	1.81	[10 <sup>3</sup> /uL]	32.0	[%]
MONO	0.21	[10 <sup>3</sup> /uL]	3.7	[%]
EO	0.07	[10 <sup>3</sup> /uL]	1.2	[%]
BASO	0.03	[10 <sup>3</sup> /uL]	0.5	[%]

WBC  
 RBC  
 PLT  
 DIFF




Sample No. P04 Rack: Tube: 2021/08/06 21:02:18


Ward: Dr.:

WBC	5.87	[10 <sup>3</sup> /uL]		
RBC	5.22	[10 <sup>6</sup> /uL]		
HGB	12.9	[g/dL]		
HCT	39.8	[%]		
MCV	76.2	[fL]		
MCH	24.7	[pg]		
MCHC	32.4	[g/dL]		
PLT	285	[10 <sup>3</sup> /uL]		
RDW-SD	33.8	[fL]		
RDW-CV	12.3	[%]		
PDW	13.3	[fL]		
MPV	11.2	[fL]		
P-LCR	34.2	[%]		
PCT	0.32	[%]		
NEUT	3.60	[10 <sup>3</sup> /uL]	61.3	[%]
LYMPH	1.92	[10 <sup>3</sup> /uL]	32.7	[%]
MONO	0.22	[10 <sup>3</sup> /uL]	3.7	[%]
EO	0.08	[10 <sup>3</sup> /uL]	1.4	[%]
BASO	0.05	[10 <sup>3</sup> /uL]	0.9	[%]


WBC



RBC



PLT



DIFF

Sample No.

P05

Rack:

Tube:


2021/08/06 21:03:26

Ward:

Dr.:

WBC	5.56	[10 <sup>3</sup> /uL]		
RBC	5.19	[10 <sup>6</sup> /uL]		
HGB	12.8	[g/dL]		
HCT	39.6	[%]		
MCV	76.3	[fL]		
MCH	24.7	[pg]		
MCHC	32.3	[g/dL]		
PLT	287	[10 <sup>3</sup> /uL]		
RDW-SD	34.1	[fL]		
RDW-CV	12.4	[%]		
PDW	13.1	[fL]		
MPV	11.2	[fL]		
P-LCR	35.6	[%]		
PCT	0.32	[%]		
NEUT	3.38	[10 <sup>3</sup> /uL]	60.8	[%]
LYMPH	1.85	[10 <sup>3</sup> /uL]	33.3	[%]
MONO	0.22	[10 <sup>3</sup> /uL]	4.0	[%]
EO	0.08	[10 <sup>3</sup> /uL]	1.4	[%]
BASO	0.03	[10 <sup>3</sup> /uL]	0.5	[%]

WBC



RBC



PLT



DIFF




Sample No. P06 Rack: Tube: 2021/08/06 21:04:29  
 Ward: Dr.:

WBC	5.72	[10 <sup>3</sup> /uL]		
RBC	5.24	[10 <sup>6</sup> /uL]		
HGB	12.8	[g/dL]		
HCT	40.0	[%]		
MCV	76.3	[fL]		
MCH	24.4	[pg]		
MCHC	32.0	[g/dL]		
PLT	281	[10 <sup>3</sup> /uL]		
RDW-SD	34.4	[fL]		
RDW-CV	12.5	[%]		
PDW	12.9	[fL]		
MPV	11.2	[fL]		
P-LCR	35.4	[%]		
PCT	0.31	[%]		
NEUT	3.55	[10 <sup>3</sup> /uL]	62.2	[%]
LYMPH	1.82	[10 <sup>3</sup> /uL]	31.8	[%]
MONO	0.22	[10 <sup>3</sup> /uL]	3.8	[%]
EO	0.10	[10 <sup>3</sup> /uL]	1.7	[%]
BASO	0.02	[10 <sup>3</sup> /uL]	0.5	[%]

WBC  
 RBC  
 PLT  
 DIFF

Sample No. P07 Rack: Tube: 2021/08/06 21:05:40  
 Ward: Dr.:

WBC	5.58	[10 <sup>3</sup> /uL]		
RBC	5.26	[10 <sup>6</sup> /uL]		
HGB	12.8	[g/dL]		
HCT	40.0	[%]		
MCV	76.0	[fL]		
MCH	24.3	[pg]		
MCHC	32.0	[g/dL]		
PLT	295	[10 <sup>3</sup> /uL]		
RDW-SD	34.2	[fL]		
RDW-CV	12.4	[%]		
PDW	13.9	[fL]		
MPV	11.4	[fL]		
P-LCR	36.1	[%]		
PCT	0.33	[%]		
NEUT	3.44	[10 <sup>3</sup> /uL]	61.7	[%]
LYMPH	1.78	[10 <sup>3</sup> /uL]	31.9	[%]
MONO	0.24	[10 <sup>3</sup> /uL]	4.3	[%]
EO	0.08	[10 <sup>3</sup> /uL]	1.4	[%]
BASO	0.04	[10 <sup>3</sup> /uL]	0.7	[%]



DIFF



Sample No.

P08

Rack:

Tube:


2021/08/06 21:06:44

Ward:


Dr.:

WBC	5.62	[10 <sup>3</sup> /uL]		
RBC	5.22	[10 <sup>6</sup> /uL]		
HGB	12.8	[g/dL]		
HCT	39.7	[%]		
MCV	76.1	[fL]		
MCH	24.5	[pg]		
MCHC	32.2	[g/dL]		
PLT	291	[10 <sup>3</sup> /uL]		
RDW-SD	34.0	[fL]		
RDW-CV	12.4	[%]		
PDW	13.4	[fL]		
MPV	11.3	[fL]		
P-LCR	35.6	[%]		
PCT	0.33	[%]		
NEUT	3.48	[10 <sup>3</sup> /uL]	62.0	[%]
LYMPH	1.75	[10 <sup>3</sup> /uL]	31.1	[%]
MONO	0.25	[10 <sup>3</sup> /uL]	4.4	[%]
EO	0.10	[10 <sup>3</sup> /uL]	1.8	[%]
PLAS	0.04	[10 <sup>3</sup> /uL]	0.7	[%]


WBC



RBC



PLT



DIFF

Sample No. P09 Rack: Tube: 2021/08/06 21:07:57  
 Ward: Dr.:

WBC	5.65	[10 <sup>3</sup> /uL]		
RBC	5.30	[10 <sup>6</sup> /uL]		
HGB	13.0	[g/dL]		
HCT	40.2	[%]		
MCV	75.8	[fL]		
MCH	24.5	[pg]		
MCHC	32.3	[g/dL]		
PLT	290	[10 <sup>3</sup> /uL]		
RDW-SD	34.0	[fL]		
RDW-CV	12.5	[%]		
PDW	13.3	[fL]		
MPV	11.2	[fL]		
P-LCR	34.9	[%]		
PCT	0.33	[%]		
NEUT	3.44	[10 <sup>3</sup> /uL]	60.9	[%]
LYMPH	1.84	[10 <sup>3</sup> /uL]	32.6	[%]
MONO	0.24	[10 <sup>3</sup> /uL]	4.2	[%]
EO	0.09	[10 <sup>3</sup> /uL]	1.6	[%]
BASO	0.04	[10 <sup>3</sup> /uL]	0.7	[%]

WBC  
 RBC  
 PLT  
 DIFF

Sample No.

P10

Rack:

Tube:


2021/08/06 21:08:59

Ward:

Dr.:

WBC	5.67	[10 <sup>3</sup> /uL]		
RBC	5.29	[10 <sup>6</sup> /uL]		
HGB	13.0	[g/dL]		
HCT	40.1	[%]		
MCV	75.8	[fL]		
MCH	24.6	[pg]		
MCHC	32.4	[g/dL]		
PLT	296	[10 <sup>3</sup> /uL]		
RDW-SD	34.3	[fL]		
RDW-CV	12.5	[%]		
PDW	13.0	[fL]		
MPV	11.4	[fL]		
P-LCR	35.5	[%]		
PCT	0.34	[%]		
NEUT	3.46	[10 <sup>3</sup> /uL]	61.0	[%]
LYMPH	1.81	[10 <sup>3</sup> /uL]	31.9	[%]
MONO	0.23	[10 <sup>3</sup> /uL]	4.1	[%]
EO	0.12	[10 <sup>3</sup> /uL]	2.1	[%]
BASO	0.05	[10 <sup>3</sup> /uL]	0.9	[%]

WBC



RBC



PLT



DIFF