## Sysmex India Pvt Ltd - West

1002, Damji Shamji Business Galleria, 10th Floor, LBS Marg, Kanjurmarg (West), India, 400078 Tel. (+91) 22 6112 6666 Fax. (+91) 22 2577 6790

Email: sysmex@sysmex.co.in
Web: https://www.sysmex.co.in
GST Reg. No: 27AADCS1551J1ZC Co. Reg. No:



CSR No: WO-100080313 **Revision No:** 

## **Customer Service Report**

	(	Customer Serv	rice Report					
Customer Name	Institute of Post Gr	raduate Teaching	and Research in Ayurv	eda				
Contact Person	Agnihotram Gopina	Agnihotram Gopinath						
Installation Date				Time				
Equipment ID	10109604	10109604			Updesh Dwivedi			
Model	BX-3010	BX-3010		ime	08/09/2019 11:34:00 AM			
Serial No.	7237-0555	7237-0555		e	08/09/2019 06:07:00 PM			
Total Cycle Count								
Defect Analysis			<u>'</u>		1			
Issue Code		C	Cause Code					
Issue			Cause					
Astion Tales								
Action Taken  1. Unpack instrument and remove	fiving motorials							
2. Check all the accessories are s	upplied							
3. Check instrument for leakage o	r damage							
4. Place instrument on flat table at	nd perform Level Adjustm	nent of the analyz	zer					
5. Perform level adjustment of the	analyzer							
6. Connect analyzer to peripherals	(PC, printer, LIS etc)							
7. Connect tubings, external tanks	and tank sensors							
8. Install the silencer in High Conc	and Low Conc. waste to	ubes						
9. Wash the containers (Wash So	lution 1 and 2) and fill with	h respective was	h solutions					
10. Wash the System Water conta	niner and fill with distilled	water						
11. Install BX-3010 Analyzer softw	rare on PC							
13. Prepare sample and reagent for	or Dye Test							
14. Check instrument operation af	ter Power ON							
15. Execute Prime sequence for 1	0 cycles followed by Was	sh Cuvette						
16. Perform Automatic Gain and v	erify results are within ac	cceptable range						
17. Perform Water Blank and verif	y that Absorbance for ea	ach cuvette is with	nin acceptable range					
18. Register reagents and sample	s for Dye test							
19. Enter chemistry parameter and	d calibration settings for [	Dye Test						
19. Load samples and reagents fo	or Dye Test							
20. Execute the Run and verify that	at Dye test results are wit	thin acceptable ra	ange					

21. Verify accuracy and precision (within run) using controls and pooled serum

Items Consumed								
Product ID	VPN	Product Description	QTY UOM					
Name:- Updesh Dwivedi		<u> </u>	Acknowledged By : Agnihotram Gopinath					
Sign:		poleth	Sign:					
<b>Date:</b> 08/09/2019 06:14:55 PM		M	Date: 08/09/2019 06:14:55 PM					

Remarks