a PerkinElmer company

CIN: U85195GA1988PTC005717

Thrissur Office:

Thattil Nadakalan Complex, Kuriachira P.O.,

Thrissur - 680 006, INDIA.

Tel.: (0487) 2250793, 6544048. Mobile No.: 9645084048.

E-mail: tulipkerala@gmail.com

Coralyzer Series Fully Automatic Biochemistry Analyzer, Inspection/calibration Report

Test QC refer to: <<Coralyzer Series Fully Automatic Analyser>> Registered Standard Model: Coralyzer Junior, S/N: AES2KOKH12056

SRI SAMRAJ LAB'S , CUDDALORE

	used for calibration /			I Date:- 30.11.21	Due Date:-29.1	7		
	Inspection Item	· · · · · · · · · · · · · · · · · · ·			Result	Conclusion		
		a) Analytical method may analyse 40 pos at the same time.			OK			
		b) Programming item			OK			
		c) Automatic priority inspect of emergency sample.			OK	-		
		d) Calibration linearity/non-linearity multi point calibration.			OK			
		e) QC at may operate most of QC at same item. g)Storage			OK			
					OK			
		h) Printer multiple report formats available, user can self edit report formats.			OK	PASS		
		i) STAT			OK	1733		
		j) Sample rerun/rech	mple rerun/recheck. OK					
		k) Sample /reagent p	probe liquid check.		ОК			
		I) Mixing function.		OK .				
		m)Results storage, query, quality control chart, the whole process			OK			
		n) cuvette condetion			ОК	1		
		o) Crash gard functi	on.		ОК			
2	Dark Abs	Absorbance >= 3.400) A		5.2	ОК		
3	Linearity range of	The maximum absor	bance should be not less than 3.7	00A within the	340nm: 5.1	OK		
	absorbance.	relative range of ± 59	elative range of ± 5%			OK		
4	Accuracy of absorbance	Less than ± 0.040A at 1.00A			0.012	OK		
		Less than ± 0.060A at 1500A			0.027	OK		
		Less than ± 0.020A within 0.500A			0.048	OK		
5	Accuracy	Absorbance change should be less than 0.005A/min		340nm: 0.003A	OK			
	· ·	,			505nm: 0.002A	OK		
6	Reproducibility of absorbance	cofficient of variation cv ≤ 3%			2.73%	ОК		
7	Temperature of	Accuracy less than than ± 0.1° C			37°C	OK		
,	reaction cuvette 37	Fluctuate of temperature less than ±0.1°C			0.005°C	OK		
8	Teperature of reagent storage	Under normal work condition, reagent storage should be 2-8°C			6.0 °C	ОК		
9	Carryover	Sample carryover les	ss than 0.05%		PASS			
		Items	Concentration	CV				
10	Accuracy of parameter	UREA	76.3-109.9	≦ 3%	1.30%	OK		
10		GLUC	81.0-116.6	≦ 3%	1.20%	OK		
		SGOT	55-87	≦ 3%	2.70%	OK		
		a) Whole appearance of the surface should not have obvious mechanical abrasion,				,		
11		corrosion and surface coating spalling. b) All connections are stable and fixed well. c) Moving parts running smooth and CPU working well.			DAGG			
	Apperance				PASS	-		
						1		
	Protect ground	The impedance between ground outlet from power supply net and ground						
12	connect	protected touchable metal should not more than 0.1 Ω			PASS			
13	Power grid, power supply, circuit	Basic insulation:50Hz the sine 230V alternating voltage test. should not apper brakdown or occur electric arcing.			PASS			

Inspection/Calibration date: 30.11.21

Check: S Nistianth

re-check: S Nishanth

Corporate Office: Gitanjali, Tulip Block, Dr. Antonio No Bogo Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA. Tel.: (0832) 2458546-51, Fax: (0832) 245854 Pimail sales@tulipgroup.com, Website: www.tulipgroup.com