UNMATCHED SERVICE  
SINCE 1979...

**TRANSASIA BIOMEDICALS LIMITED**  
**PERFORMANCE QUALIFICATION**

**TRANSASIA®**  
Bio-Medicals Ltd.

Customer Name: UJALA CYGNUS RAINBOW HOSPITAL, AGRA  
 Model: EC 90 Serial Number: 4457-2-202209

**ANNEXURE-A**

**TABLE OF CONTENTS**

| Sr.No. | Title               | Page. No. |
|--------|---------------------|-----------|
| 1.     | Approval            | 1         |
| 2.     | Objectives          | 2         |
| 3.     | Pre-Requisites      | 2         |
| 4.     | Test Plan           | 2         |
| 5.     | Acceptance Criteria | 2         |
| 6.     | Summary             | 2         |

**1. APPROVAL:**

**1.1 Prepared By:**

| Name              | Designation            | Signature         | Date |
|-------------------|------------------------|-------------------|------|
| VED PRAKASH SINGH | APPLICATION SPECIALIST | <i>Vedprakash</i> |      |

**1.2 Checked By:**

| Name | Designation | Signature | Date |
|------|-------------|-----------|------|
|      |             |           |      |
|      |             |           |      |
|      |             |           |      |

**1.3 Approved By:**

| Name | Designation | Signature | Date |
|------|-------------|-----------|------|
|      |             |           |      |

UNMATCHED SERVICE  
 SINCE 1979...

## TRANSASIA BIOMEDICALS LIMITED

### PERFORMANCE QUALIFICATION

|  |                     |                |               |
|--|---------------------|----------------|---------------|
| Customer Name: UJALA CYGNUS RAINBOW HOSPITAL, AGRA | Model: <u>EC 90</u> | Serial Number: | 4457-2-202209 |
|--|---------------------|----------------|---------------|

**2. OBJECTIVE:**

The objective of this protocol is to establish documented evidence for the Performance Qualification of EC 90 and to ensure that the results obtained are within the pre-determined Acceptance Criteria.

**3. PRE-REQUISITES:**

Following Pre-requisites are required before the execution of Performance Qualification.

- Completion of Installation Qualification (IQ) prior to PQ.
- Completion of Operational Qualification (OQ) prior to PQ.

**4. TEST PLAN:**

Perform precision study for both Erba Norm & Erba Path controls and check the recovery (for both level QC) & CV% (for both level QC) during the Performance Qualification of EC 90.

- Follow the EC 90 Quick Start Guide for Test procedure.

**5. ACCEPTANCE CRITERIA:**


Acceptance criteria are based on the precision studies conducted by manufacturer for EC 90.

**6. SUMMARY:**

| Level                           | Norm QC | Path QC |
|---------------------------------|---------|---------|
| Lot #                           |         |         |
| Mean                            |         |         |
| SD                              |         |         |
| CV%                             |         |         |
| Range Limit Low                 |         |         |
| Range Limit High                |         |         |
| Acceptable Criteria (Max.value) |         |         |
| Status                          | PASS    | PASS    |

Traceability: Instrument Raw Data printouts.

Note: This report is effective from the date of approval.

UNMATCHED SERVICE  
SINCE 1979...

## TRANSASIA BIOMEDICALS LIMITED

### PERFORMANCE QUALIFICATION

**TRANSASIA®**  
Bio-Medicals Ltd.

Customer Name: UJALA CYGNUS RAINBOW  
HOSPITAL, AGRA

Model: EC 90

Serial  
Number:

4457-2-  
202209

#### Annexure-B

Instrument Raw Data to be attached as printout attachments in Annexure-B.

UNMATCHED SERVICE  
SINCE 1979...

## TRANSASIA BIOMEDICALS LIMITED

### OPERATIONAL QUALIFICATION

Customer Name: UJALA CYGNUS RAINBOW  
HOSPITAL, AGRA

Model: EC 90

Serial  
Number:

4457-2-  
202209

TRANSASIA®  
Bio-Medicals Ltd.

#### ANNEXURE-A

#### TABLE OF CONTENTS

| Sr.No. | Title               | Page. No. |
|--------|---------------------|-----------|
| 1.     | Approval            | 1         |
| 2.     | Objectives          | 2         |
| 3.     | Pre-Requisites      | 2         |
| 4.     | Test Plan           | 2         |
| 5.     | Acceptance Criteria | 2         |

#### 1. APPROVAL:

##### 1.1 Prepared By:

| Name             | Designation      | Signature | Date |
|------------------|------------------|-----------|------|
| PUSHPENDRA SINGH | SERVICE ENGINEER |           |      |

##### 1.2 Checked By:

| Name | Designation | Signature | Date |
|------|-------------|-----------|------|
|      |             |           |      |
|      |             |           |      |
|      |             |           |      |

##### 1.3 Approved By:

| Name | Designation | Signature | Date |
|------|-------------|-----------|------|
|      |             |           |      |

UNMATCHED SERVICE  
SINCE 1979...

## TRANSASIA BIOMEDICALS LIMITED

### OPERATIONAL QUALIFICATION

|   |                     |                   |                   |  |
|---|---------------------|-------------------|-------------------|--|
| Customer Name: UJALA CYGNUS RAINBOW<br>HOSPITAL, AGRA | Model: <u>EC 90</u> | Serial<br>Number: | 4457-2-<br>202209 | <b>TRANSASIA®</b><br>Bio-Medicals Ltd. |
|   |                     |                   |                   |  |

#### 2. OBJECTIVE:

The objective of this protocol is to establish documented evidence for the Operational Qualification of EC 90

#### 3. PRE-REQUISITES:

Following Pre-requisites are required before the execution of Performance Qualification.

- Completion of Installation Qualification (IQ) prior to PQ.


#### 4. TEST PLAN:

- Setting up the System settings by selecting date, time and language.
- Checked Tubes inserted, found OK.
- Checked EC Cartridge inserted with Biosensor found OK.

#### 5. ACCEPTANCE CRITERIA:

Prime Up and 2P Calibration are resulted with NO ERROR.

Note: This report is effective from the date of approval.

UNMATCHED SERVICE  
SINCE 1979...

## TRANSASIA BIOMEDICALS LIMITED

### INSTALLATION QUALIFICATION

Customer Name: UJALA CYGNUS RAINBOW  
HOSPITAL, AGRA

Model: EC 90

Serial  
Number:

4457-2-  
202209

**TRANSASIA®**  
Bio-Medicals Ltd.

#### ANNEXURE-A

#### TABLE OF CONTENTS

| Sr.No. | Title                   | Page. No. |
|--------|-------------------------|-----------|
| 1.     | Approval                | 1         |
| 2.     | Objectives              | 2         |
| 3.     | Pre-Requisites          | 2         |
| 4.     | Unpacking and placement | 2         |
| 5      | Installation            | 2         |

#### 1. APPROVAL:

##### 1.1 Prepared By:

| Name             | Designation      | Signature | Date |
|------------------|------------------|-----------|------|
| PUSHPENDRA SINGH | SERVICE ENGINEER |           |      |

##### 1.2 Checked By:

| Name | Designation | Signature | Date |
|------|-------------|-----------|------|
|      |             |           |      |
|      |             |           |      |
|      |             |           |      |

##### 1.3 Approved By:

| Name | Designation | Signature | Date |
|------|-------------|-----------|------|
|      |             |           |      |

UNMATCHED SERVICE  
SINCE 1979...

## TRANSASIA BIOMEDICALS LIMITED

### INSTALLATION QUALIFICATION

TRANSASIA®  
Bio-Medicals Ltd.

Customer Name: UJALA CYGNUS RAINBOW  
HOSPITAL, AGRA

Model: EC 90

Serial  
Number:

4457-2-  
202209

#### 2 OBJECTIVE:

The objective of this protocol is to establish documented evidence for the Installation Qualification of EC-90.

#### 3 PRE-REQUISITES:

Initial Inspection of the unit carried out and the details are as follows:

- Analyzer Condition Report:  
Found the Analyzer along with accessories to have been delivered in satisfactory condition and no external physical damaged observed.
- External Requirements for Installation:
- Input voltage of 100-240VAC 50/60Hz 60W
- Recommended Operational Conditions 15-35 °C / 15-90% Humidity
- In order to ensure the safety of the analyzer and the operator, following conditions must be guaranteed.
- The power supply system (Protection Class B & Category II) must be compatible with the voltage and current specifications indicated on the serial number plate attached to the rear of the analyzer; and it is recommended to check the efficiency of the electrical system from time to time.
- The material for the operator's safety (gloves, containers for disposing of used expendable items, surface disinfectants) must always be available.


#### 4 UNPACKING AND PLACEMENT :

- The unpacking & positioning of the analyzer must comply with the instructions indicated in user manual.

#### 5 INSTALLATION :

- The installation of the analyzer must be in accordance to the instructions indicated in the EC 90\_Quick Start Guide (Available in Accessory box along with the analyzer).
- After verification of the installation of the analyzer, flip toggle the power switch to ON (I).
- Once switched on, the analyzer guides the operator to perform 11 relevant steps needed to install the analyzer.
- On confirming the final step, the analyzer performs a Start Up Completion (Prime Up followed by 2P Calibration)
- The analyzer returns to Home Screen after successful Start Up completion process indicating readiness for sample measurement.

Note: This report is effective from the date of approval.

Page 2 of 2

## Certificate of Analysis

Product name: EC 90

Serial No.: 4457-2-202209


Catalog No.: INS00079

Manufacturing Date: 09/2022

### Results of testing

| Testing             | Limit/criterion        | Test results | Remarks |
|---------------------|------------------------|--------------|---------|
| Instrument QC       | pass                   | pass         | ok      |
| QC measurements     | more than 91,7 % valid | 100 %        | ok      |
| Conclusive SW check | pass                   | pass         | ok      |
| Completeness        | pass                   | pass         | ok      |

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## System QC protocol

### EC 90 Analyzer

|                          |               |                    |           |
|--------------------------|---------------|--------------------|-----------|
| Serial number of product | 4457-2-202209 | Manufacturing Date | 09/2022   |
| SW version               | 1.11.3        | Test performed by  | Burianová |
| EC cartridge lot         | 2111139       | Test performed at  | 29.9.2022 |

| No.                            | Performance content        | Limit/Criterion   | Sample    | Notice   | Pass/Fail |
|--------------------------------|----------------------------|---|-----------|----------|-----------|
| 1.                             | Final QC instrument        | valid protocol  | -         | -        | ✓         |
| <b>Basic QC measurements</b>   |                            |   |           |          |           |
| 2.                             | QC measurement performance | less than 8,3% of all values are out of target range      | Erba Norm | 100 %    | ✓         |
|                                |                            |   | Lot:      | 2105150B |           |
| Erba Path                      |                            |   | 100 %     | ✓        |           |
| Lot:                           |                            |   | 2203101   |          |           |
| 3.                             |                            |   |           |          |           |
| <b>Conclusive system check</b> |                            |   |           |          |           |
| 4.                             | SW check                   | all tasks could be done per SW, no SW-freeze              | -         | -        | ✓         |
| 5.                             | Mechanical check           | no mechanical damage detectable                           | -         | -        | ✓         |
| 6.                             | Put into operation         | no additional maintenance for put into operation required | -         | -        | ✓         |
| 7.                             | Closing activities         | Cleaning, visual tube check, shutdown                     | -         | -        | ✓         |

| Disposition | Name      | Date      | Signature        |
|-------------|-----------|-----------|------------------|
| Pass/Fail   | Burianová | 29.9.2022 | <i>Burianová</i> |

|                |               |
|----------------|---------------|
| Serial number: | 4457-2-202209 |
| Date:          | 29.09.2022    |
| Technician     | Burianová     |

|  |                     |
|--|---------------------|
| Materials                              | Lot Number          |
| EC cartridge plus iCa<br>Elite H Clean | 2111139<br>20210906 |

*buianova*

|           |   | QC      |       |       |        | Target values |      |      |      |
|-----------|---|---------|-------|-------|--------|---------------|------|------|------|
|           |   | Na      | K     | iCa   | Cl     | Na            | K    | iCa  | Cl   |
| Erba Norm | 1 | 121,260 | 4,028 | 0,574 | 87,990 | 121           | 3,97 | 0,53 | 88,7 |
|           | 2 | 122,310 | 4,201 | 0,586 | 88,650 | 133           | 4,36 | 0,62 | 97,7 |
|           | 3 | 121,870 | 4,080 | 0,579 | 89,670 | 156           | 6,64 | 1,11 | 130  |
|           | 4 |         |       |       |        | 192           | 7,30 | 1,29 | 142  |
|           | 5 |         |       |       |        |               |      |      |      |
|           | 6 |         |       |       |        |               |      |      |      |
|           | 7 |         |       |       |        |               |      |      |      |
|           | 8 |         |       |       |        |               |      |      |      |
|           | 9 |         |       |       |        |               |      |      |      |
| Mean      |   | 121,870 | 4,080 | 0,579 | 88,650 |               |      |      |      |

|                 |  |
|-----------------|--|
| Requirement     | 8,3% of all QC measurement can be out of the range |
| Material        | Erba Norm  |
| Values OutR     | 0  |
| Total measur. % | 12   |
| Result          | OK   |

|           |   | Na      | K     | iCa   | Cl      | Na  | K    | iCa  | Cl  |
|-----------|---|---------|-------|-------|---------|-----|------|------|-----|
| Erba Path | 1 | 170,750 | 6,565 | 1,070 | 127,750 | 174 | 6,64 | 1,11 | 130 |
|           | 2 | 170,980 | 6,719 | 1,085 | 127,060 | 192 | 7,30 | 1,29 | 142 |
|           | 3 | 171,240 | 6,522 | 1,078 | 127,640 |     |      |      |     |
|           | 4 |         |       |       |         |     |      |      |     |
|           | 5 |         |       |       |         |     |      |      |     |
|           | 6 |         |       |       |         |     |      |      |     |
|           | 7 |         |       |       |         |     |      |      |     |
|           | 8 |         |       |       |         |     |      |      |     |
|           | 9 |         |       |       |         |     |      |      |     |
| Mean      |   | 170,980 | 6,522 | 1,078 | 127,640 |     |      |      |     |

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## Konečné operace (Closing activities)

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*[Handwritten signature]*  
*[Handwritten date: 25/04/23]*

Erba<sup>®</sup>  
Mannheim

EC 90

QC PROTOCOL

SN 04457 SW: v1.13.3

26/04/2023 07:55

TEST#: 000233  
TYPE: Plasma/Serum  
LEVEL: 2  
LOT: S062244

Na: 169.4 mmol/l OK  
[ 152.0 - 188.0 ]  
K: 6.35 mmol/l OK  
[ 5.81 - 7.13 ]  
Cl: 130.3 mmol/l OK  
[ 116.0 - 140.0 ]

*[Handwritten signature]*  
*[Handwritten date: 26/04/23]*

*[Handwritten signature]*  
*[Handwritten date: 26/04/23]*