



PT EQAS RESULT REVIEW FORM

Issued by QA:

Department: Hematology CBC	Month / Year: May-2023
PT Provider: MHL EQAS	Distribution / Lot No: 230104
Analyte(Test): TLC	
Machine used: Advia 560	Kit used: Advia 560
Test done by: Ritika	Done on: 12-5-2023
Report uploaded by: Amit Kumar	Report uploaded on: 12-5-2023
Software used: MHL PORTAL	
EQAS Report received on 22-5-2023	
Reviewed by: Amit Kumar.	Approved by: DR. Shweta.
Observations: ALL RESULT Received is satisfactory, TLC is shows $\pm 3sd$ deviation due to unit issue after the unit correction value is acceptable range.	

*If results are ungraded, please attach self-assessment. If satisfactory no corrective and preventive action required. But if unsatisfactory, corrective and preventive action is required with optimum evidence.

Self-Assessment Result: Satisfactory Unsatisfactory

Attachment: Yes No

Please investigate failure for unsatisfactory performance and submit a report to QA and Medical Director in 2 Days.

Approved by

Signature
(Lab Director/Designee)



INVESTIGATION FORM FOR UNACCEPTABLE PT/EQAS RESULTS

Issued by QUA:

Department	Haematology		
PT provider	MHL EQAS		
Survey Name:	MHL EQAS 230104		
Date Survey received	07-5-2023	Date Analysis Performed:	12-5-2023
Date Survey Results Submitted	12-5-2023	Date Results received:	22-5-2023

Investigation performed by	AMIT KUMAR	
Lab director	DR. SHWETA	Date 22-5-2023

Unacceptable Result 1			
Specimen No.	03	Analyte	TLC
Reported result	1.5	Intended result/range	1500
Acceptable Lower Limit	775.62	Acceptable Higher Limit	1672.9
SD	224.32	CV%	20.50

Unacceptable Result 2			
Specimen No.	04	Analyte	TLC
Reported result	8.48	Intended result/range	8480
Acceptable Lower Limit	7824.96	Acceptable Higher Limit	11615.04
SD	947.52	CV%	9.75

Unacceptable Result 3			
Specimen No.		Analyte	
Reported result		Intended result/range	
Acceptable Lower Limit		Acceptable Higher Limit	
SD		CV%	

EVALUATION OF POSSIBLE SOURCES OF ERROR

Clerical	Yes	No	N/A
Incorrect transcription of the result from the instrument read-out or report? (Check the raw data.)			✓
Incorrect instrument/method/reagent reported on the result form (Check instrument log book)			✓
Mis-match of the units of measure between the result form and the instrument results			✓
Incorrect decimal placement			✓
Errors in calculations.			✓
Mis-match between result reported on the 'result form' and on the 'proficiency testing evaluation report'			✓
Mis-match between result reported on the 'result form' and on the 'proficiency testing evaluation report'			✓

Any other clerical problem?			✓
<p>A response of "Yes" to any of these questions may indicate a clerical error.</p> <p>Although reporting of PT results is unlike those for patient results, clerical errors may indicate a need for additional staff training, review of instructions provided with PT or investigation of the reporting format provided by the testing device. If results reported on the result form do not match the results found on the evaluation report, please contact PT provider.</p>			

Procedural	Yes	No	N/A
Written procedure not followed			✓
Improper preparation of reagents			✓
Reagents open stability not within acceptable range			✓
Faulty standards run			✓
Unacceptable IQC results (Run accepted in non-linear range/though controls were out of range)			✓
Incorrect interpretation of microscopic examinations			✓
Incorrect dilution or pipetting error			✓
Incorrect staining or interpretation			✓
Time delay between reconstitution and analysis			✓
Media preparation related			✓
Antibiotic disc potency			✓
Any other procedural problem			✓

A response of "Yes" to any of these questions may indicate a procedural error. These errors indicate inappropriate operation of equipment or performance of a method.

A review of the instructions provided with the proficiency testing material and/or review of laboratory procedures may be required.

Analytical	Yes	No	N/A
Most recent calibration unacceptable or not within established stability limits at the time of PT			✓
Instrument repaired or replaced recently at time of PT run			✓
Review of past PT results indicate unevenly distributed data or a bias			✓
Intended result not within the measuring range for the instrument			✓
Was instrument maintenance performed on schedule?	✓	✗	
Review of records indicate there was related instrument/test problems noted prior to or after the PT was performed			✓
Any other analytical problem?			✓

A response of "Yes" to any of these questions may indicate an analytical error. These types of errors could indicate a failure to follow recommended instrument maintenance and calibration.

Specimen handling	Yes	No	N/A
Survey specimens not reconstituted as indicated in the Kit Instructions			✓
Survey specimens not stored as indicated in the Kit Instructions			✓
Special instructions provided in the Kit Instructions not performed as indicated			✓
Correct tests not performed on the correct vial of proficiency testing material			✓
Survey specimen mix-up			✓
Any other problem related to specimen handling?			✓

A response of "Yes" to any of these questions may indicate a specimen handling error. These types of errors could indicate a failure to read the material provided with the Surveys materials.

PT Material	Yes	No	N/A
Late shipment			✓
Hemolysed sample			✓
Bacterial contamination			✓
Perceived survey bias			✓
Poor growth in culture			✓
Unstable PT material			✓
Matrix effect incompatible with method			✓
No comparable peer group			✓
Inappropriate peer group based on method reported on result form			✓
Acceptable range too low			✓
Any other problem?			✓

A response of "Yes" to any of these questions may indicate a problem with the PT material. If a delay in receipt of material in the laboratory is an issue, ensure timely receipt of Surveys after arrival in the institution. If you believe your result was compared to an inappropriate peer group, verify the method reported on the result form. Contact PT provider for additional information if needed.

Miscellaneous	Yes	No	N/A
Was analysis performed by technically competent personnel? (Check technical competence records)	✓		
Were all patient reports satisfactory on day of run? (i.e. there were no patient complaints)	✓		
Were the results of rerun (where possible) the same?	✓		
Any other actions (please specify)			✓

SUMMARY REPORT

Root-cause analysis

UNIT issue found.

Corrective Action proposed

unit changed as required.

Corrective Action taken (To be filled within 30 days of receipt of EQAS result)

Unit changed ~~done~~ after that EQAS TLC result found within the $\pm 2SD$ limits.

Evidence that problem is solved after corrective action (To be attached)

Attaching instrument printout, EQAS Report.

Preventive Action (if any) proposed:

Result will review after or before submitting on the portal by the Lab incharge, Lab head or Lab director.

Conclusion and future plan if no evidence found:

Result will review after or before submitting on the portal by the Lab incharge, manager and Lab head or Director.

Recorded by (Sign / Date)	 22-5-2023
Reviewed by Lab/Medical Director (Sign / Date)	 22-5-2023
Reviewed by Head - QA (Sign / Date)	 22-5-2023

Lab : 880

Lab Name : LIFEWELL DIAGNOSTICS PVT LTD (CHANDIGARH)

Report Date : 19/05/2023

Cycle No : 230104

Sample No : 03

Sample Date : 14/04/2023

Haematology - CBC, Total Leucocytes (WBC) count, /c.mm

Comparative Statistics


	N	Mean	SD	CV%	Excluded
All Methods	167	1,099.66	217.27	19.76	48
Impedance - (Your Method)	142	1,094.34	224.32	20.50	36
ADVIA 560 - (Your Analyzer)	7	1,280.00	0.00	0.00	6

Evaluation Based On - Equipment wise Peer Grouping


Your Result	1.50	SDI / z-Score	*
Assigned Value	1,224.26	Average SDI	0.00
SDPA	87.400	Uncertainty of AV	14.96
Standard Unit Result Val	1.50	%DEV	-99.88

Note : # - Not evaluated / Bimodal distribution
 X : SDI more than 4


Frequency Histogram




Z-Score Trend - Within Round



Frequency Histogram Your Method




Frequency Histogram Your Analyzer



Lab : 880
 Lab Name : LIFEWELL DIAGNOSTICS PVT LTD (CHANDIGARH)
 Report Date : 19/05/2023

Frequency Histogram Your Method



Note : SDI with * marked are not considered for evaluation

Comments :

Haematology - CBC, Total Leucocytes (WBC) count		N	Mean	SD	CV%	Excluded
	All Analyzer	167	1,099.66	217.27	19.76	48
Your Analyzer :	ADVIA 560	7	1,280.00	0.00	0.00	6
	Beckman Coulter DxH 560	24	1,159.33	35.95	3.10	9
	Beckman Coulter DxH 520	19	1,140.63	61.26	5.37	3
	Horiba Yumizen H500	19	671.54	53.05	7.90	6
	Beckman Coulter DxH 500	18	1,153.85	61.17	5.30	4
	Beckman Coulter DxH 500	17	1,153.85	61.17	5.30	4
	Horiba Yumizen H550	12	656.25	109.54	16.69	4
	Sysmex XP-100	9	1,320.00	83.67	6.34	4

Remarks: Interpretation of your reports

- a) If your SDI lies between +/- 2.0 SDI then your results are well within limits.
- b) If your SDI lies between +/-2.0 to 2.99 SDI then it is a warning alert flag.
- c) If your SDI is +/- 3.0 SDI or more then it is an action alert flag.

Interpretation when an analyte is reported for multiple samples in a Round:

a) Reported on 2 samples:

- Two out of two pass - Acceptable performance
- One out of two pass - Review for random error & take action if appropriate
- Both failed - Unacceptable performance, Review for systematic error

b) Reported on 3 samples:

- All pass - Acceptable performance
- Two out of three pass - Review for random error & take action if appropriate
- One out of three pass - Unacceptable performance, Review for systematic error
- All failed - Unacceptable performance, Review for systematic error

Dr. Puneet Kumar Nigam

Dr Puneet Kumar Nigam
 PT coordinator & Technical Manager, MHL EQAS
 Unit No. 409-416.
 Commercial Building - 1A
 Kohinoor Mall, Kiroi Road, Kurla (W),
 Mumbai - 400070

This is a computer generated report hence signature is not required



PC-1022

Lab : 880
Lab Name : LIFEWELL DIAGNOSTICS PVT LTD (CHANDIGARH)
Report Date : 19/05/2023

Cycle No : 230104
Sample No : 04
Sample Date : 14/04/2023

Haematology - CBC, Total Leucocytes (WBC) count, /c.mm

Comparative Statistics

- All Methods
- Impedance - (Your Method)
- ADVIA 560 - (Your Analyzer)

N	Mean	SD	CV%	Excluded
166	8,233.83	900.93	10.94	38
142	8,212.39	914.85	11.14	29
7	9,720.00	947.52	9.75	5


Evaluation Based On - Equipment wise Peer Grouping

Your Result	8.48
Assigned Value	8,617.43
SDPA	841.000
Standard Unit Result Val	8.48


SDI / z-Score	.
Average SDI	0.00
Uncertainty of AV	104.38
%DEV	-99.91

Note : # - Not evaluated / Bimodal distribution
X : SDI more than 4


Frequency Histogram




Z-Score Trend - Within Round



Frequency Histogram Your Method



Frequency Histogram Your Analyzer



Lab : 880
Lab Name : LIFEWELL DIAGNOSTICS PVT LTD (CHANDIGARH)
Report Date : 19/05/2023

Cycle No : 230104
Sample No : 04
Sample Date : 14/04/2023

Frequency Histogram Your Method



Note : SDI with * marked are not considered for evaluation

Comments :

Haematology - CBC, Total Leucocytes (WBC) count

	N	Mean	SD	CV%	Excluded
All Analyzer	166	8,233.83	900.93	10.94	38
Your Analyzer : ADVIA 560	7	9,720.00	947.52	9.75	5
Beckman Coulter DxH 560	23	8,508.67	201.95	2.37	8
Beckman Coulter DxH 520	19	8,315.63	520.99	6.27	3
Horiba Yumizen H500	19	6,641.33	859.95	12.95	4
Beckman Coulter DxH 500	18	8,166.43	229.03	2.81	4
Beckman Coulter DxH 500	17	8,144.62	222.73	2.74	4
Horiba Yumizen H550	12	7,541.11	1,091.00	14.47	3
Sysmex XP-100	9	8,671.43	390.36	4.50	2

Remarks: Interpretation of your reports

- If your SDI lies between +/- 2.0 SDI then your results are well within limits.
- If your SDI lies between +/-2.0 to 2.99 SDI then it is a warning alert flag.
- If your SDI is +/- 3.0 SDI or more then it is an action alert flag.

Interpretation when an analyte is reported for multiple samples in a Round:

a) Reported on 2 samples:

- Two out of two pass - Acceptable performance
- One out of two pass - Review for random error & take action if appropriate
- Both failed - Unacceptable performance, Review for systematic error

b) Reported on 3 samples:

- All pass - Acceptable performance
- Two out of three pass - Review for random error & take action if appropriate
- One out of three pass - Unacceptable performance, Review for systematic error
- All failed - Unacceptable performance, Review for systematic error

Dr Puneet Kumar Nigam

Dr Puneet Kumar Nigam
PT coordinator & Technical Manager, MHL EQAS
Unit No. 409-416.
Commercial Building - 1A
Kohinoor Mall, Kirod Road, Kurla (W),
Mumbai - 400070

This is a computer generated report hence signature is not required



PC-1022



Measure



AS



Database



Print

Sample ID:
Date:

HCBC SAMPLE 3 Patient ID
10/05/2023 09:47 Mode:

WBC:	1.50	LL	10 ³ / μ L	[5.00	-	10.00]	DIFF
NEU:	0.00	LL	10 ³ / μ L	[2.00	-	7.50]	
LYM:	0.69	L	10 ³ / μ L	[1.30	-	4.00]	
MON:	0.06	LL	10 ³ / μ L	[0.15	-	0.75]	
EO:	0.07		10 ³ / μ L	[0.00	-	0.50]	
BAS:	1.33	HH	10 ³ / μ L	[0.00	-	0.75]	
NEU%:	0.0	LL	%	[40.0	-	75.0]	
LYM%:	45.7	H	%	[21.0	-	40.0]	
MON%:	3.7		%	[3.0	-	7.0]	
EO%:	4.6		%	[0.0	-	5.0]	
BAS%:	88.6	HH	%	[0.0	-	1.5]	BAS
RBC:	1.03	LL	10 ⁵ / μ L	[4.00	-	5.50]	
HGB:	3.1	LL	g/dL	[12.0	-	17.0]	
HCT:	12.0	LL	%	[36.0	-	52.0]	
MCV:	116.7	H	fL	[78.0	-	95.0]	
MCH:	29.7		pg	[27.0	-	32.0]	
MCHC:	25.4	L	g/dL	[30.0	-	35.0]	
RDWsd:	100.1	H	fL	[46.0	-	59.0]	
RDWcv:	---	**	%	[0.0	-	16.0]	
PLT:	27	LL	10 ³ / μ L	[150	-	400]	
PCT:	0.02		%		-	15.0]	
MPV:	6.0	L	fL	[8.0	-	15.0]	
PDWsd:	15.7		fL		-		
PDWcv:	40.0		%		-		
PLCR:	19.13		%		-		
PLCC:	5		10 ³ / μ L		-		

Warning flags

ADT

Interpretive flags

Leukopenia?
Neutropenia? Bas

Morphology at 400x

Down

Up

Measure

AS

Database

Print

Sample ID:

HCBC SAMPLE 04 Patient ID:

Date:

10/05/2023 09:50 Mode:

WBC:	8.48		$10^3/\mu\text{L}$	[5.00 - 10.00]
NEU:	0.00	LL	$10^3/\mu\text{L}$	[2.00 - 7.50]
LYM:	1.72		$10^3/\mu\text{L}$	[1.30 - 4.00]
MON:	0.47		$10^3/\mu\text{L}$	[0.15 - 0.70]
EO:	0.25		$10^3/\mu\text{L}$	[0.00 - 0.50]
BAS:	---	**	$10^3/\mu\text{L}$	[0.00 - 0.15]
NEU%:	0.0	LL	%	[40.0 - 75.0]
LYM%:	20.3	L	%	[21.0 - 40.0]
MON%:	5.6		%	[3.0 - 7.0]
EO%:	3.0		%	[0.0 - 5.0]
BAS%:	---	**	%	[0.0 - 1.5]
RBC:	2.15	L	$10^3/\mu\text{L}$	[4.00 - 5.50]
HGB:	7.7	L	g/dL	[12.0 - 17.0]
HCT:	28.3	L	%	[36.0 - 52.0]
MDV:	131.5	H	fL	[76.0 - 96.0]
MCH:	35.9	H	pg	[27.0 - 32.0]
MCHC:	27.3	L	g/dL	[30.0 - 36.0]
RDWsd:	100.2	H	fL	[48.0 - 59.0]
RDWcv:	---	**	%	[0.0 - 16.0]
PLT:	211		$10^3/\mu\text{L}$	[150 - 400]
PCT:	0.16		%	
MPV:	7.5	L	fL	[8.0 - 15.0]
PDWsd:	15.2		fL	
PDWcv:	38.8		%	
PLCR:	19.38		%	
PLCC:	41		$10^3/\mu\text{L}$	

DIFF

BAS

RBC

PLT

Warning flags

Interpretive flags

Morphological flags

uADT

Neutropenia?
Anemia? Macrocy

G

Down

Up

Abundance

Flow



42%



10%