

Sysmex India Pvt. Ltd.
 1002, Danji Shamji Business Galleria,
 10th Floor, LBS Marg,
 Kanjurmarg (West),
 India, 400078
 Tel. (+91) 22 6112 6666 Fax. (+91) 22 2577 6790
 Email: sysmex@sysmex.co.in
 Web: <https://www.sysmex.co.in>
 GST Reg. No. 27AADCS1551J1ZC Co. Reg. No: ABC



CSR No.: WO-100263467
 Revision No.:

Customer Service Report

Customer Name	Sri Vishnu Healthcare Unit		
Contact Person	E.Mahipal		
Installation Date	30/01/2021	Login Date/Time	
Equipment ID	10133465	FSR Name	M Prathik Goud
Model	XN-350	From Date/Time	28/01/2023 05:01:02 PM
Serial No.	16464	To Date/Time	28/01/2023 07:54:00 PM
Total Cycle Count	18,181	Case No.	

Defect Analysis	
Issue Code	Cause Code
Issue	Cause
Calibration	Calibration

Action Taken	Action Code	Value Before Adjustment	Value After Adjustment
Test Calibration Task	Completed		

Remarks

Items Consumed				
Product ID	VPN	Product Description	QTY	UOM

Name:- M Prathik Goud	Acknowledged By : E.Mahipal
Sign: <i>M. Prathik</i>	Sign : <i>malj</i>
Date: 29/01/2023 05:45:53 PM	Date: 29/01/2023 05:45:53 PM

Customer Confidential

Sysmex India Pvt. Ltd.

1002, Damji Shamji Business Galleria,
10th Floor, LBS Marg,
Kanjurmarg (West),
India 400078
Tel: (+91) 22 6112 6666 Fax: (+91) 22 2577 6790
Email: sysmex@sysmex.co.in
Web: <https://www.sysmex.co.in>
GST Reg No: 27AADCS1551J1ZC Co. Reg. No: ABC



CSR No.: WO-100263470
Revision No.:

Customer Service Report

Customer Name	Sri Vishnu Healthcare Unit		
Contact Person	E.Mahipal		
Installation Date	30/01/2021	Login Date/Time	
Equipment ID	10133465	FSR Name	M Prathik Goud
Model	XN-350	From Date/Time	28/01/2023 02:20:05 PM
Serial No.	16464	To Date/Time	28/01/2023 04:59:00 PM
Total Cycle Count	18,181	Case No.	

Defect Analysis

Issue Code	Cause Code
Issue	Cause

Action Taken	Action Code	Value Before Adjustment	Value After Adjustment
1. Pre-check: Perform QC analysis with XN-L Check or XN Check at least 2 times	Completed		
2. RBC Detector Block: Open the metal cover of the RBC detector block and clean the internal surfaces with moist paper towel	Completed		
3. RBC Detector Block: Perform "Drain Isolation Chamber" in Maintenance Menu of IPU program	Completed		
4. RBC Detector Block: Clean transducer aperture using aperture brush dipped in cellclean	Cleaned		
5. RBC/HGB Mix Chamber: Activate "Drain Reaction Chamber". Dispense 2ml of cellclean auto into the RBC chamber and activate sequence again to drain cellclean. Repeat twice	Completed		
6. WBC Chamber: Activate "Drain Reaction Chamber". Dispense 2ml of cellclean auto into the reaction chambers and activate sequence again to drain cellclean. Repeat twice	Completed		
7. FCM Unit: Check laser current value in "Service Menu" of IPU program [Range: 100mA to 140 mA]	Checked	129.5	
8. Vacuum Line: Clean vacuum line using ethanol or methanol (Approx 30ml)	Cleaned		
9. Shut Down: Shut down analyzer before proceeding	Completed		
10. Vacuum Line: Rinse vacuum trap chamber using DI or distilled water	Completed		
11. RBC/HGB Mix Chamber: Remove Filter No. 13 below RBC/HGB chamber and flush it with cellclean using syringe	Completed		
12. WBC Chamber: Remove Filter No. 13 below WBC reaction chamber and flush it with cellclean using syringe	Completed		
13. Waste Chamber: Remove waste chamber and clean the internal surface and the waste sensor with gauze dipped in diluted cellclean.	Cleaned		
14. Waste Chamber: Rinse waste chamber trap chamber with DI or distilled water	Completed		
15. Waste Chamber: Replace the motorized pinch valve pharmed tubing (Total Length : 150mm) (3 Years or 60K Cycle)	NA		
16. FCM Unit: Clean flowcell using syringe and 3ml of cellclean in a back and forth motion at least 20 times	Cleaned		

Systemx India Pvt. Ltd.

1002, Danji Shamji Business Galleria,
10th Floor, LBS Marg,
Kanjurmarg (West),
India, 400078
Tel. (+91) 22 6112 6666 Fax. (+91) 22 2577 6790
Email: systemx@systemx.co.in
Web: https://www.systemx.co.in
GST Reg. No: 27AADCS1551J1ZC Co. Reg. No: ABC



CSR No.: WO-100263470
Revision No.:

Customer Service Report

Customer Name	Sri Vishnu Healthcare Unit		
Contact Person	E.Mahipal		
Installation Date	30/01/2021	Login Date/Time	
Equipment ID	10133465	FSR Name	M Prathik Goud
Model	XN-350	From Date/Time	28/01/2023 02:20:05 PM
Serial No.	16464	To Date/Time	28/01/2023 04:59:00 PM
Total Cycle Count	18,181	Case No.	
Action Taken	Action Code	Value Before Adjustment	Value After Adjustment
17. Piercer: Clean exterior of pipette/piercer using gauze dipped in cellclean	Cleaned		
18. Piercer: Check piercer, replace if it is blunt. Replace piercer. Reset Piercer cycle count if it is replaced (120K Cycles)	Checked		
19. Piercer: Clean rinse cup using cellclean then rinse with DI or distilled water	Cleaned		
20. Piercer: Clean pipette/piercer protector using cellclean then rinse with DI or distilled water	Cleaned		
21. Sheath syringe: Check sheath syringe for leakage. Replace if there is a leak	Checked		
22. Sheath syringe: Remove old grease and lubricate the sheath motor spindle using Kluber L60 (6 Months)	Lubricated		
23. WB pump: Check WB pump for leakage. Replace WB pump if there is a leak	Checked		
24. Sampler: Remove old grease and lubricate the sampler drive unit (XN-530/550) with Grease Molykote (6 Months)	NA		
25. Pneumatic unit: Replace air pump. Reset airpump cycle count if it is replaced (2 Years or 30K Cycles)	NA		
26. Power On: Power on analyzer in service mode and ensure there are no errors. Run auto rinse once and ensure background check pass	Completed		
27. Pneumatic unit: Check 0.06 Mpa pressure from IPU software [Target: ≥ 0.06MPa]	Checked	0.0612	
28. Pneumatic unit: Check -0.03 Mpa pressure from IPU software [Target: ≥ -0.03MPa]	Checked	-0.0303	
29. Mechanical Adjustment: Check pipette/piercer position. Adjust if it is out of alignment	Checked		
30. Mechanical Adjustment: Check sampler hand and tube holder position for XN-550. Adjust if it is out of alignment	Checked		
31. Mechanical Adjustment: Check Barcode Reader functionality. Adjust barcode reading position if it is out of alignment. (XN-530/550)	NA		
32. Electronics Adjustment: Check blood aspiration sensor blank. Adjust if it is out of range [Range: 4800 to 5200]	Adjusted	7167	4941
33. Electronics Adjustment: Check blood aspiration sensor span using XN-L Check or XN Check Level 2. Adjust if it is out of range [Range: 12,500 to 13,500]	Adjusted		13057
34. Electronics Adjustment: Check RBC Clog Level. Adjust if it is out of range [Range: 99 to 101]	Adjusted	100	
35. Electronics Adjustment: Check HGB Blank. Adjust if it is out of range [Range: 4800 to 5200]	Adjusted	7638	5000
36. Operation check: Test the host communication	Completed		



Customer Service Report

Customer Name	Sri Vishnu Healthcare Unit		
Contact Person	E.Mahipal		
Installation Date	30/01/2021	Login Date/Time	
Equipment ID	10133465	FSR Name	M Prathik Goud
Model	XN-350	From Date/Time	28/01/2023 02:20:05 PM
Serial No.	16464	To Date/Time	28/01/2023 04:59:00 PM
Total Cycle Count	18,181	Case No.	

Action Taken	Action Code	Value Before Adjustment	Value After Adjustment
37. Operation check: Perform QC analysis with XN-L Check or XN Check at least 3 times	Completed		
38. Operation check: Perform full backup and copy mmb files to desktop	Completed		
39. General Cleaning: Clean Main Unit covers	Cleaned		
40. General Cleaning: Clean IPU touch panel	Cleaned		
41. General Cleaning: Clean peripheral device (printer, etc.)	Cleaned		
42. General Cleaning: Clean the interior of sampler for XN-550	Cleaned		

Remarks

Items Consumed				
Product ID	VPN	Product Description	Qty	UOM

Name:- M Prathik Goud	Acknowledged By : E.Mahipal
Sign: 	Sign : 
Date: 29/01/2023 05:24:09 PM	Date: 29/01/2023 05:24:09 PM

Customer Confidential

Calibrator Calibration History

Instrument Nickname: XN-L
Material: XN CALCalibration Date: 2023/01/28 18:42:00
Lot No.: 30232101

Logon Name: sysmex

	WBC	RBC	HGB	HCT	PLT
Target	7.535	4.372	11.83	34.71	247.1
No. 2	7.53	4.40	11.8	34.8	236
No. 3	7.59	4.39	11.9	34.9	240
No. 4	7.52	4.37	11.9	34.7	236
No. 5	7.37	4.42	11.8	35.0	237
No. 6	7.58	4.39	11.9	34.9	235
No. 7	7.73	4.36	11.8	34.6	238
No. 8	7.48	4.39	11.9	34.9	243
No. 9	7.45	4.32	11.8	34.2	242
No. 10	7.60	4.36	11.8	34.7	244
No. 11	7.32	4.43	12.0	35.1	242
Range Value	0.41	0.11	0.2	0.9	9
Max Range	0.58	0.13	0.2	1.0	26
Mean Value	7.517	4.383	11.86	34.78	239.3
Delta Percent (%)	0.24	0.25	0.25	0.20	3.26
Acceptable Limit (%)	2.27	1.25	0.78	2.64	4.16
Service Limit (%)	14.00	4.00	5.00	5.00	10.00
Current Rate (%)	102.3	98.9	99.0	101.2	102.9
New Rate (%)	102.5	98.7	98.7	101.2	106.3


Sample No.: QC-23501401
 Patient ID:
 Name:
 Sample Comment:

Ward: Adapter:


Pos.: 2023/01/28 19:40:30 WB
 Doctor:
 Birth: Sex:
 Nickname: XN-L

WBC	2.36	[10 ³ /uL]		
RBC	2.32	[10 ⁶ /uL]		
HGB	5.8	[g/dL]		
HCT	16.9	[%]		
MCV	72.8	[fL]		
MCH	25.0	[pg]		
MCHC	34.3	[g/dL]		
PLT	57	[10 ³ /uL]		
RDW-SD	49.5	[fL]		
RDW-CV	18.5	[%]		
PDW	9.0	[fL]		
MPV	10.5	[fL]		
P-LCR	26.1	[%]		
PCT	0.06	[%]		
NEUT	1.07	[10 ³ /uL]	45.4	[%]
LYMPH	0.64	[10 ³ /uL]	27.1	[%]
MONO	0.24	[10 ³ /uL]	10.2	[%]
EO	0.27	[10 ³ /uL]	11.4	[%]
BASO	0.14	[10 ³ /uL]	5.9	[%]
IG	0.29	[10 ³ /uL]	12.3	[%]


WDF




WDF-CBC



RBC



PLT



WBC IP Message

RBC IP Message

PLT IP Message


Sample No.: QC-23501402
 Patient ID:
 Name:
 Sample Comment:

Ward: Adapter:


Pos.: 2023/01/28 19:42:49 WB
 Doctor:
 Birth: Sex:
 Nickname: XN-L

WBC	7.19	[10 ³ /uL]		
RBC	4.31	[10 ⁶ /uL]		
HGB	12.7	[g/dL]		
HCT	35.9	[%]		
MCV	83.3	[fL]		
MCH	29.5	[pg]		
MCHC	35.4	[g/dL]		
PLT	266	[10 ³ /uL]		
RDW-SD	45.6	[fL]		
RDW-CV	14.9	[%]		
PDW	7.8	[fL]		
MPV	9.2	[fL]		
P-LCR	14.1	[%]		
PCT	0.25	[%]		
NEUT	3.58	[10 ³ /uL]	49.8	[%]
LYMPH	1.61	[10 ³ /uL]	22.4	[%]
MONO	0.61	[10 ³ /uL]	8.5	[%]
EO	0.92	[10 ³ /uL]	12.8	[%]
BASO	0.47	[10 ³ /uL]	6.5	[%]
IG	1.01	[10 ³ /uL]	14.0	[%]


WDF




WDF-CBC



RBC



PLT



WBC IP Message

RBC IP Message

PLT IP Message


Sample No.: QC-23501403
 Patient ID:
 Name:
 Sample Comment:

Adapter:
 Ward:


Pos.: 2023/01/28 19:44:39 WB
 Doctor:
 Birth: Sex:
 Nickname: XN-L

WBC	16.47	[10 ³ /uL]		
RBC	5.36	[10 ⁶ /uL]		
HGB	16.8	[g/dL]		
HCT	46.9	[%]		
MCV	87.5	[fL]		
MCH	31.3	[pg]		
MCHC	35.8	[g/dL]		
PLT	644	[10 ³ /uL]		
RDW-SD	47.0	[fL]		
RDW-CV	15.1	[%]		
PDW	8.1	[fL]		
MPV	8.9	[fL]		
P-LCR	12.2	[%]		
PCT	0.58	[%]		
NEUT	8.51	[10 ³ /uL]	51.7	[%]
LYMPH	3.11	[10 ³ /uL]	18.9	[%]
MONO	1.43	[10 ³ /uL]	8.7	[%]
EO	2.23	[10 ³ /uL]	13.5	[%]
BASO	1.19	[10 ³ /uL]	7.2	[%]
IG	2.24	[10 ³ /uL]	13.6	[%]


WDF




WDF - CBC



RBC



PLT



WBC IP Message

RBC IP Message

PLT IP Message