# ERBA- H560 AUTOMATED HEMATOLOGY ANALYZER

# PERFORMANCE QUALIFICATION

For

**ERBA H-560** 

66

Marketed by:
Transasia Bio-Medicals Ltd.,
(ISO 9002 CERTIFIED)
Transasia House,
Chandivali Studio road,
Andheri (E),
MUMBAI – 400 072

## **Table of Contents:**

Sr. No.	Contents	Page No.
I	Approval of the PQ procedure	3
II	Instructions	4
III	Scope	5
IV	Performance Qualification	6
V	System Certification	10

## I. Approval of the PQ procedure

Both RARI and Transasia are jointly responsible for conducting the Performance Check of the Hematology Analyzer, Model: ERBA – H560, Serial No. K1104B2211056 in the clinical lab of the Hospital as per the attached protocol.

Date:

0 0

Protocol Performed By: Transasia Representative

Name: MR. AJAY SINGH

Title: PERFORMANCE QUALIFICATION

Company: TRANSASIA BIO-MEDICALS LTD.

Signature:

Validation Team from

Name

Designation : Department :

Name

Designation
Department

**Customer Authorizations:** 

Name: DR. MEENAKSHI SURI

Title: PERFORMANCE QUALIFICATION

Site : RARI Jammy.

Name : Title: PERFORMANCE QUALIFICATION

Site :

Signature:

Signature:

Date:

Date:

#### II. Instructions

1. An authorized TRANSASIA representative will check for the performance of the instrument and enter the specific data as outlined in the Performance Qualification. Each result will be noted and dated.

2. Performance checks on a regular basis described in the Further Performance

Checks (vide-infra) will be the responsibility of the customer's personnel.

3. Employee of (Customer) RARI, Jammu Laboratory will verify each result and sign on the last page. The members of the validation team will carry this out.

4. ALL deviations from the acceptance criteria detailed in this document will be noted in the COMMENTS section at the end of each PQ protocol. All resolutions to such problems will also be noted in the COMMENTS section and must be resolved prior to the issuance of a SYSTEM CERTIFICATION. These will be an additional cost to the purchasing institution (CUSTOMER). However, this additional cost will be waived when this test is conducted at the time of the initial performance check of new instruments.

5. Any test data that does not meet the specified acceptance criteria will be submitted to the appropriate laboratory personnel for the solution. All steps taken

subsequently will be documented.

6. This document contains proprietary information and is in no way to be copied, photographed, or duplicated in any way without expressed written authorization by the Production Manager at Transasia Bio-Medicals Ltd., Transasia House, Mumbai.

Validation Team:

Name

Mr. Afay Singh.
Sv. Service Engg.

Signature

#### III. Scope

This Performance Qualification protocol will be performed on the Hematology Analyzer, Model ERBA -H560, Serial No K1104B2211056 located in the clinical laboratory of RARI Jammu. This Protocol will define the documentation that will be used to evaluate the instrument's installation in accordance with the manufacturer's specifications and intended use. Successful completion of this protocol will verify that the instrument identified is performing in accordance with the intended usage.

Trained, knowledgeable personnel will perform qualification studies.

Any exceptional conditions encountered during the qualification studies will be identified for review. Exceptional conditions will be investigated and the appropriate course of action determined. All documents will be initialed and dated.

Validation Team:

Name

Mr. Afgy Sings.
Sv. Service Engg

Designation

Signature

# IV. Performance Qualification

a. Instrument Identification		Verified Date		
Model Name     Serial Number	ERBA – H560 K1104B2211056	9/1/23		
b. Following is a list of test	s to be performed and veri	ned:		
Test No. Test Nam		W. I Dodo		
02 Sample Processing	Ability to proce	ess samples		
03 Further Performance	e Checks Regular Main	ntenance NA		
Signature	Agay Single Service Ex Busles.	78		
Name  Designation	Agay Single Service Cy Dougles.	7 °		

**Verified Date** 

## c. Performance Testing

Test 1

Test Name:

Sample Processing

Purpose:

**Ability to Process Samples** 

Method:

## 1. Run the control samples five times consecutively

Acceptance Criteria: Each of the results obtained above should be within the range as specified in the control chart.

## **Parameters Values for Verification:**

#### **RBC Count:**

PD .	General Values	Results Obtained	Pass	Fail
Test	Control Values	2.23	PASS	
1.	2.31	2.23	PASS	
2.	2.31		PASS	
3.	2.31	2.22	PASS	_
4.	2.31	2.21		
5	2.31	2.24	PASS	

Validation Team:

Name

Designation

Mr. Afay Single. Sv. Service Evgg.

Signature.

#### **WBC Count:**

	1 77 1	Results Obtained	Pass	Fail
Test	Control Values		PASS	
1.	3.42	3.42	PASS	
2	3.42	3.41	PASS	
3	3.42	3.38		
3.	3.42	3.32	PASS	
4.		3.60	PASS	
5.	3.42			

Hemoglobin:

Hemog	globin:	D 14 Obtained	Pass	Fail
Test	Control Values	Results Obtained		
1	5.7	5.7	PASS	
1.		5.7	PASS	
2.	5.7		PASS	
3.	5.7	5.6	PASS	
4.	5.7	5.7		
7.		5.7	PASS	
5.	5.7		MANAGEMENT OF THE STATE OF THE	

HCT

TG . 137.1	Results Obtained	Pass	Fail
17.3			
17.3			
17.3	18.1		
	17.9		
	18.2	PASS	
	Control Values 17.3 17.3 17.3 17.3 17.3	17.3     18.1       17.3     18.1       17.3     18.1       17.3     17.9	Control varies       18.1       PASS         17.3       18.1       PASS         17.3       18.1       PASS         17.3       17.9       PASS         17.3       17.9       PASS

Platelet Count:

Platele	et Count:	7 1 01/1 1	Pass	Fail
Test	Control Values	Results Obtained		1 an
1	67	71	PASS	
1.	67	72	PASS	
2.	67	68	PASS	
3.	67		PASS	
4.	67	66	PASS	
5.	67	67	ITADO	

Validation Team: TRANSASIA BIO-MEDICALS LTD

Name.

MR. AJAY SINGH

Designation

SR. SERVICE ENGG.

Signature

#### Test 2

#### Test Name:

- 1. Tests for checking the performance of the instruments during analysis
- 2. Tests for checking the long-term performance of the instrument

#### Purpose:

The purpose of the above checks is to ensure the reliability of the results being obtained.

#### Method:

1. During Sample analysis:

To run control samples each time the instrument is used for sample analysis and verification of the results of the controls to be within the reference range to be established by the performance of the precision experiments.

2. Long-term Performance

This is to be checked by Levy Jennings plots to be updated once in six months

Mv. Afgy Singh.
Sv. Service Enge.

Validation Team:

Name

Designation

Signature

## V. System Certification

Study data has determined that the system described in this document either meets all criteria outlined in this Performance Qualification Protocol, or exceptional conditions have been identified and documentation included. Exceptional conditions, if any, have been addressed. The system is ready for specified usage.

## Report Performed By: Transasia Representative

Name: MR. AJAY SINGH

Title: PERFORMANCE QUALIFICATION

Company: TRANSASIA BIO-MEDICALS LTD.

Signature:

Date:

### Customer Authorizations:

Name: DR. MEENAKSHI SURI

Title: PERFORMANCE QUALIFICATION

Site: JAMMU

Signature:

Date:

Name:

Title: PERFORMANCE QUALIFICATION Signature:

Site: Date: