Customer: Vioma Lifesciences Pvt Ltd

Page 1 of 8

Equipment: BC 760 Auto Hematology Analyzer

Title: IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

## Installation Qualification for Mindray BC 760 Sr. No.:- DG7-32000377

Customer Name : Vioma Lifesciences Pvt Ltd

Contact Person : Mrs Shilpa Surpur

Instrument Model : BC 760

Serial No. : DG7-32000377

Date Of installation : 22/06/2023

The instrument was installed And was found to be working satisfactory. Preliminary Customer Training was provided And standardization of the parameters Were done. The results were found to be within the expected range and system found to be working satisfactorily.

**Customer: Vioma Lifesciences Pvt Ltd**

Service Engineer		Customer	
Name:	Nandkumar Zagade	Name:	Dr Niraj Gujar
Signature:		Signature:	

Customer: Vioma Lifesciences Pvt Ltd

Page 2 of 8

Equipment: BC 760 Auto Hematology Analyzer

Title: IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

## Installation Certificate For BC 760 Sr. No.- DG7-32000377

This is to certify that the instrument BC 760 Sr. No:- DG7-32000377 is successfully installed and commissioned at :

**Vioma Lifesciences Pvt Ltd**

Indofill Industries Ltd , Off Swami Vivekanand Road , Azad Nagar , Sandoz Baug, P.O. Thane West , Maharashtra-400607

and the Installation Protocol/ checklist has been Successfully completed for the above Instrument.

Date of Installation : 22/06/2023

MBQ , Technical Services Department

Service Engineer		Customer	
Name:	Nandkumar Zagade	Name:	Dr Niraj Gujar
Signature:		Signature:	

Customer: Vioma Lifesciences Pvt Ltd

Page 3 of 8

Equipment: BC 760 Auto Hematology Analyzer

Title: IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

## Installation Qualification for Mindray BC 760 Sr. No.:- DE5-230000074

Customer Name : Vioma Lifesciences Pvt Ltd

Contact Person : Mrs Shilpa Chati

Instrument Model : BC 760

Serial No. : DG7-32000377

Date Of installation : 22/06/2023

**Initial Inspection of the unit carried out and the details are as follows:**

### System Condition Report:

Found the system to have been delivered in satisfactory condition and no external physical damaged observed on the same, Package was kept in good Condition as per the directional indicators like not tilt, indicating the system has not been subjected to mechanical shocks or stored in any manner, so as to cause any damage to the same.

Customer: Vioma Lifesciences Pvt Ltd

Page 4 of 8

Equipment: BC 760 Auto Hematology Analyzer

Title: IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

**Found all the required accessories as required.**

**Installation Procedure And Checklist Attached for the records.**

**External Requirement for Installation:**

1. Input voltage of 220V-240V/50Hz or 60 Hz.
2. Recommended operating Temperature is 15-30 degree Celsius, with in Relative Humidity 30-85% and Atmospheric pressure 70-106kPa.

**Installation Qualification for BC 760 Sr. No.:- DG7-32000377**

**Carried out all the installation procedures as per the installation procedure and checklists.**

**Carried out all the necessary checks and alignments.**

**Carried out all the necessary system checks and tests.**

**Handover the Instrument for operators Training And Qualifications**

**Technical Services Department**

Service Engineer		Customer	
Name:	Nandkumar Zagade	Name:	Dr Niraj Gujar
Signature:		Signature:	

**Customer:** Vioma Lifesciences Pvt Ltd

Page 5 of 8

**Equipment:** BC 760 Auto Hematology Analyzer

**Title:** IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

## **Performance Qualification for Mindray BC 760 Sr. No.: -DG7-32000377**

### **Calibration Parameters**

Checked and found all the control [parameters to be within the acceptable CV limits and in range.

Checked and found all controls to be within the acceptable SD.

### **System Certification:**

Study data has determined that the system described in this document either meets the necessary criteria outlined in this Performance Qualification Protocol, or exceptional conditions have been identified and documentation included.

**The system is ready for Specific Usage.**

Customer: Vioma Lifesciences Pvt Ltd

Page 6 of 8

Equipment: BC 760 Auto Hematology Analyzer

Title: IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

## Operational Qualification:

### System Certification:

Study data has determined that the system described in this document either meets all criteria outlined in this Operational Protocol, or exceptional conditions have been identified and documentation included.

Exceptional Conditions, If any Have been Addressed.

**The System is ready for specific usage.**

Customer Authorization: **Vioma Lifesciences Pvt Ltd**

**Indofill Industries Ltd , Off Swami Vivekanand Road , Azad Nagar ,  
Sandoz Baug,P.O.Thane West , Maharashtra-400607**

Service Engineer		Customer	
Name:	Nandkumar Zagade	Name:	Dr Niraj Gujar
Signature:		Signature:	

Customer: Vioma Lifesciences Pvt Ltd

Page 7 of 8

Equipment: BC 760 Auto Hematology Analyzer

Title: IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

## Operations Qualifications for BC 760 Sr. No.- DG7-32000377

1. Verified all the Mechanical Movements : Done
2. Verified Hydraulics : Done
3. Verified Electrical Systems : Done
4. Verified the all reagents Systems : Done

### Technical Services Department

Service Engineer		Customer	
Name:	Nandkumar Zagade	Name:	Dr Niraj Gujar
Signature:		Signature:	

Customer: Vioma Lifesciences Pvt Ltd

Page 8 of 8

Equipment: BC 760 Auto Hematology Analyzer

Title: IQ, OQ & PQ DOCUMENT

MBQ/NKZ/031072023/BC760

## Instrument Setup

1. Assembled the instrument accessories.
2. Removed the shipping Clamps.
3. Connected the LD LYSE REAGENT .
4. Connected the LH LYSE REAGENT
5. Connected the FD DYE REAGENT
6. Connected the DILUENT REAGENT.
7. Connected the Waste Tubing
8. Connected the Power cord and connection cord.
9. Initialise the machine and follow the installation procedure.

## Operational Instection

1. Checked and found Mechanical Movements OK.
2. Checked and found Hydraulics OK.
3. Checked and found Electricals OK.
4. Checked with Controls And Samples, Results Are found OK.

Service Engineer		Customer	
Name:	Nandkumar Zagade	Name:	Dr Niraj Gujar
Signature:		Signature:	

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## Directions

- (1) The operators fill the results in the appropriate tables refer to the template.
- (2) In the result row, fill the form with test data in digital or OK/NG, fill the conclusion row in PASS or FAIL.
- (3) For inspection items not applicable to certain configurations, mark checkbox with "No", and leave the "Conclusion" cell blank.
- (4) the item which carry "\*" are just for reference, not the criterion for the final result.