

To Whom It May Concern

For ISO 15189:2012 and ISO 15189:2014 accredited Laboratories — requirements regarding "Calibration & Verification Procedures" [1]

All In vitro Diagnostics Products which are manufactured and distributed by Roche Diagnostics GmbH and for which a Free-Sales-Certificate is issued, are CE-marked.

The In-Vitro-Diagnostics Directive of the European Union [2A.] which is currently switching to IVD Regulation 2017/746/EU (final timeline: May 26, 2022) [2B.] requires for all CE marked products that the manufacturer assures compliance of the products with the requirements of the mentioned directive or regulation. This means that all processes in development and manufacturing of Roche Diagnostics GmbH products are guided by a Quality Management System. Our Quality Management System is in compliance with the requirements from ISO 13485:2016 [3] and 21 CFR Part 820 [4].

The mentioned regulations and standards require that the production systems and measuring devices used are qualified and the manufacturing and test procedures are validated. This status has to be assured by scheduled maintenance and by regular qualification resp. validation reviews and updates.

All physical quantities, calibrators and controls used in Roche Diagnostic systems are fully traceable to certified standards or reference materials. The performance of all In-vitro diagnostics systems of Roche Diagnostics GmbH at the customer site is assured if regular Quality Control measurements, cleaning and maintenance procedures as described in the instructions for use or service documentation are performed. By having controlled internal procedures and by running the tasks required in the respective user documentation, all In-vitro diagnostics systems of Roche Diagnostics GmbH will be performed as specified during their defined lifetime.

Additional calibration or verification procedures are NOT required by the user in order to assure the specified performance of every system of Roche Diagnostics GmbH. Only if a user deviates from these manufacturer's recommendations, the user have to establish site-specific calibration and verification procedures as part of his accreditation process.

- [1] ISO 15189:2012/ ISO 15189:2014 Medical laboratories — Requirements for quality and competence
- [2] A. Directive 98/79/EC of the European Parliament and of the Council of the 27 October 1998 on vitro diagnostics medical devices;
B. IVD Regulation 2017/746/EU of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU
- [3] EN ISO 13485:2016 Medical devices — Quality management systems-Requirements for regulatory purposes
- [4] CFR Part 820, Quality System regulations 21 Regulations on medical devices


Mannheim, 10. August 2021

Sincerely,

Roche Diagnostics GmbH

i.V./on behalf of the company

ppa/on behalf of the company


ECA5294AC4E94AF...

Andrea Weber
Manager Global Regulatory Affairs
Centralised and Point of Care Solutions


A7F0BA9FE91A46A...

Ralf Zielenski
Head Q&R Compliance, PRRC RDG
Centralised and Point of Care Solutions

Roche Diagnostics GmbH
Sandhofer Straße 116
D-68305 Mannheim



Operational Qualification for cobas® c311

Notice:

The steps described in OQ.1 have to be carried out after a new system installation and after any repair action which requires additional calibration.

If the service action does not affect the measurement performance, perform only steps OQ.2 and OQ.3 of the Operation Qualification.

Description

OQ.1 Calibration

Calibration of all photometric parameters successful (attached printout)

Calibration of all ISE parameters successful (attached printout)

OQ.2 Quality Control

Specify the type of control used:
 PCCC1

QC of all photometric parameters within acceptable range (see attached results)

QC of ISE parameters within acceptable range (see attached results)

OQ.3.1 Accuracy check for ISE

Perform test with analytical reagents

			Number of det.
Na	ACN	989	15
K	ACN	990	15
Cl	ACN	991	15

Sample solution: PNU (code 300 / Cat. No. 10171735, 10171743, 10651257).
Fill 21 Hitachi cups with PNU (code: 300) and perform Na, K and Cl tests. Calculate the CV.

Accuracy check for ISE was within acceptable range



OQ.3.2 Accuracy check for Photometric Assays

Perform test with analytical reagents

	Number of det.
2-point/end-point Assay	21
Rate A Assay	15

Sample solution: PNU (code 300 / Cat. No. 10171735, 10171743, 10651257).
Fill 21 Hitachi cups with PNU (code: 300) and perform 21 determinations of each parameter.

Accuracy check for Photometric Assays was within acceptable range yes

Deviation Report: Any discrepancies found during the installation must be documented in the space below. Roche personnel will then investigate the deviation and decide upon the most appropriate action to be taken.

Deviation #1	
Investigation	
Action taken	
Deviation resolved satisfactorily?	specify



Attachments

Precision Check

17/12/22

19:39

TEST	N	MEAN	UNIT	RANGE	MAX.	MIN.	SD	CV(%)
Cl	00015	116.41	mmol/L	1.0	116.9	115.9	0.27	0.23
CREJ2	00015	4.143	mg/dL	0.14	4.23	4.09	0.036	0.88
K	00015	7.721	mmol/L	0.08	7.75	7.67	0.019	0.24
Na	00015	149.3	mmol/L	1	150	149	0.5	0.31

Precision Check

17/12/22

19:42

TEST	N	MEAN	UNIT	RANGE	MAX.	MIN.	SD	CV(%)
ASTL	00021	162.24	U/L	2.9	163.8	160.9	0.69	0.43
GLUC3	00021	252.94	mg/dL	2.5	254.1	251.6	0.64	0.25

Data Monitor

17/12/22 18:07

Ser/P1 N000014 001 ID Pcheck12
17/12/22 ASTL GLUC3
17:50:27 162.2 252.3
bmserv

Ser/P1 N000015 001 ID Pcheck13
17/12/22 ASTL GLUC3
17:50:51 163.4 253.2
bmserv

Ser/P1 N000016 001 ID Pcheck14
17/12/22 ASTL GLUC3
17:51:15 162.6 253.3
bmserv

Ser/P1 N000017 001 ID Pcheck15
17/12/22 ASTL GLUC3
17:51:39 162.4 253.6
bmserv

Ser/P1 N000018 001 ID Pcheck16
17/12/22 ASTL GLUC3
17:52:03 162.8 252.9
bmserv

Ser/P1 N000019 001 ID Pcheck17
17/12/22 ASTL GLUC3
17:52:27 163.8 253.2
bmserv

Ser/P1 N000020 001 ID Pcheck18
17/12/22 ASTL GLUC3
17:52:51 162.8 254.1
bmserv

Ser/P1 N000021 001 ID Pcheck19
17/12/22 ASTL GLUC3
17:53:15 162.6 253.7
bmserv

Ser/P1 N000022 001 ID Pcheck20
17/12/22 ASTL GLUC3
17:53:39 162.9 252.6
bmserv

Ser/P1 N000023 001 ID Pcheck21
17/12/22 ASTL GLUC3
17:54:03 162.2 253.4
bmserv

Data Monitor

17/12/22

18:07

Ser/Pl N000003 001 ID Pcheck1
17/12/22 ASTL GLUC3
17:46:51 160.9 252.8
bmserv

Ser/Pl N000004 001 ID Pcheck2
17/12/22 ASTL GLUC3
17:47:09 161.6 252.3
bmserv

Ser/Pl N000005 001 ID Pcheck3
17/12/22 ASTL GLUC3
17:47:28 161.7 253.1
bmserv

Ser/Pl N000006 001 ID Pcheck4
17/12/22 ASTL GLUC3
17:47:47 161.9 251.6
bmserv

Ser/Pl N000007 001 ID Pcheck5
17/12/22 ASTL GLUC3
17:48:06 161.4 252.1
bmserv

Ser/Pl N000008 001 ID Pcheck6
17/12/22 ASTL GLUC3
17:48:25 161.3 252.6
bmserv

Ser/Pl N000009 001 ID Pcheck7
17/12/22 ASTL GLUC3
17:48:45 162.1 252.4
bmserv

Ser/Pl N000010 001 ID Pcheck8
17/12/22 ASTL GLUC3
17:49:03 162.1 253.3
bmserv

Ser/Pl N000011 001 ID Pcheck9
17/12/22 ASTL GLUC3
17:49:22 161.8 252.2
bmserv

Ser/Pl N000012 001 ID Pcheck10
17/12/22 ASTL GLUC3
17:49:41 162.2 253.6
bmserv

Ser/Pl N000013 001 ID Pcheck11
17/12/22 ASTL GLUC3
17:50:03 162.3 253.5
bmserv

Data Monitor

19/09/23

08:49

C004026 097
19/09/23 RF-II
08:35:20 51.4
bmserv

LOT 62690000

RFCO2

Data Monitor

19/09/23

08:54

C009026 094	LOT 52520500	PCCC1			
19/09/23	ALB2	ALP2L	AMYL2	ASTL	DBIL2
08:35:20	3.14	89	83	46.0	0.820
bmserv			ReagEX		
	BILT3	CA2	CHO21	CREJ2	GGT12
	1.008	8.81	94.4	0.98	58
	GLUC3	LIP	LDHI2	TP2	TRIGL
	106.6	46.7	172	4.91	117.5
	UA2	UREAL	PHOS2	LDLC3	IRON2
	4.8	41.5	3.65	54.0	110.41
	CRP4	UIBCI	ALTL	HDLC4	
	5.93	205.0	47.9	27.0	

Data Monitor

19/09/23

09:00

C010026 095	LOT 53571900	PCCC2			
19/09/23	ALB2	ALP2L	AMYL2	ASTL	DBIL2
08:35:20	4.96	259	270	148.3	2.198
bmserv			ReagEX		
	BILT3	CA2	CHO2I	CREJ2	GGT12
	3.654	14.17	169.3	3.73	222
	GLUC3	LIP	LDHI2	TP2	TRIGL
	252.3	101.8	315	7.66	222.6
	UA2	UREAL	PHOS2	LDLC3	IRON2
	10.1	121.4	8.35	92.4	245.63
	CRP4	UIBCI	ALTL	HDLC4	
	53.72	254.3	127.6	55.0	

Data Monitor

19/09/23

09:11

C010027 095
19/09/23 Na
09:05:20 143
bmserv

LOT 53571900 PCCC2
K C1
7.31 106.0

Data Monitor

19/09/23

10:19

C010029 095
19/09/23 Na
10:13:28 142
bmserv

LOT 53571900 PCCC2
K C1
7.44 107.7