

GLUCOSE

METHOD : GOD / POD Method

PRINCIPLE : Glucose is Oxidised to Gluconic Acid and Hydrogen Peroxide in the presence of Glucose Oxidase. Hydrogen Peroxide further reacts with Phenol and 4 – Amino Antipyrine by the catalytic action of Peroxidase to form a Red Coloured Quinimine Dye Complex. Intensity of the Colour formed is directly proportional to the amount of Glucose present in the Sample.

NORMAL RANGE :

Fasting : 70 to 110 mgs/dl

Post Prandial: 80 to 140 mgs/dl

SAMPLE : Serum Or Sodium Fluoride Plasma.

MATERIAL REQUIRED : Glucose Reagent, Clean and Dry Test tubes, Micro Pipettes, Semi Auto Analyzer.

SYSTEM PARAMETERS :

Reaction : End Point Sample Volume : 0.01 ml

Wavelength : 505 nm Reagent Volume : 1.00 ml

Incubation Time : 37 ° C / 10 Min Standard Conc. : 100 mg/dl

Incubation Time : R. T. / 30 Min Linearity : 500 mgs/dl

PROCEDURE :

Bring the Glucose Reagent to room temperature.

Take Glucose Reagent 1.00 ml in a test tube.

Add 0.01 ml Serum / Plasma and Wait..

Take reading in Semi Auto Analyzer.

Prepared by : Lab Technician	Authorized by : Medical Officer