

# UNITY GASTRO CARE AND DIAGNOSTIC

DOC. NO.-UGCD/BIO/SOP/01

STANDARD OPERATING PROCEDURE FOR  
BIOCHEMISTRY

## PHOSPHOROUS

- a) **PURPOSE OF THE EXAMINATION:** Estimation of blood phosphorous level.  
b) **PRINCIPLE AND METHOD OF THE PROCEDURE USED FOR EXAMINATIONS:** Phosphomolybdate

### REACTION PRINCIPLE

Ammonium molybdate + Sulphuric acid + Phosphate  $\longrightarrow$  Inorg. Phosphorus molybdate complex  
The complex absorption is maximal at 340 nm.

### c) PERFORMANCE CHARACTERISTICS:

As per Kit Instruction. Measuring range 0.5 – 13 mg/dl.

d) **TYPE OF SAMPLE:** Serum.

e) **PATIENT PREPARATION:** As mentioned in QMSP 17.

### f) TYPE OF CONTAINER AND ADDITIVES:

Serum vial / vacutainer. No additives required.

### g) REQUIRED EQUIPMENT AND REAGENTS:

- ✓ TRANSASIA ERBA.
- ✓ Centrifuge.
- ✓ Sample (plasma).
- ✓ Micropipettes of variable volume from 0 – 1000  $\mu$ l.
- ✓ Isotonic Saline or Reagent grade water FS Diluent Pack 2 or FS Diluent Pack 3.

### h) ENVIRONMENTAL AND SAFETY CONTROLS:

### i) CALIBRATION PROCEDURES:

We perform Full calibration as per plan and when required with calibration kit (CAL KIT 1)

### PROCEDURAL STEPS:

- Blood is collected using standard laboratory procedures
- 3.0 ml of whole blood for serum preparation in vacutainer.
- Centrifuge whole blood at 3000 rpm for 15 minutes for serum separation.
- Put the vacutainer having clear serum into the sample rack of TRANSASIA ERBA

**STORAGE :** Separated serum should be stored in stoppered containers at -18 to -20 degree C for 1 week

**PROCEDURE FOR CONFIRMATION / CORRECTION:** 1. Dilute the sample with specialty diluent or isotonic range .

2. Reanalyze

3. Multiply the results by the dilution factor to obtain an estimate of the original sample amylase concentration.

Issue No. 01	Issue Date: 28.11.2023	Prepared By:	Copy No. 01	Page 1 of 2
Rev. No.: 00	Rev. Date: Nil	Approved By:	Issued By:	



# UNITY GASTRO CARE AND DIAGNOSTIC

DOC. NO.-UGCD/BIO/SOP/01

STANDARD OPERATING PROCEDURE FOR  
BIOCHEMISTRY

## ENVIRONMENTAL CONDITIONS:

1. Temp: 20-25°C.
2. Humidity: not more than 70%.
3. Dust free.

**k)QUALITY CONTROL PROCEDURES:** We run controls (level 1 & level 2) BIORAD every day,

**l)INTERFERENCES:** No significant interference up to 60 mg/dl of conjugated bilirubin and 40 mg/dl of unconjugated bilirubin.. And triglycerides concentration of 2000 mg/dl. Haemoglobin upto 1000.

**m)PRINCIPLE OF PROCEDURE FOR CALCULATING RESULTS INCLUDING, WE RELEVANT, THE MEASUREMENT UNCERTAINTY OF MEASURED QUANTITY VALUES:**

**n)BIOLOGICAL REFERENCE INTERVALS:** 2.5-4.5 mg/dl

**o)REPORTABLE INTERVAL OF EXAMINATION RESULTS:** Normally after 8 hours in case of emergency 3 hour .

**p)INSTRUCTIONS FOR DETERMINING QUANTITATIVE RESULTS WHEN A RESULT IS NOT WITHIN THE MEASUREMENT INTERVAL:** NA

**q)ALERT/CRITICAL VALUES: :**

**r)LABORATORY CLINICAL INTERPRETATION:** Not done

## s)POTENTIAL SOURCES OF VARIATION:

Increase in enzyme activity is observed in various hepatobiliary diseases and pancreatitis, acute myocardial infarction, heavy use of alcohol, carcinoma of breast and lung, neoplasms and carcinoma of prostate.

## t)REFERENCES:

- ✓ Tietz textbook of clinical chemistry. 4<sup>th</sup> ed. Philadelphia: W.B. Saunders Company; 1996.p.351-374.
- ✓ Kit Literature

Issue No. 01	Issue Date: 28.11.2023	Prepared By:	Copy No. 01	Page <b>2</b> of <b>2</b>
Rev. No.: 00	Rev. Date: Nil	Approved By:	Issued By:	