

# Urinalysis Reagent Strips (Urine) Package Insert

For rapid detection of multiple analytes in human urine.  
For *in vitro* diagnostic use only.

## INTENDED USE

The Urinalysis Reagent Strips (Urine) are five plastic strips onto which several separate reagent areas are affixed. The test is for the detection of one or more of the following analytes in urine: Ascorbic acid, Glucose, Bilirubin, Ketone (Acetoacetic acid), Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite and Leukocytes.

## SUMMARY

Urine undergoes many changes during states of disease or body dysfunction before being voided. Urinalysis is a significant part of a patient's health care. Urinalysis Reagent Strips (Urine) can be used in a part of routine health screening. The diagnosis and monitoring of metabolic or systemic diseases that affect kidney function, endocrine disorders and diseases or disorders of the urinary tract.<sup>1,2</sup>

## PRINCIPLE AND EXPECTED VALUES

**Ascorbic Acid:** This test involves decolorization of Timm's reagent. The presence of ascorbic acid causes the color of the test field to change from blue-green to orange.

**Glucose:** This test is based on the enzymatic reaction that occurs between glucose oxidase, peroxidase and chromogen. Glucose is first oxidized to produce gluconic acid and hydrogen peroxide in the presence of glucose oxidase. The hydrogen peroxide then reacts with potassium tetraoxygen in the presence of peroxidase. The event is green chromophore. Low amounts of glucose are normally excreted in urine.<sup>3</sup> Glucose excretion of less than 100 mg/dL is considered normal. Glucose excretion of 100 mg/dL or more is considered abnormal. For semi-quantitative results, read at 30 seconds only.

**Bilirubin:** This test is based on azo-coupling reaction of bilirubin with diazotized dihydroxime in a strongly acidic medium. Varying bilirubin levels will produce a pinkish-purple color for positive results. Ketones are normally not present in urine. Detectable levels may occur in urine during physiological stress conditions such as fasting, pregnancy and frequent strenuous exercise.<sup>4,5</sup> In starvation diets or in other abnormal carbohydrate metabolism situations, ketones appear in the urine in excessively high concentration before serum ketones are elevated.

**Specific Gravity:** This test is based on the apparent pKa change of certain pre-treated polyelectrolytes in relation to ionic concentration. In the presence of an indicator color in urine of increasing ionic concentration. Randomly collected urine may vary in specific gravity from 1.003-1.040. Twenty-four hour urine from healthy adults with normal diets and fluid intake will have a specific gravity of 1.015-1.022.<sup>6</sup> In cases of severe renal damage, the specific gravity is fixed at 1.010, the value of the glomerular filtrate.

**Blood:** This test is based on the peroxidase-like activity of hemoglobin which catalyzes the reaction of current-hydroperoxide and 3,3',5,5'-tetramethylbenzidine. The resulting color ranges from orange to green to dark blue. Any green spots or green color development on the reagent area within 60 seconds is significant and the urine specimen must be examined further. Blood is often, but not invariably, found in the urine of menstruating females.

**pH:** This test is based on a double indicator system which gives a broad range of colors covering the entire urinary pH range. Colors range from yellow and green to red/blue. The expected range for normal urine specimens is pH 4.5-8, with an average result of pH 6.

**Protein:** This reaction is based on the phenenon known as the "protein error" of pH indicators where an indicator that is highly buffered will change color in the presence of proteins (anions) as the indicator releases hydrogen ions to the protein. As a constant pH, the development of any green color is due to the presence of protein. Colors range from yellow to yellow-green for negative results and green to green-blue for positive results. 1-14 mg/dL of protein may be excreted by a normal kidney.<sup>7</sup> A color matching urinary block greater than this indicates significant proteinuria. For urine with high specific gravity, the test area may most closely match the trace color block even though only normal concentrations of protein are present. Clinical judgment is required to evaluate the significance of trace results.

**Urobilinogen:** This test is based on a modified Ehrlich reaction between 5-dihydroxyacetophenone and urobilinogen acid in strongly acidic medium to produce a pink color. Urobilinogen is one of the major compounds produced in heme synthesis and is a normal substance in urine. The expected range for normal urine with this test is 0.2-1.0 mg/dL (3.5-17 µmol/L). A result of 2.0 mg/dL (35 µmol/L) may be of clinical significance, and the patient specimen should be further evaluated.

**Nitrite:** This test depends upon the conversion of nitrate to nitrite by the action of Gram negative bacteria in the urine. In an acidic medium, nitrite in the urine reacts with p-aminobenzenesulfonamide to produce a pink color. Nitrite in the urine is dependent on the normal urinary flora. Urine specimens with nitrite are more likely to be contaminated with bacteria. Urine specimens with nitrite test ranges from as low as 45% in cases where

bladder inoculation occurred, to as high as approximately 80% in cases where bladder inoculation took place for at least 4 hours.

**Leukocytes:** This test reveals the presence of granulocyte esterase. The esterase cleave a derivatized pyrazole amino acid ester to liberate derivatized hydroxy pyrazole. This pyrazole then reacts with a diazonium salt to produce a large pink to purple color. Normal urine specimens generally yield negative results. Trace results may be of questionable clinical significance. When trace results occur, it is recommended to retest using a fresh specimen from the same patient. Repeated trace and positive results are of clinical significance.

## REAGENTS AND PERFORMANCE CHARACTERISTICS

Based on the dry weight at the time of impregnation, the concentrations given may vary with manufacturing tolerances. The following table below indicates read times and performance characteristics for each parameter:

Reagent	Read Time	Composition	Description
Ascorbic Acid (ASC)	30 seconds	0.3% w/w 2,6-dichlorophenolindophenol, 99.7% w/w buffer and non-reactive ingredients	Detects ascorbic acid as low as 0.25-0.50 mg/dL (0.25-0.50 mmol/L)
Glucose (GLU)	30 seconds	0.3% w/w glucose oxidase, 10.0% w/w potassium tetraoxygen, 13.0% w/w non-reactive ingredients	Detects glucose as low as 50-100 mg/dL (2.5-5.0 mmol/L). Results may be read at 10 seconds for qualitative semi-quantitative results or at 30 seconds for semi-quantitative results.
Bilirubin (BIL)	30 seconds	2,4-dichloroaniline, 0.5% w/w sodium sulfite, 99.5% w/w non-reactive ingredients	Detects bilirubin as low as 0.4-0.8 mg/dL (6.8-13.6 µmol/L)
Ketone (KET)	40 seconds	5% w/w sodium nitroprusside, 95% w/w buffer	Detects acetoacetic acid as low as 2.5-5.0 mg/dL (0.5 mmol/L)
Specific Gravity (SG)	45 seconds	17.5% w/w buffer and non-reactive ingredients; 55% w/w polyelectrolyte (polyvinylidene sulfonate sodium salt)	Determines urine specific gravity. Between 1.000 and 1.030. Results correlate with values obtained by refractometer (within 0.005).
Blood (BLO)	60 seconds	4% w/w hemolyzed diaminotetrahydroxy-TMB, 6% w/w uricemic, 90% w/w buffer and non-reactive ingredients	Detects free hemoglobin as low as 0.015-0.062 mg/dL or 5-10 EryUL in urine. Specimens with ascorbic acid will give a false negative result (within 450 mg/dL sodium salt).
pH	60 seconds	Sodium salt, 94.5% w/w non-reactive ingredients	Permits the quantitative determination of pH values within the range of 5-10.
Protein (PRO)	60 seconds	0.3% w/w tetrabromophenol blue, 99.7% w/w buffer and non-reactive ingredients	Detects albumin as low as 30 mg/dL (0.3 g/L) and non-albumin protein as low as 0.2-1.0 mg/dL (0.2-10 µmol/L).
Urobilinogen (URO)	60 seconds	2.5% w/w p-aminobenzenesulfonamide, 91.5% w/w buffer and non-reactive ingredients	Detects urobilinogen as low as 0.2-1.0 mg/dL (3.5-17 µmol/L).
Nitrite (NIT)	60 seconds	95.5% w/w non-reactive acid, 4.5% w/w p-aminobenzenesulfonamide	Detects urinary nitrite as low as 0.05-0.1 mg/dL in urine with a low specific gravity and less than 30 mg/dL ascorbic acid.
Leukocytes (LEU)	120 seconds	0.3% w/w derivatized pyrazole amino acid ester, 0.4% w/w diazonium salt, 92.5% w/w buffer, 0.1% w/w non-reactive ingredients	Detects leukocytes as low as 10-25 white blood cells per high power field in clinical urine.

The performance characteristics of the Urinalysis Reagent Strips (Urine) have been determined in both laboratory and clinical tests. Parameters of importance to the user are sensitivity, specificity, accuracy and precision. Generally, this test has been developed to be specific for the parameters to be measured with the exception of the interferences listed. Please refer to the Limitations section in this package insert.

Interpretation of visual results is dependent on several factors, the variability of color perception, the presence or absence of inhibitory factors, and the lighting conditions when the strip is read. Each color block on the chart corresponds to a range of analyte concentrations.

## PRECAUTIONS

- For *in vitro* diagnostic use only. Do not use after the expiration date.
- Do not use if the cap is damaged or the cap is missing.
- Do not touch the reagent areas of the strip.
- Do not use any damaged strips that may have deteriorated.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used strip should be discarded according to local regulations after testing.

## STORAGE AND STABILITY

Strips are packaged in the closed cap and stored at room temperature or refrigerated (2-30°C). Keep out of direct sunlight. The strip is stable through the expiration date when the cap is closed. Do not remove the cap until ready to use. Do not use strips that have been opened for more than 30 days. Do not use strips that have been opened for more than 30 days. Do not use strips that have been opened for more than 30 days. Do not use strips that have been opened for more than 30 days.

## SPECIMEN COLLECTION AND PREPARATION

A urine specimen must be collected in a clean and dry container and used as soon as possible. Do not centrifuge. The use of urine preservatives is not recommended. If preservatives are used, the urine should be refrigerated and the specimen immediately analyzed. Do not use urine that has been refrigerated for more than 24 hours. Do not use urine that has been refrigerated for more than 24 hours. Do not use urine that has been refrigerated for more than 24 hours.

## MATERIALS

- Specimen container
- Reagent Strips
- Materials Required But Not Provided
- Specimen collection container
- Timer

## DIRECTIONS FOR USE

- Remove the top, urine specimen, and/or control to reach room temperature (15-30°C) prior to testing.
- Remove the cap from the closed cap and use it as soon as possible. Do not use the cap for more than one specimen. The cap should be discarded after use. Do not use the cap for more than one specimen. The cap should be discarded after use. Do not use the cap for more than one specimen. The cap should be discarded after use.
- Remove the cap from the closed cap and use it as soon as possible. Do not use the cap for more than one specimen. The cap should be discarded after use. Do not use the cap for more than one specimen. The cap should be discarded after use.
- Remove the cap from the closed cap and use it as soon as possible. Do not use the cap for more than one specimen. The cap should be discarded after use. Do not use the cap for more than one specimen. The cap should be discarded after use.
- Remove the cap from the closed cap and use it as soon as possible. Do not use the cap for more than one specimen. The cap should be discarded after use. Do not use the cap for more than one specimen. The cap should be discarded after use.

## INTERPRETATION OF RESULTS

Results are obtained by direct comparison of the color blocks printed on the device to the color blocks on the chart. The color blocks represent normal values. Results may vary due to the normal values. In the event of unexplained or questionable results, the following tests are recommended: confirm that the specimen has been tested within the expiration date; retest on the control strip; compare the results of the problem patient specimen with the control strip; repeat the test using a fresh specimen; if the problem persists, describe any other test results and/or control panel test results.

## QUALITY CONTROL

For best results, performance of reagent strips should be confirmed by using known positive and negative specimens. Each secondary should include a new test, a control, and a reagent strip. Each secondary should include a new test, a control, and a reagent strip. Each secondary should include a new test, a control, and a reagent strip. Each secondary should include a new test, a control, and a reagent strip.

## LIMITATIONS

Note: As with all diagnostic and therapeutic tests, all results must be considered with other clinical information available to the physician.

**Ascorbic Acid:** This test is highly specific for glucose. No substance except in urine other than glucose is known to give a positive result. The reagent area does not react with ketones, lactates, glycerol, fructose or other metabolic substances. It may be depressed in specimens with high specific gravity (1.025) and with ascorbic acid concentrations of 210 mg/dL.

## ESBIOGRAPHY

1. *Am J Physiol* 1961; 201: 1-12.

2. *Am J Physiol* 1961; 201: 1-12.

3. *Am J Physiol* 1961; 201: 1-12.

4. *Am J Physiol* 1961; 201: 1-12.

5. *Am J Physiol* 1961; 201: 1-12.

6. *Am J Physiol* 1961; 201: 1-12.

7. *Am J Physiol* 1961; 201: 1-12.

8. *Am J Physiol* 1961; 201: 1-12.

9. *Am J Physiol* 1961; 201: 1-12.

10. *Am J Physiol* 1961; 201: 1-12.

## REFERENCES

1. *Am J Physiol* 1961; 201: 1-12.
2. *Am J Physiol* 1961; 201: 1-12.
3. *Am J Physiol* 1961; 201: 1-12.
4. *Am J Physiol* 1961; 201: 1-12.
5. *Am J Physiol* 1961; 201: 1-12.
6. *Am J Physiol* 1961; 201: 1-12.
7. *Am J Physiol* 1961; 201: 1-12.
8. *Am J Physiol* 1961; 201: 1-12.
9. *Am J Physiol* 1961; 201: 1-12.
10. *Am J Physiol* 1961; 201: 1-12.

## CONTACT INFORMATION

For more information, contact the following:

1. *Am J Physiol* 1961; 201: 1-12.

2. *Am J Physiol* 1961; 201: 1-12.

3. *Am J Physiol* 1961; 201: 1-12.

4. *Am J Physiol* 1961; 201: 1-12.

5. *Am J Physiol* 1961; 201: 1-12.

6. *Am J Physiol* 1961; 201: 1-12.

7. *Am J Physiol* 1961; 201: 1-12.

8. *Am J Physiol* 1961; 201: 1-12.

9. *Am J Physiol* 1961; 201: 1-12.

10. *Am J Physiol* 1961; 201: 1-12.

## NOTES

1. *Am J Physiol* 1961; 201: 1-12.

2. *Am J Physiol* 1961; 201: 1-12.

3. *Am J Physiol* 1961; 201: 1-12.

4. *Am J Physiol* 1961; 201: 1-12.

5. *Am J Physiol* 1961; 201: 1-12.


6. *Am J Physiol* 1961; 201: 1-12.

7. *Am J Physiol* 1961; 201: 1-12.

8. *Am J Physiol* 1961; 201: 1-12.

9. *Am J Physiol* 1961; 201: 1-12.

10. *Am J Physiol* 1961; 201: 1-12.


**RPO DIAGNOSTIC REAGENT CO. INC.**  
 2300 RIVINGTON STREET, NEW YORK, N.Y. 10021  
 (212) 261-1000  
 FAX: (212) 261-1001  
 Telex: 270 200 RPO  
 Cable: RPO  
 © 1983 RPO Diagnostic Reagent Co. Inc.



ACON Biotech (Hangzhou) Co., Ltd.  
 No.210 Zhenzhong Road, West Lake District  
 Hangzhou, P.R. China, 310030

ACON

Mission

Urinalysis Reagent Strips

REF U031-101

10U

IND

ISO 13485

Σ 100

LOT URS3080061

2025-10-09

CE

LEU 120 s	-	15 ±	75 +
NIT 60 s	-	+	pk nsa nsa/s
URO 60 s	0.2(3.5)	1(17)	25
PRO 60 s	-	15(0.15) ±	20-30
pH 60 s	5.0	6.0	6.5
BLO 60 s	-	±	2+
SG 45 s	1.000	1.025	1.050
KET 40 s	-	5(0.5) ±	50-80
BIL 30 s	-	1(17) +	20-30
GLU 30 s	-	10(0.5) ±	20-30

