



results of the strip test and those from the refractive index method is only 0.005. To make it more accurate, 0.006 may be added to readings from urines with pH equal or greater than 5. Urine reading instrument can automatically make these adjustments in strip-readings. The urine reagent constituents such as glucose or radiopaque dye won't make any changes in the test. Highly buffered alkaline urines may cause the low readings comparing with the other methods. Elevated specific gravity readings may occur in the presence of moderate quantities of protein (9.7-7.5g/L).

#### Blood

Trace reaction may vary among the patients. Clinical judgments are required for individual cases. The presence of green spots (band erythrocytes) or green colour (haemoglobin/myoglobin) on the reagent area within 60 seconds indicates for further diagnostic check. Blood is often found in the urine of the menstruating females. Haemoglobin 150 µg/L-520 µg/L is approximately equivalent to 5-15 cells/L intact erythrocytes.

The reagent strip is highly sensitive to haemoglobin and this can be used as a supplementary to the microscopic examination. The sensitivity of the strip might be reduced in urine with a large amount of specific gravity. The strips are equally sensitive to myoglobin as to haemoglobin. Certain oxidizing contaminants, such as hypochlorite, may lead to false positive results. Microbial peroxidase associated with urinary tract infection may also produce a false positive result. Ascorbic acid less than 5.0mmol/L in urine may not influence the result of the test.

#### pH

The strip tests for pH values are generally in the range of 5.0-8.5 visually and 5.0-9.0 instrumentally.

#### Protein

The reagent area is more sensitive to albumin than to globulins, haemoglobin, Bence-Jones protein and mucoprotein. So a 'Negative Result' is not good enough to indicate that these proteins don't exist in urine. Normally no protein is detectable in urine with conventional methods, although a minute amount of protein is excreted through a normal kidney; though a minute amount when the colour is darker than mark on the chart. False positive results may be obtained in highly buffered alkaline urines. Urine specimens contaminated with quaternary ammonium compounds and cleansers containing phenothiazine may also produce false positive results.

#### Urobilinogen

The reagent strips can detect urobilinogen in low amount as 3 µmol/L (approximately 0.2 Ehrlich unit/L) in urine. A result of 33 µmol/L in urine indicates the critical value, representing the transition from normal to abnormal, which requires further check on patients and specimens. The negative results are not final to determine the absence of urobilinogen.

#### Nitrite

Gram-negative bacteria in urine converts nitrate, derived from food/drink in urine. The reagent strip is essential to nitrite and won't react with the other substances in urine. Pink spots or edges on the strip should not be interpreted as positive result, but any degrees of uniform pink colour development should be taken as positive result. The degree of colour development (the numbers of bacteria are not in direct proportion). The negative result doesn't mean the absence of bacteria in a large amount. Negative result may occur (1) when urine doesn't contain organism that caused the conversion from nitrate to nitrite. (2) when urine has not remained in the bladder long enough (four hours up to let the nitrate convert into nitrite) (3) the nitrate in the foods is absent. Large High volume of specific gravity in urine may reduce the sensitivity of the test. 1.4mmol/L ascorbic acid or less won't interfere the test result.

#### LEUKOCYTES

Test area read with esterase in leucocytes (granulocytic leucocytes). Normal urine specimens generally yield negative result, positive result (+ or greater) are clinically significant. Individually observed 'Trace' results may be of questionable clinical significance, however 'Trace' results observed repeatedly may be clinically significant. 'Positive' results may occasionally be found with random specimens from females due to contamination of the specimen by vaginal discharge. Elevated glucose concentrations (160mg/dL) or high specific gravity may cause decreased test results.

#### SPECIFIC PERFORMANCE CHARACTERISTICS

Specific performance characteristics are based on clinical and analytical studies. In clinical specimens, the sensitivity depends upon several factors, the variability of colour perception, the presence or absence of inhibitory factors typically found in urine, specific gravity, pH, and the lighting conditions when the product is read visually. Each colour block or instrumental display value represents a range of values. Because of specimen and reading variability, specimens with analyte concentrations that fall between nominal levels may give results at either level. Results at levels greater than the

second positive level for the Protein, Glucose, Ketone, and Urobilinogen tests will usually be within one level of the true concentration. Equal agreement between visual results and instrumental results might not be found because of the inherent differences between the perception of the human eye and the optical system of the instrument.

#### Notes

This product is used for in vitro diagnosis. The strips must be kept in the original bottle. Never use the products after the expiration date. Each strip can be used only once. Do not remove the desiccating. If strips are removed and tightly after taking out the strips. The strips should be stored in a dry place at the temperature between 2°C-30°C. Do not store the strips in refrigerator and keep them away from direct sunlight. Do not touch the reagent area of the strip. Protection against ambient moisture, light and heat is essential to guard against altered reagent reactivity. Deterioration may result in discoloration or darkening of the reagent area of the strip. If this happens, or the test results are questionable or inconsistent with the expected results, check and make sure the strips are within the expiration date and also run a control. Please dispose the reagent strips as waste according to Treatment Regulation of Lab Biological Materials.

#### SENSITIVITY AND TEST RANGE OF URINALYSIS STRIPS

Item	Sensitivity	Instrumental test range	Visual test range
Glucose (mmol/L)	2.8-5.5	Neg-55	Neg-110
Protein (g/L)	0.15-0.3	Neg-3.0	Neg-20.0
Ketone (acetone acid) (mmol/L)	0.5-1.0	Neg-7.8	Neg-16
Blood (cells/L)	5-15	Neg-200	
Bilirubin (µmol/L)	3.3-17	Neg-Trace	
Nitrite (µmol/L)	13-22	Neg-Trace	
Leukocytes (cells/L)	5-15	Neg-500	
Urobilinogen (µmol/L)	3.3-16	3.3-131	3.3-218
pH	5.0-9.0	5.0-9.0	5.0-8.5
Specific gravity		1.005-1.200	1.000-1.010

#### REACTIVE INGREDIENTS (based on dry weight at time of impregnation)

Protein	ethylmercaptan blue butter sensitive blue pigments azobenzene sensitive blue pigments sensitive blue pigments	0.1% w/w 97.48% w/w 2.5% w/w 26.0% w/w 1.2% w/w 1.2% w/w 37.2% w/w
Blood	diisopropylbenzene alkyon peroxide butyl sensitive blue pigments sensitive blue pigments	1.7% w/w 0.2% w/w 11.8% w/w 26.2% w/w
Glucose	glucose oxidase (Glucooxidase 12311) peroxidase butter sensitive blue pigments	0.2% w/w 11.8% w/w 26.2% w/w
Ketone	sodium nitroprusside butter sensitive blue pigments	5.1% w/w 14.4% w/w 2.7% w/w
Leukocytes	pyrrole amino acid ester diagonium salt butter sensitive blue pigments	4.1% w/w 0.4% w/w 1.2% w/w 2.7% w/w
Nitrite	p-aminosalicylic acid (N-Nitrosation) 4-dihydroxyquinone butter sensitive blue pigments	1.3% w/w 60.6% w/w 8.2% w/w 4.8% w/w
Specific Gravity	benzoin butter sensitive blue pigments sensitive blue pigments sensitive blue pigments	95.2% w/w 5.0% w/w 3.5% w/w 3.5% w/w 3.5% w/w
pH	amylal and blue decolorative ingredients	41.1% w/w 41.1% w/w
Bilirubin	2,4-dichlorophenyl dimethyl amine salt butter decolorative ingredients	6.6% w/w 44.1% w/w 47.1% w/w
Urobilinogen	p-toluidino benzaldehyde butter decolorative ingredients	42.1% w/w 42.1% w/w 1.8% w/w

#### Notes on symbols and marks

LOT	Batch code	Use by
Single use	Single use	In Vitro Diagnostic Use
Manufactured by	CE	Shown at
Please read package insert	CE	These labels conform to the
Adapted Representative	CE	directive 90/269/EEC
	CE	Catalogue number

**DIRUI**

DIRUI Industrial Co., Ltd  
95 Yunhe Street, New & High Tech,  
Development Zone  
Changshun, Jilin 130012 P.R. China

Tel: +86(431)85100409  
Fax: +86(431)85172581  
E-mail: dirui@dirui.com.cn  
Http://www.dirui.com.cn