



# YUMIZEN H550

(Serial no: 909YAXH02625)

Fully Automated Hematology Analyzer

**Performance Qualification** 

For

# LUPIN HEALTHCARE LTD, LATUR

#246, Okhla Industrial Estate, Phase III, New Delhi 110020, India, Tel: 011 4646 5000.

Visit us: http://www.horiba.com/in/



### A. Performance Qualification

#### A. Instrument Identification:

**Instrument Name** 

: YUMIZEN H550

Serial Number

: 909YAXH02625

# B. Following is the list of test to be performed and verified

Blank Reference cycle: To verify the Startup Cycle of the instrument.

Serial No: 909YAXH02625

Parameters	Acceptable Range	Observed Value
WBC 10 <sup>3</sup> /mm <sup>3</sup>	$\leq 0.3 \times 10^3 / \text{ mm}^3$	0.07
RBC 10 <sup>6</sup> /mm <sup>3</sup>	$\leq 0.02 \times 10^6 / \text{mm}^3$	0.00
HGB g/Dl	≤ 0.3 g/dl	0.0
PLT 10 <sup>3</sup> /mm <sup>3</sup>	≤ 5×10 <sup>3</sup> /mm <sup>3</sup>	0

Conducted By:



Precision Study: Precision is checked by running blood sample in 10 replicates & getting CV% in within acceptance.

Serial No: 909YAXH02625

Parameters	CV % Acceptance	CV % Observed	Comments
RBC 10 <sup>6</sup> /mm <sup>3</sup>	< 2.0	0.78	Ok
HGB g/dL	< 1.5	0.37	Ok
HCT %	< 2.0	0.88	Ok
PLT 10 <sup>3</sup> /mm <sup>3</sup>	< 5.0	2.32	Ok
WBC 10 <sup>3</sup> /mm <sup>3</sup>	< 2.5	2.14	Ok

Conducted By:



<u>Calibration</u>: To calibrate the Instrument using calibrator (ABX Minocal) and verify the same.

**Procedure:** Go to Quality Assurance icon on main screen and then Calibration icon. Run Calibrator (ABX Minocal) 10 times, without taking the values of first run, calibrate the instrument using average of the last 10 runs. Run Calibrator as a sample to verify the calibration. **Lot:** CX471; Expiry: 05/08/2022.

Serial No: 909YAXH02625

Parameter	Target Value (As per Kit Insert)	Mean Value	Observed CV%	Acceptance CV%	Comments
WBC	8.90	8.82	1.12	<2%	PASSED
RBC	4.57	3.98	0.44	<2%	PASSED
HGB	13.2	13.2	0.37	<1.5%	PASSED
НСТ	38.6	33.7	0.48	<2%	PASSED
PLT	256	215	3.04	<5%	PASSED
MPV	10.9	10.8	1.73	<3%	PASSED

Conducted By:



- <u>Control Runs:</u> The quality of the analyzer is checked by running three levels of Controls & getting the values in the range as per the kit insert.
  - Lot: ABX Difftrol PX436; Exp: 05/09/2022.

Serial No: 909YAXH02625

STATUS
PASSED
PASSED
PASSED

Conducted By:



- <u>Carryover Study:</u> Carry over is checked by running quality controls /samples (high & low) in 3 replicates & getting CV% in within acceptance.
  - Carry Over %= (L1-L3)\*100/(H3-L3).

Serial No: 909YAXH02625

Parameters	WBC 10 <sup>3</sup> /mm <sup>3</sup>	RBC 10 <sup>6</sup> /mm <sup>3</sup>	HGB g/dL	нст %	PLT 10 <sup>3</sup> /mm <sup>3</sup>
Carry Over (%)	0.69	0.97	0	0.91	0.89
Manufacturer acceptable CV%	<1%	<1%	<1%	<1%	<1%
Status	Passed	Passed	Passed	Passed	Passed

Conducted By:



## **B. PERFORMANCE CERTIFICATE:**

**Instrument Name** 

: YUMIZEN H550

Serial Number

: 909YAXH02625

**Customer Details** 

: LUPIN HEALTHCARE, LATUR

with complete address

Barmade Hospital And Life Advanced Test Tube Baby

Centre, Opposite Datta Mandir, Adarsh Colony, Ausa Road,

Latur, Maharashtra 413531.

**Installation Date** 

:14/07/2022

Warranty expires on

: RR (Reagent Rental) - 5yrs

Prepared by:	HORIBA Medical - HORIBA India Pvt. Ltd.			
Name:	Yogesh Gawade			
Title: Application S	pecialist	Sign:	Date:14/07/2022	
Approved by:	LUPIN HEALTHCARE, LATUR			
Name:	Mr. Mustkim Shaikh			
Title: Head of Department		Sign:	Date: 14/07/2022	

**Deviation:** No

**Conclusion:** Instrument has been qualified for Performance.



# **ATTACHMENTS**

Attachment No	<b>Document Title</b>	Page Nos.	Verified By Sign/Date
1	Blank Cycle	1	
2	Repeatability	2	
3	Calibration	4	
4	Calibrator as a sample	6	
5	<b>Quality Control</b>	7	
6	Carry over	10	