



CALIBRATION REPORT

Customer Name : BM Lab and Digital X-ray
Equipment Name: EM 200 Biochemistry Analyzer
Make :ERBA
Serial Number : NA
Date of Calibration :11-03-2024

Asset ID :BM/BA/01
Calibration Due date :11-03-2025

Item to be checked/replace		Check	Comments
DI water Filter	Replace .	<input checked="" type="checkbox"/>	
Tubing	Replace Laundry dispense tubings	<input checked="" type="checkbox"/>	
	Replace the Laundry aspiration tubings.	<input checked="" type="checkbox"/>	
	Replace the detergent tube	<input checked="" type="checkbox"/>	
	Clean the waste tubings.	<input checked="" type="checkbox"/>	
Arm assemblies	Clean probes internally/externally	<input checked="" type="checkbox"/>	
	Check the centering	<input checked="" type="checkbox"/>	
	Clean troughs / drains/opto coupler assemblies.	<input checked="" type="checkbox"/>	
Tray assemblies	Clean optointerruptor& coupler assemblies	<input checked="" type="checkbox"/>	
	Clean the container assemblies.	<input checked="" type="checkbox"/>	
Detergent Filter	Replace filter		
Lamp	Replace	<input type="checkbox"/>	
Fuses	Check replace if blown	<input checked="" type="checkbox"/>	Checked, good
Laundry Assembly	Clean probes	<input checked="" type="checkbox"/>	
	Clean the optos	<input checked="" type="checkbox"/>	
	Check opto coupler connections	<input checked="" type="checkbox"/>	
	Clean the dust on shafts	<input type="checkbox"/>	
Stirrer assembly	Clean stirrer paddles externally	<input checked="" type="checkbox"/>	
	Clean stirrer troughs & Drains	<input checked="" type="checkbox"/>	
	Clean the opto coupler assembly.	<input checked="" type="checkbox"/>	
	Clean the dust on Shaft	<input checked="" type="checkbox"/>	
	Check the electrical motor connections.	<input checked="" type="checkbox"/>	
Syringe assemblies	Clean the syringe path	<input checked="" type="checkbox"/>	
	Clean the opto assembly.	<input checked="" type="checkbox"/>	
Sampler Rack assembly(XL 1000)	Clean the transport assembly	<input checked="" type="checkbox"/>	
	Check the motor connections	<input checked="" type="checkbox"/>	
	Clean opto coupler assemblies.	<input checked="" type="checkbox"/>	
RCT cuvette ring set	Clean cuvettes.	<input checked="" type="checkbox"/>	
	Check the heater/sensor connections.	<input checked="" type="checkbox"/>	
	Clean the opto coupler assembly.	<input checked="" type="checkbox"/>	
Reagent cooling	Clean the heat sink	<input checked="" type="checkbox"/>	





Item to be checked/replace		Check	Comments
Reagent assembly	cooling		
	Clean the heat sink	<input checked="" type="checkbox"/>	
	Clean the Fans	<input checked="" type="checkbox"/>	
Purging	Perform Cal A & Cal B purging	<input checked="" type="checkbox"/>	
Calibration	Call	<input checked="" type="checkbox"/>	Calibration done AST and ALT for Precision and Linearity verification
	brate ISE	<input checked="" type="checkbox"/>	
Probe wash	Perform probe wash	<input checked="" type="checkbox"/>	

Item to be checked/replace		Check	Comments
Auto wash	Perform Autowash	<input checked="" type="checkbox"/>	
Software Adjustment	Check the stirrer centering	<input checked="" type="checkbox"/>	
	Check the R1 probe centering	<input checked="" type="checkbox"/>	
	Check the R2 probe centering	<input checked="" type="checkbox"/>	
	Check the laundry centering	<input checked="" type="checkbox"/>	
R1/R2/ SPT/probes/troughs	Check the flow of DI water from probes . Clean path if not proper.	<input checked="" type="checkbox"/>	
	Check the flow of DI water from troughs adjust if required.	<input type="checkbox"/>	
CPU	Execute scan disk if required	<input checked="" type="checkbox"/>	
	Defragmentation HDD if required	<input type="checkbox"/>	not required
	Check program version	<input checked="" type="checkbox"/>	
Data verification after PM	Check the 1 bar pressure	<input checked="" type="checkbox"/>	
	Check Rct temperature(37°C)	<input checked="" type="checkbox"/>	
	Check control material data or sample data	<input type="checkbox"/>	
	Check RGT temperature(below 12°C)	<input checked="" type="checkbox"/>	
	Check filter gains/cell blank	<input checked="" type="checkbox"/>	
Exterior Instrument	Clean instrument exterior panels.	<input checked="" type="checkbox"/>	
	Clean peripheral equipment (Printer monitor trolley etc)	<input checked="" type="checkbox"/>	
Operation Review	Operation testing	<input checked="" type="checkbox"/>	

- Backup done
- Test Laniary Standardization Checked - Passed
- Checked all the settings
- Linearity and SD is checked as per calibration of the Equipment

Engineer Name: Pradeep Kumar
Signature:

Customer Name: Jithendra
Sign & Stamp :



Calibration Raw Data

BM Laboratory

EM 200 calibration History

Test: SGOTD

Calibration date: 11 March -2024 12 :23:20 Acceptable limit(%) :15

Pos	Consumable	Lot No	Conc	Delta Abs /Min	Factor
S1	BLANK	1234	0.000	-0.0009	0.0000
S2	MULTI CAL	2310072	96.400	-0.0381	-2587.9197

Calibration date:11 March -2024 12 :23:20

R1 Lot No :S022371-309

R2 Lot No : S022371-309

Indication: Successful Calibration



Calibration Raw data

BM Laboratory

EM 200 calibration History

Test: SGPTD

Calibration date: 11 March -2024 12 :22:21 Acceptable limit(%) :15

Pos	Consumable	Lot No	Conc	Delta Abs /Min	Factor
S1	BLANK	1234	0.000	-0.0009	0.0000
S2	MULTI CAL	2310072	91.300	-0.0326	-2871.0691

Calibration date:11 March -2024 12 :22:21


R1 Lot No :S032310-311

R2 Lot No : S032310-311

Indication: Successful Calibration



Calibration Raw Data (SGOT)



Calibrator Multical Kit Insert

XL MULTICAL



Lot No.: 2310072

Expiry: 2025-04



ASSAYED VALUES USING ERBA REAGENTS
HODNOTY PRO SOUPRAVY ERBA

Abstr. Zkratka	Parameter Analyt	Methodology Metoda	Kit Souprava	Unit Jednotka	Value Hodnota	Unit Jednotka	Value Hodnota
ALB	Albumin	BCG	ALB 501 ALB 250504	g/l	3.27	g/l	33.7
ALP	Alkaline phosphatase Alkalická fosfatáza	IFCC with AMP buffer IFCC s AMP pufrém	ALP 113 ALP AMP 102505	U/l	320	µkat/l	5.34
		GGCK with MEG buffer GGCK s MEG pufrém	ALP 102 500	U/l	368	µkat/l	6.14
		IFCC without PDP IFCC bez PDP	ALTOPT 330 ALTOPT 250505	U/l	91.0	µkat/l	1.52
ALTOPT	Alanine aminotransferase Alaninaminotransferáza	IFCC with PDP IFCC s PDP	ALTOPT 330 ALTOPT 102505 + PDP 50	U/l	91.3	µkat/l	1.52
AMY	Amylase Amyláza	CMPO3	AMY 110 AMY 5000 5 100	U/l	187	µkat/l	3.12
		IFCC - maltoheptasid	AMY 200500	U/l	171	µkat/l	2.86
AMY-P	Pancreatic amylase Pancreatická Amyláza	EPS - maltoheptasid	PAMY 400	U/l	171	µkat/l	2.86
AST/GOPT	Aspartate aminotransferase Aspartátaminotransferáza	IFCC without PDP IFCC bez PDP	AST/GOPT 100 AST/GOPT 400000	U/l	96.4	µkat/l	1.61
		IFCC with PDP IFCC s PDP	AST/GOPT 330 AST/GOPT 250505 + PDP 50	U/l	56.4	µkat/l	1.61
BL-D	Bilirubin Direct Bilirubin přímý	Diazo with sulphuric acid Diazo s kyselou sulfurovou	BL D 200500 BL D 30 300	mg/dl	2.32	µmol/l	39.6
		Diazo with dichloroaniline (DCA) Diazo s dichloroanilinem (DCA)	BL D DCA 320 BL D DCA 500	mg/dl	2.38	µmol/l	40.9
BL-T	Bilirubin Total Bilirubin celkový	Diazo with sulphuric acid Diazo s kyselou sulfurovou	BL T 300500 BL T 30 300	mg/dl	4.83	µmol/l	82.5
		Diazo with dichloroaniline (DCA) Diazo s dichloroanilinem (DCA)	BL T DCA 330 BL T DCA 500	mg/dl	4.49	µmol/l	76.8
CA	Calcium Vápník	Arenasso III	CA 100 CA 100050	mg/dl	10.6	mmol/l	2.64
CHE	Cholinesterase Cholinesteráza	Butyrylthiocholine Butyrylthiocholin	CHE 130	U/l	6160	µkat/l	84.2
CHOL	Cholesterol	CHOD-PAP	CHOL 440 CHOL 1000500100	mg/dl	146	mmol/l	3.77
CK	Creatine kinase NAC Kreatinkináza NAC	IFCC	CK 110 CK 100050	U/l	373	µkat/l	6.23
CKMB	Creatine kinase MB Kreatinkináza MB	Immunochemical / IFCC Imunochémická / IFCC	CK MB 110 CK MB 100	U/l	426	µkat/l	7.11
CL	Chloride Chlór	Mercuric Thiocyanate Thiokyanátová metoda	CL 130	mmol/l	101	mmol/l	101
CREA	Creatinine Kreatinin	Jaffe's kinetic Jaffe kinetický	CREA 210 CREA 500	mg/dl	3.96	µmol/l	362
		Enzymatic Enzymatický	CREA 410 370 CREA 210 284	mg/dl	3.85	µmol/l	341
GGT	Gamma-glutamyl transferase Gamma-glutamyltransferáza	Gluta-C, standardized against IFCC Gluta-C, standardizováno vůči IFCC	GGT 110 GGT 100050	U/l	104	µkat/l	1.74
GLU	Glucose Glukóza	GOD-POD GOD-POD	GLU 410 GLU 1000500100	mg/dl	200	mmol/l	11.1
		Hockinase Hexokináza	GLU 100 200	mg/dl	200	mmol/l	11.1
LDH	Lactate dehydrogenase Laktátdehydrogenáza	DEKIC	LDH 110 LDH 100	U/l	500	µkat/l	8.34
		IFCC	LDH 100050	U/l	263	µkat/l	4.39
LIP	Lipase Lipáza	Enzymatic colorimetric test Enzymatické kolorimetrické stanovení	LIP 110	U/l	92.9	µkat/l	1.58
MG	Magnesium Hořčík	Xylyl Blue Metoda s xylylovou modří	MG 30 MG 200	mg/dl	2.72	mmol/l	1.12
PHOS	Phosphate Fosfor	Ammonium Molybdate UV Molybdat UV	PHOS 100 PHOS 210	mg/dl	5.62	mmol/l	1.78
TG	Triglycerids Triglyceridy	GPO-Tinder	TG 40 TG 1000500100	mg/dl	138	mmol/l	1.56
TP	Total Protein Celková bílkovina	Buret	TP 40 TP 200050	g/l	4.88	g/l	48.8
UA	Uric Acid Kyselina močová	Uricase-Tinder	UA 375440	mg/dl	5.48	µmol/l	325
		Uricase-colorimetric with AOD Uricase-kolorimetrická s AOD	UA 500	mg/dl	5.48	µmol/l	325
UREA	Urea Močovina	Uricase-GLDH, kinetic Uricase-GLDH, kinetický	UREA 175 UREA 200050	mg/dl	101	mmol/l	16.8



12000000
20230000



Erba Lachema s.p., Katschek 2219/Id, 621 09 Brno, CZ
e-mail: @prognosa@erba.com, www.erbatestprogn.com

A-234010001

Date of manufacturing order: 10.03.2023

Scanned with CamScanner





KOMMUNJE HEALTHCARE

Mob : 8088314654

Linearity and SD and CV checked

Control: Biorad Lot Number :89730 Expiry :31-01-2026 (Control Level -1)

Parameter: AST Mean value : 41.1U/L Range :33 -49.2 U/L

Sl No	Values
1	41.1
2	40.9
3	39.7
4	40.5
5	40.2
6	40.8
7	40.5
8	40.9
9	41.8
10	41.3
11	39.5
12	39.7
13	39.8
14	40.7
15	41.7
16	40.8
17	40.8
18	40.7
19	40.8
20	41.2

Calculated Mean: 40.67 SD : 0.64

CV:1.56

Laniary Check Passed

Acceptable Limit for CV is :2.7



Kommunje House Manila Village, Post Murva, Bantwal Taluk, Dakshina Kannada - 574243



AnyScanner



Linearity and SD and CV checked

Control: Bio-Rad Lot Number: 89730 Expiry :31-01-2026 (Control Level 1)

Parameter: ALT Mean value : 26.5 Range :21.1 -31.9 U/L

Sl No	Values
1	26.9
2	26.5
3	27.1
4	27.3
5	27.4
6	26.8
7	26.5
8	26.4
9	26.8
10	27.2
11	27.1
12	27.3
13	26.6
14	26.8
15	26.7
16	26.9
17	26.9
18	26.8
19	27.1
20	27.0

Calculated Mean:26.91 SD :0.28 CV:1.05

Laniary Check Passed



Acceptable Limit for CV is :2.7



Lyphochek Assayed Chemistry Control Levels 1 and 2



C-310-5 Level 1 12 x 5 mL
C-315-5 Level 2 12 x 5 mL
313X MiniPak 2 x 5 mL



Level 1 89731
Level 2 89732



<https://www.myinserts.com/89730>

INTENDED USE

Lyphochek Assayed Chemistry Control is intended for use as an assayed quality control serum to monitor the precision of laboratory testing procedures for the analytes listed in this package insert.

SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Multiple levels of control are available to allow monitoring of the test system's reliability.

For customers in Germany: Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (Rili-BÄK regulation).

REAGENT

This product is prepared from human serum with added chemicals, purified biochemical material (tissue extracts of human and animal origin), therapeutic drugs, stabilizers and preservatives. This product is provided in lyophilized form for increased stability.

STORAGE AND STABILITY



This product will be stable until the expiration date when stored unopened at 2 to 8°C.

Reconstituted and Refrigerated: After reconstituting and storing tightly capped at 2 to 8°C, this product will be stable as follows:

- All analytes: 7 days

Except:

- T3 (Free), Acid Phosphatase (Total) and Prostatic Acid Phosphatase (PAP): 3 days

Reconstituted and Frozen: When reconstituted and stored tightly capped at -10 to -20°C, this product will be stable as follows:

- All analytes: 30 days

Except:

- Tobramycin: 20 days

- T3 (Free): 10 days

Once thawed, do not refreeze this product. Discard the remaining material.

This product is shipped under refrigerated conditions.

PROCEDURE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

Using a volumetric pipet or equivalent, reconstitute each vial with 5.0 mL of distilled or deionized water. Replace the stopper and allow this product to stand for approximately 20 minutes swirling occasionally.

Before sampling, gently swirl the vial several times to ensure homogeneity. If performing trace metal analysis, do not mix by inversion. After each use, promptly replace the stopper and return to the appropriate storage condition.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

LIMITATIONS

1. This product should not be used past the expiration date.
2. If there is evidence of microbial contamination or excessive turbidity in the reconstituted product, discard the vial.
3. This product is not intended for use as a standard.

WARNING



Biological source material. Treat as potentially infectious.

Each human donor unit used to manufacture this product was tested as required by FDA accepted methods. Tests results were non-reactive or negative for evidence of infection due to Human Immunodeficiency Virus (HIV), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV). This product may also contain other human source materials for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

Hazard (H) and Precautionary (P) Statements

EJH208 Contains Gentamicin, sulfate (salt) May produce an allergic reaction.

Safety Data Sheet (SDS) available for professional users on www.bio-rad.com.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a freeze-dried product manufactured under rigid quality control standards. To obtain consistent assay values, the control requires proper storage and handling as described.

ASSIGNMENT OF VALUES

The mean values and corresponding $\pm 3SD$ ranges in the Assignment of Values Data Charts (available separately) were derived from replicate analyses and are specific for this lot of product. Data from Unity™ Interlaboratory Program are included in the determination of some ranges. The tests listed were performed by the manufacturer and/or independent laboratories using manufacturer supported reagents and a representative sampling of this lot of product. It is recommended that each laboratory establish its own acceptable ranges and use those provided only as guides. Laboratory established ranges may vary from those listed during the life of this control. [Customers in Germany have to follow the requirements as described in the Rili-BÄK regulation.] Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation and reagents, or by manufacturer test method modifications.

Unity™ Interlaboratory Program is a proprietary Bio-Rad software program holding more than 2 billion QC data points from thousands of laboratories.

INSTRUCTIONS FOR OBTAINING THE DATA CHARTS

The Data Charts are available through the Internet, at www.myinserts.com/89730. Follow the directions at the website to receive email notifications of insert updates. Alternate methods for receiving data charts are available by contacting your local Bio-Rad Laboratories Office.

IMPORTANT PRODUCT INFORMATION



Significant changes are highlighted!



INSTRUMENT (1)

	Units	Level 1 - 89731		Level 2 - 89732	
		Mean	Range	Mean	Range
BIOSYSTEMS B400/A25/A15 (5) (continued)					
Magnesium (Xylidyl blue)	mg/dL	2.24	1.79 - 2.89		3.63 - >4.00
Phosphorus (Phosphomolybdate method)	mg/dL	3.54	2.90 - 4.18	7.45	6.11 - 8.79
Protein Serum (Total) (Biuret)	g/dL	6.94	6.11 - 7.77	4.54	4.00 - 5.08
Triglycerides (Glycerol phosphate oxydase/peroxydase)	mg/dL	156	133 - 179	87.3	74.2 - 100
Urea (Urease, UV) (4)	mg/dL	32.4	27.5 - 37.3	92.2	78.4 - 106
Uric Acid (Uricase, colorimetric)	mg/dL	4.52	3.84 - 5.20	8.82	7.50 - 10.1
ERBA XL SERIES (5)					
Albumin	g/dL	4.14	3.51 - 4.77	2.71	2.29 - 3.13
Albumin (India Market Only)	g/dL	4.25	3.62 - 4.88	2.72	2.30 - 3.14
Alkaline Phosphatase (IFCC Optimized) (2)	U/L	121	97.0 - 145	488	389 - 587
Alkaline Phosphatase (India Market Only) (2)	U/L	143	113 - 173	571	457 - 685
ALT/SGPT (2)	U/L	25.8	20.7 - 30.9	94.2	75.3 - 113
ALT/SGPT (Modified IFCC) (India Market Only) (2)	U/L	26.5	21.1 - 31.9	95.1	76.2 - 114
Amylase (2)	U/L	58.3	46.6 - 70.0	367	295 - 439
AST/SGOT (2)	U/L	41.1	33.0 - 49.2	204	162 - 246
AST/SGOT (Modified IFCC) (India Market Only) (2)	U/L	41.9	33.5 - 50.3	204	162 - 246
Bilirubin (Direct)	mg/dL	0.510	0.390 - 0.630	1.04	0.770 - 1.31
Bilirubin (Direct) (DCA)	mg/dL	0.520	0.400 - 0.640	1.37	1.04 - 1.70
Bilirubin (Total)	mg/dL	1.22	0.920 - 1.52	4.94	3.71 - 6.17
Bilirubin (Total) (DCA)	mg/dL	1.28	0.950 - 1.61	4.90	3.67 - 6.13
Calcium	mg/dL	9.28	7.90 - 10.7	12.3	10.5 - 14.1
Calcium (India Market Only)	mg/dL	9.35	7.94 - 10.8	12.2	10.4 - 14.0
Carbon Dioxide (CO2)	mEq/L	36.4	29.2 - 43.6	18.5	14.9 - 22.1
Chloride	mEq/L	97.5	87.6 - 107	84.7	76.3 - 93.1
Chloride (Mercuric Thiocyanate)	mEq/L	97.3	82.6 - 112	90.3	76.8 - 104
Cholesterol (HDL) (7)	mg/dL	59.1	47.4 - 70.8	24.6	19.8 - 29.4
Cholesterol (HDL) (India Market Only)	mg/dL	79.2	63.3 - 95.0	27.4	22.0 - 32.8
Cholesterol (LDL)	mg/dL	124	100 - 148	59.4	47.4 - 71.4
Cholesterol (LDL) (India Market Only)	mg/dL	145	115 - 175	70.7	56.6 - 84.8
Cholesterol (Total)	mg/dL	237	201 - 273	104	89.0 - 119
Cholinesterase (2)	U/L	7680	6144 - 9216	2125	1699 - 2551
Creatine Kinase (CK) (2)	U/L	138	111 - 165	484	388 - 580
Creatinine (Alkaline picrate method)	mg/dL	2.50	1.99 - 3.01	5.60	4.49 - 6.71
Creatinine (Enzymatic)	mg/dL	1.65	1.32 - 1.98	5.23	4.18 - 6.28
Creatinine (Enzymatic) (India Market Only)	mg/dL	1.63	1.30 - 1.96	5.42	4.34 - 6.50
Gamma Glutamyltransferase (GGT) (2)	U/L	60.8	48.5 - 73.1	165	132 - 198
Glucose	mg/dL	87.5	74.3 - 101	303	258 - 348
Glucose (Hexokinase) (India Market Only)	mg/dL	85.2	72.3 - 98.1	294	249 - 339
Glucose (HK)	mg/dL	84.2	71.6 - 96.8	297	252 - 342
Iron	µg/dL	205	163 - 247	56.7	45.3 - 68.1
Iron (IIBC)	µg/dL	115	91.0 - 139	136	109 - 163
Lactate Dehydrogenase (LDH) (2)	U/L	336	270 - 402	832	667 - 997
Lipase (2)	U/L	57.4	46.0 - 68.8	73.3	58.6 - 88.0
Lithium	mEq/L	0.660	0.600 - 0.720	2.10	1.89 - 2.31
Magnesium	mg/dL	2.15	1.82 - 2.48	4.18	3.55 - 4.81
Phosphorus (India Market Only)	mg/dL	3.68	3.14 - 4.22	6.90	5.85 - 7.95
Phosphorus (New Formulation)	mg/dL	3.61	3.07 - 4.15	7.09	6.04 - 8.14
Potassium	mEq/L	3.74	3.38 - 4.10	5.58	5.01 - 6.15



INSTRUMENT (1)

	Units	Level 1 - 89731		Level 2 - 89732	
		Mean	Range	Mean	Range
ERBA XL SERIES (5) (continued)					
Protein Serum (Total)	g/dL	6.53	5.54 - 7.52	4.23	3.60 - 4.96
Sodium	mEq/L	149	134 - 164	127	115 - 139
Triglycerides	mg/dL	153	129 - 177	93.1	79.0 - 107
Triglycerides (India Market Only)	mg/dL	160	136 - 184	87.3	74.1 - 101
Urea (4)	mg/dL	34.6	29.5 - 39.7	102	87.0 - 117
Uric Acid	mg/dL	5.21	4.43 - 5.99	10.7	9.20 - 12.2
Uric Acid (Single Reagent)	mg/dL	5.88	5.01 - 6.75	10.0	8.50 - 11.5
FUJI DRI-CHEM Series (5)					
Albumin	g/dL	6.00	5.60 - >6.00	4.20	3.60 - 4.80
Alkaline Phosphatase (ALP-PH) (IFCC) (2)	U/L	52.0	41.0 - 63.0	334	289 - 379
Alkaline Phosphatase (ALP-PH) (JSCC) (2)	U/L	164	132 - 196	996	864 - 1127
ALT/SGPT (2)	U/L	19.0	13.0 - 26.0	76.0	67.0 - 85.0
Amylase (2)	U/L	80.0	68.0 - 92.0	315	279 - 351
AST/SGOT (2)	U/L	39.0	33.0 - 46.0	210	188 - 232
Bilirubin (Direct)	mg/dL	0.200	<0.100 - 0.400	1.40	1.20 - 1.60
Bilirubin (Total)	mg/dL	1.10	0.800 - 1.40	4.00	3.20 - 4.80
Calcium	mg/dL	9.40	8.10 - 10.7	12.8	11.2 - 14.3
Carbon Dioxide (CO2)	mEq/L	38.4	29.0 - >40.0	17.9	13.0 - 23.0
Chloride	mEq/L	89.0	83.0 - 96.0	85.0	78.0 - 91.0
Cholesterol (HDL) (7)	mg/dL	93.0	75.0 - >110	32.0	23.0 - 41.0
Cholesterol (Total)	mg/dL	209	187 - 231	98.0	86.0 - 110
Cholinesterase (2)	U/L	253	220 - 287	67.0	58.0 - 76.0
Creatine Kinase (CK) (2)	U/L	140	119 - 161	515	448 - 583
Creatinine	mg/dL	1.30	0.900 - 1.70	4.70	4.20 - 5.20
Gamma Glutamyltransferase (GGT) (2)	U/L	48.0	37.0 - 58.0	184	159 - 209
Glucose	mg/dL	71.0	62.0 - 81.0	258	234 - 283
Lactate Dehydrogenase (LDH) (2)	U/L	142	123 - 161	298	259 - 337
LAP Amylase (2)	U/L	16.0	<10.0 - 31.0	18.0	<10.0 - 33.0
Lipase (2)	U/L	45.0	29.0 - 61.0	103	87.0 - 119
Magnesium	mg/dL	3.60	3.10 - 4.10	5.40	4.80 - 6.00
Phosphorus	mg/dL	3.10	2.70 - 3.60	7.50	6.60 - 8.40
Potassium	mEq/L	3.90	3.20 - 4.60	5.90	5.40 - 6.40
Protein Serum (Total)	g/dL	5.10	4.60 - 5.70	3.50	3.00 - 4.00
Sodium	mEq/L	151	145 - 157	130	125 - 135
Triglycerides	mg/dL	104	86.0 - 123	68.0	54.0 - 82.0
Urea Nitrogen (BUN) (4)	mg/dL	16.2	14.2 - 18.2	47.0	42.6 - 51.4
Uric Acid	mg/dL	4.20	3.70 - 4.80	10.1	9.20 - 11.1
ILa6 Taurus (5)					
Albumin	g/dL	4.30	3.66 - 4.95	2.70	2.43 - 2.97
Alkaline Phosphatase (2)	U/L	125	100 - 150	493	404 - 582
ALT/SGPT (2)	U/L	24.6	17.6 - 31.6	88.5	72.6 - 104
Amylase (2)	U/L	68.6	54.9 - 82.3	399	339 - 459
AST/SGOT (2)	U/L	39.2	31.4 - 47.0	185	152 - 218
Bilirubin (Direct)	mg/dL	0.390	0.190 - 0.590	1.07	0.800 - 1.33
Bilirubin (Total)	mg/dL	1.20	0.960 - 1.44	4.64	3.80 - 5.48
Calcium (o-cresolphthalein complexone)	mg/dL	9.13	8.22 - 10.0	12.5	11.3 - 13.8
Chloride	mEq/L	92.8	83.5 - 102	79.1	71.2 - 87.0
Cholesterol (HDL) (7)	mg/dL	62.7	50.2 - 75.2	24.5	19.6 - 29.4

