



XN-L Series Calibration Report

Customer : Unipath Specialty Laboratory (Gujarat) LLP
Model: XN-350

Date: 2/12/2023
Serial No: 17249

1 MECHANICAL ADJUSTMENTS

Positions	Perform
Piercer/Pipette Position	Performed
Tube Holder Position	NA
Hand Position	NA
BCR Position	NA

2 BACKGROUND VERIFICATION

Parameters	Results	Acceptable Range	Status
WBC	0.00	$\leq 0.10 \times 10^3/\mu\text{L}$	PASS
RBC	0.00	$\leq 0.02 \times 10^6/\mu\text{L}$	PASS
HGB	0.0	$\leq 0.1 \text{ g/dL}$	PASS
PLT-I	0	$\leq 10 \times 10^3/\mu\text{L}$	PASS
RBC-O*	NA	$\leq 0.02 \times 10^6/\mu\text{L}$	NA
PLT-O*	NA	$\leq 10 \times 10^3/\mu\text{L}$	NA

* Only applicable when RET license is activated

3 OPTICAL AXIS VERIFICATION

3.1 CELLPACK DFL 10ml + DUKE 4207A 5 Drops. The diluted latex was analyzed to determine the Optical Axis of the FSC (Rough). The results are as below :

Parameters	Results	Acceptable Range	Status
FSC (X) Rough	NA	70 - 130	NA
FSC (W) Rough	NA	Less than 0.15	NA

3.2 CELLPACK DCL 10ml + DUKE 4207A 5 Drops. The diluted latex was analyzed to determine the Optical Axis of the FSC (Fine). The results are as below :

Parameters	Results	Acceptable Range	Status
FSC (X) Fine	92.3	70 - 130	PASS
FSC (W) Fine	0.119	Less than 0.35	PASS

3.3 CELLPACK DCL 10ml + PS-FLUORED-L2830-4 2 Drops. The diluted latex was analyzed to determine the Optical Axis of the SFL(Fine) & SSC(Fine). The results are as below :

Parameters	Results	Acceptable Range	Status
SFL (X)	118.3	For Record Purpose	
SFL (W)	0.127	lower than 0.20	PASS
SSC (X)	60.3	50 - 100	PASS
SSC (W)	0.166	For Record Purpose	

NOTE: Please attach FSC(X) rough, FSC(X) fine, SFL(X) fine, SSC(X) fine peaks screenshots.



4 RBC CLOG LEVEL ADJUSTMENT

4.1 The RBC clog adjustment is monitored after the background check. The result is as below :

Parameters	Results	Acceptable Range	Status
RBC CLOG	100	100 ± 1	PASS

5 HGB BLANK ADJUSTMENT

Parameters	Results	Acceptable Range	Status
HB Blank Gain	5054	5000 ± 200	PASS

6 ASPIRATION SENSOR BLANK ADJUSTMENT

Parameters	Results	Acceptable Range	Status
Aspiration Sensor	5101	5000 ± 200	PASS

7 ASPIRATION SENSOR SPAN ADJUSTMENT

7.1 *Material : XN-CHECK / XN CAL

Lot :- 32381402

Expiry :- 05/12/2023

* Cross whichever not applicable

Parameters	Results	Acceptable Range	Status
Aspiration Sensor Span	13064	13000 ± 500	PASS

8 SENSITIVITY ADJUSTMENTS

8.1 Material :XN-CAL

Lot :- 33032101

Expiry :- 03/12/2023

Parameters	Results	Assay Value	Acceptable Range	Status	Scattergram normal?
WDF-X	163.5	163.5	± 3 ch	PASS	Yes
WDF-Y	95.4	94.4	± 3 ch	PASS	
WDF-Z / WDF-FSC	113.7	114.5	± 2.5 ch	PASS	

Parameters	Results	Assay Value	Acceptable Range	Status	Scattergram normal?
RBC-X	NA	NA	± 1 ch	#VALUE!	NA
RBC-Y	NA	NA	± 3 ch	#VALUE!	
RBC-Z	NA	NA	± 2 ch	#VALUE!	

Parameters	Results	Assay Value	Acceptable Range	Status
MCV	80.6	81.3	± 1 ch	PASS
MPV	9.8	9.9	± 0.3 ch	PASS
HGB	12.3	12.30	± 0.2g/dl	PASS

NOTE: Please attach sensitivity adjustment screenshots.



9 Calibration

9.1 Please refer to Whole Blood Calibration Report Attached

9.2 PD Mode Calibration (NOTE: If Applicable)

Material: XN CAL (1:7 Dilution)
Lot:- 33032101
Expiry:- 03/12/2023

Parameters	1	2	3	4	5	Average	Assay Value	Old Cal	New Cal	% Diff	Status
RBC_PD_CAL	4.38	4.35	4.32	4.40	4.30	4.35	4.402	1030	1042	-1.20	PASS
PLT_PD_CAL	252	239	242	248	246	245	247.9	980	990	-1.02	PASS
HGB_PD_CAL	12.3	12.2	12.2	12.4	12.3	12.3	12.30	1010	1012	-0.16	PASS
WBC_PD_CAL	7.27	7.32	7.41	7.50	7.59	7.42	7.593	970	993	-2.36	PASS
WBC-D/WDF_PD	6.87	6.85	6.94	7.04	7.15	6.97	7.100	1005	1024	-1.87	PASS
RET_PD_CAL*	NA	NA	NA	NA	NA	#DIV/0!	NA	NA	#VALUE!	#VALUE!	#VALUE!
RBCO_PD_CAL*	NA	NA	NA	NA	NA	#DIV/0!	NA	NA	#VALUE!	#VALUE!	#VALUE!
PLTO_PD_CAL*	NA	NA	NA	NA	NA	#DIV/0!	NA	NA	#VALUE!	#VALUE!	#VALUE!

**NOTE: Only Applicable when RET license has been activated*
NOTE: Please attach PD Mode Calibration results screenshots

9.3 Body Fluid Calibration (NOTE: If Applicable)

Material: XN CAL
Lot:- 33032101
Expiry:- 03/12/2023

Parameters	1	2	3	4	5	Average	Assay Value	Old Cal	New Cal	% Diff	Status
RBC_BF_CAL*	4.292	4.36	4.36	4.29	4.23	4.307	4.402	1010	1032	-2.20	PASS
WBC_BF_CAL*	7.093	7.04	7.06	7.12	7.04	7.070	7.100	978	982	-0.43	PASS

**NOTE: Only Applicable when BF license has been activated*
NOTE: Please attach BF Mode Calibration results screenshots

10 QC Verification

10.1 Please attach QC radar charts print outs for QC runs after the calibration.

11 Pipetors/Dilutors reproducibility and accuracy checked.



12 Certification

We certify that the **XN-350** Automated Hematology Analyzer S/N: **17249** has been successfully calibrated in accordance with the manufacturer's recommendations.

Report and Commissioning Performed By :

Signature (Engineer 1)

Name: Mehul Ekavadiya
Date: 2/12/2023

Signature (Engineer 2)

Name: _____
Date: _____

Report Reviewed and Accepted By :

Signature (Customer)

Name: Ms. Pinkal Patel
Date: 2/12/2023