LAB MONTHLY SUMMARY

Lab Name

BALAJI MEDICAL CENTRE

Lab No

3280

Month

December

Year

2019

Constituent Group

HbA1c

Date of Result Entered: 18/12/2019

Date of Report Published: 04/01/2020

SINO Analyte Method

HbA1c

OTHER IMMUNO ASSAYS

357

Interpretation

7.30

11.06

0.81

6.9 %

-0.50

0.09

SDI Range Within -1.0 to +1.0

Excellent. Between ± 1.0 to ± 2.0 Good.

Between ± 2.0 to ± 3.0 Accept with caution. Warning Signal.

Beyond ±3.0 Unacceptable performance. Action Signal.

Page 1 of 1

Homogeneity and Stability of the sample is passed.

Panela Christudos

Dr. Pamela Christudoss CMC EQAS Co-Ordinator

Christian Medical College, Vellore

****** End of Report ******