

## LAB MONTHLY SUMMARY

Lab Name BALAJI MEDICAL CENTRE Lab No 8227  
 Month January Year 2020  
 Constituent Group Chemistry II

Date of Result Entered : 13/01/2020

Date of Report Published : 11/02/2020

Sl.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants CV	SD	Your Value	SDI	U
1	GLUCOSE I	GOD-POD	OTHERS ( Any other Analyzer )	850	134.43	9.69	13.03	128.6 mg/dl	-0.45	0.89
2	UREA I	UREASE BERTHELOT	OTHERS ( Any other Analyzer )	281	39.56	14.36	5.68	35.1 mg/dl	-0.79	0.68
3	CREATININE I	JAFFE-END POINT	OTHERS ( Any other Analyzer )	350	1.77	13.35	0.24	2 mg/dl	0.97	0.03
4	T.BILIRUBIN I	DIAZONIUM SALT ( Colorimetric ) /JENDRASSIK	OTHERS ( Any other Analyzer )	597	2.04	21.38	0.44	2.6 mg/dl	1.28	0.04
5	T-PROTEIN I	BIURET - colorimetric	OTHERS ( Any other Analyzer )	690	5.22	11.48	0.60	5.2 g/dl	-0.03	0.05
6	ALBUMIN I	BCG - colorimetric	OTHERS ( Any other Analyzer )	674	3.19	10.36	0.33	3.1 g/dl	-0.27	0.03
7	URIC ACID I	ENZYMATIC / URICASE Colorimetric	OTHERS ( Any other Analyzer )	631	5.96	15.93	0.95	4.2 mg/dl	-1.85	0.08
8	CHOLESTEROL I	CHOD-PAP	OTHERS ( Any other Analyzer )	776	100.70	13.79	13.88	89.2 mg/dl	-0.83	1.00
9	TRIGLYCERIDE I	GPO-PAP / Enzymatic Colorimetric / End Point	OTHERS ( Any other Analyzer )	766	115.15	15.91	18.32	111.6 mg/dl	-0.19	1.32

SDI Range	Interpretation
Within -1.0 to +1.0	Excellent.
Between $\pm 1.0$ to $\pm 2.0$	Good.
Between $\pm 2.0$ to $\pm 3.0$	Accept with caution. Warning Signal.
Beyond $\pm 3.0$	Unacceptable performance. Action Signal.

Page 1 of 1

Homogeneity and Stability of the sample is passed.  
 Data in CMC EQAS reports is confidential

Contact details:  
 Email: clinqc@cmcvellore.ac.in  
 Contact Number: 0416-2283102

*Pamela Christudoss*  
 Dr. Pamela Christudoss  
 CMC EQAS Co-Ordinator  
 Christian Medical College, Vellore

\*\*\*\*\* End of Report \*\*\*\*\*

## VIEW LAB MONTHLY SUMMARY

<b>Lab Name</b>	BALAJI MEDICAL CENTRE		<b>Details About Robust Analysis</b>							
<b>Lab No</b>	8227		<b>Detail About Monthly Summary</b>							
<b>Month</b>	January		<b>Detail about SDI</b>							
<b>Year</b>	2020									
<b>Constituent Group</b>	HbA1c									
Click on the analyte to view Graphical Data			<input type="button" value="All Method Result"/>		<input type="button" value="Print"/>					
			Date of Result Entered : 13/01/2020							
			Date of Report Published : 11/02/2020							
Sl.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants CV	SD	Your Value	SDI	U
1	HbA1c	TURBIDIMETRIC INHIBITION IMMUNO ASSAY	OTHERS ( Any other Analyzer )	171	5.33	11.42	0.61	4 %	-2.18	0.09
<b>SDI Range</b>		<b>Interpretation</b>								
Within -1.0 to +1.0		Excellent.								
Between $\pm 1.0$ to $\pm 2.0$		Good.								
Between $\pm 2.0$ to $\pm 3.0$		Accept with caution. Warning Signal.								
Beyond $\pm 3.0$		Unacceptable performance. Action Signal.								

**Links**

<a href="#">Modify Address</a>	<a href="#">Modify Methods</a>	<a href="#">Submit Lab Results</a>	<a href="#">View Lab Monthly Summary</a>
<a href="#">Yearly Summary</a>	<a href="#">Change Password</a>	<a href="#">Home</a>	<a href="#">Sign-out</a>