



**VIEW LAB MONTHLY SUMMARY**

**Lab Name**            **SHRADHHADEEP GREENCROSS PATHOLOGY**  
**Lab No**                **5524**  
**Month**                 **January**  
**Year**                    **2020**  
**Constituent Group**   **HbA1c**

[Details About Robust Analysis](#)

[Detail About Monthly Summary](#)

[Detail about SDI](#)

[Click on the analyte to view Graphical Data](#)

[All Method Result](#)

[Print](#)

Date of Result Entered : 22/01/2020

Date of Report Published : 11/02/2020

Sl.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants		Your Value	SDI	U
						CV	SD			
1	HbA1c	TURBIDIMETRIC INHIBITION IMMUNO ASSAY	Roche C111,C311,Integra 400 plus,C501, C502 ,C702	81	4.73	9.53	0.45	4.51 %	-0.49	0.10

SDI Range	Interpretation
Within -1.0 to +1.0	Excellent.
Between ±1.0 to ±2.0	Good.
Between ±2.0 to ±3.0	Accept with caution. Warning Signal.
Beyond ±3.0	Unacceptable performance. Action Signal.

Links			
Modify Address	Modify Methods	Submit Lab Results	View Lab Monthly Summary
Yearly Summary	Change Password	Home	Sign-out



*Site Maintained by*

**Swiftware Solutions**

Best Viewed at pixels 1024 X 768