# Ortho Clinical Diagnostics

403, Leela Business Park, Andheri Kurla Road, Andheri East, Mumbai – 400059 T:+91 22 6787 9300 F:+91 22 6787 9333

### **CALIBRATION CERTIFICATE**

CERTIFICATE NO: OCD/PRIMA DIAGNOSTICS/01 CALIBRATION DATE: 20/11/2019
CALIBRATED FOR: PRIMA DIAGNOSTICS CALIBRATION DUE: 20/05/2020
LOCATION: MYSORE CALIBRATED AT: IMMUNOLOGY

#### SPECIFICATION OF INSTRUMENT:

INSTRUMENT	MODEL	SERIAL NO
Vitros	VITROS ECI	30004168
Immunodiagnostics System	J NUMBER	J30004168

*Specification of Software:* 

SOFTWARE NAME	SOFTWARE VERSION
QNX os	V3.9

The Reference of Calibration: The performance of various subsystems has been calibrated/tested by trained site engineer with pre-designed calibrated tools provided for the particular subsystems by the company.

#### **Installation Procedure Checklists:**

Subsystem Name	Adjusted/Verified	
1.Sample Supply Inspection	$\checkmark$	
2.Sample Metering	$\checkmark$	
* Sample Tray Tip Pickup & Sampling Center		
* Auxiliary Tip Supply		
* Sample Dispense		
* Outer Ring Mapping		
* Sample Metering Leak Test		
3.Reagent Supply	$\overline{\checkmark}$	
* Carousel to Well Shuttle		
* Well Shuttle – Positioning		
* Pack Opener (Home/Open/Close)		
* Outer Reagent – Alignment		
* Inner Reagent – Alignment		
* Well Top Alignment		
4.Reagent Metering	$\checkmark$	
* Wash Station		
* Dispense		
-		

Ortho Clinical Diagnostics India Private Limited

## Ortho Clinical Diagnostics

403, Leela Business Park, Andheri Kurla Road, Andheri East, Mumbai – 400059 T:+91 22 6787 9300

F: +91 22 6787 9333

_	<b>T</b>		
`	Proce	ccino	Center
$\mathcal{L}$	11000	عسرون	Contor

\* Inner Lift Pin ( Home/Shuttle/Read)

\* Inner/Outer rings Home Position

\* Incubator Shuttle

\* Incubator Temperature calibration

6. Well Wash Assembly

\* Well Wash Probe Horizontal Alignment

\* Inner Ring Mapping

\* Station Verification

7. Signal Reagent Assembly

\* Signal Reagent Probe – Horizontal Alignment

\* SR Dispense

\* SR Station Verification

8.Luminometer

| 🗸 |

\* Full Calibration (optional)

\* IRS Calibration

9.Fluid Supply

abla

\* Pressure Regulator Calibration

\* Supply Level Calibration

The results of comparison are as follows:

SUBSYSTEM NAME	CALIBRATION	ATION ACCURACY	
	RANGE	OBTAINED	
SAMPLE METERING	AUTOMATIC		
REAGENT METERING	4.90ML-5.10ML	5.0	
DISPENSE CALIBRATION			
REAGENT METERING OUTER	2.05ML-2.55ML	2.2	
PROBE WASH			
WELL WASH ASPIRATION	1.8ML-2.4ML	2.1	
WELL WASH DISPENSE	6.2ML-6.4ML	6.2	
SIGNAL REAGENT DISPENSE	9.99ML-10.01ML	10.0	
PUMP A			
SIGNAL REAGENT DISPENSE	9.99ML-10.01ML	10.0	
PUMP B			

# Ortho Clinical Diagnostics

403, Leela Business Park, Andheri Kurla Road, Andheri East, Mumbai – 400059 T:+91 22 6787 9300 F:+91 22 6787 9333

**CALIBRATED BY** 

KARTHIK K (ENGINEER ) Ortho Clinical Diagnostics

NOTE:

- 1. This certificate refers only to the particular item submitted for calibration.
- 2. The calibration result reported in the certificate is valid at the time of and under the stated condition of the measurement.