

Name Ravi Kumar K
Department Customer Service
HC LD

Telephone +91 40 30922661
Fax +91 40 23243146
Mobile +91 9642327866
E-mail ravi.rk.kumar@siemens.com

Your letter of
Our reference
Date July 15th, 2019

CERTIFICATE OF CALIBRATION

Calibration Date: 16/07/2019 Next Calibration: 15/07/2020

Name of the Instrument : Urine Analyzer

Model : CLINITEK STATUS+

Make : Siemens Healthcare Pvt Ltd,

Serial No. : 304086

Customer Address : VIJAYA DIAGNOSTIC CENTRE, NANDYAL

R. S. Kumar
Calibrated by: 16/07
2019

G. M. S. Sathak
16/7/2019
Reviewed by:

Siemens Healthcare Pvt Ltd.

4th Floor, 5-9-19,
Lakshmi Narasinh Estate,
Opp. T. S. Secretariat, Saifabad,
Hyderabad - 500004, TS, India

Tel.: +91 022-3967 7000
Fax: +91 022-3967 7075
Toll free: 1800-209-1800
Email: hc_contact.india@siemens.com
www.siemens.co.in/healthcare

Registered office: 130, Pandurang Budhkar Marg, Worli, Mumbai – 400018; Corporate Identity number: U74999MH2015PTC264859
Sales Offices: Ahmedabad, Bengaluru, Bhopal, Chandigarh, Chennai, Coimbatore, Gurgaon, Hyderabad, Jaipur, Kochi, Kolkata, Lucknow, Mumbai, Nagpur, Patna, Pune.

Calibration Parameters:

S.No.	Test/Parameter	Observation	Specification Pass/Fail
1	Display	Operating correctly	Pass
2	Touch Screen	Test Passed	Pass
3	Printer	Test Print	Pass
4	Test Table	Test Passed	Pass
5	Light Source	Test Passed	Pass
6	Electronics	Test Passed	Pass

S.No	Test/Parameters	Observation	Specification
1	Positive Control	Glu: +++	1+ to 3+
		Bil: +++	2+ to 3+
		Ket: ++	1+ to 3+
		SG: 1.015	1.000 – 1.015
		Blo: +++	2+ to 3+
		pH: 6.0	5.5 to 8.5
		Pro: +++	2+ to 3+
		Uro: >= 8.0 E.u./dl	(4.0 to >=8.0) E.u./dl
		NIT: Positive	Positive
		LEU: Positive	Small - Large
2	Negative control	Glu: Negative	Negative
		Bil: Negative	Negative
		Ket: Negative	Negative
		SG: 1.020	1.010 – 1.025
		Blo: Negative	Negative
		pH: 6.5	6.5 – 8.0
		Pro: Negative	Negative
		Uro: 0.2 E.u./dl	(0.2 – 1) E.u./dl
		NIT: Negative	Negative
		LEU: Negative	Negative

The instrument is working satisfactorily, Subsequent to calibrations of the above parameters.

Note: *Check-Stix Combo Pack* Urinalysis Control is specifically designed for urine dipstick testing approved AG kit used in performing the system calibration in all Clinitek Status Urine analyzers.

Authorized Signatory

Ravikumar Kadambari

Regional Customer Care Manager

R. S. Kadambari
16/07
2019

Siemens Healthcare Pvt Ltd.

4th Floor, 5-9-19,
Lakshmi Narasinh Estate,
Opp. T. S. Secretariat, Saifabad,
Hyderabad - 500004, TS, India

Tel.: +91 022-3967 7000
Fax: +91 022-3967 7075
Toll free: 1800-209-1800
Email: hc_contact.india@siemens.com
www.siemens.co.in/healthcare

Registered office: 130, Pandurang Budhkar Marg, Worli, Mumbai – 400018; Corporate Identity number: U74999MH2015PTC264859
Sales Offices: Ahmedabad, Bengaluru, Bhopal, Chandigarh, Chennai, Coimbatore, Gurgaon, Hyderabad, Jaipur, Kochi, Kolkata, Lucknow, Mumbai, Nagpur, Patna, Pune.