LAB MONTHLY SUMMARY

Lab Name

DR.SATAV'S PATHOLOGY LABORATORY

Lab No

2438

Month

September

Year

2020

Constituent

Group

HbA1c

Date of Result Entered: 02/09/2020

Date of Report Published: 07/10/2020

SI.No Analyte

Method / Principle Name

Analyzer Name

No of **Participants**

DV

Participants SD

Your Value

U SDI

HbA1c

HPLC-ION EXCHANGE

Biorad (D10 / D100 /

5.73 473

C٧ 7.35

5.3 % 0.42

-1.02

0.04

	Variant Turbo)
SDI Range	Interpretation
Within -1.0 to +1.0	Excellent.
Between ± 1.0 to ± 2.0	Good.
Between ± 2.0 to ± 3.0	Accept with caution. Warning Signal.
Beyond ±3.0	Unacceptable performance. Action Signal.

Page 1 of 1

Homogeneity and Stability of the sample is passed.

Data in CMC EQAS reports is confidential

Contact details:

Email:clinqc@cmcvellore.ac.in Contact Number: 0416-2283102

Panela Christudoss

Dr. Pamela Christudoss

CMC EQAS Co-Ordinator Christian Medical College, Vellore

****** End of Report *****