



PRIME AGENCIES

for your lab needs...

| | | |
|--|--|---|
| Instrument Service Report | | Ref. No: 010 |
| Customer Name : <i>Pranav Hospital</i> Address : <i>Madinaguda Hyderabad.</i> | <u>Message Received</u> By: <i>M. Sreenivas</i> | <u>Work carried out at</u> <input checked="" type="checkbox"/> Site <input type="checkbox"/> Service Centre |
| | Date: Time: | |
| Contact Person : <i>Mrs. Aruna</i> | Visit Date : <i>20/06/2020</i> | <u>Nature of Visit</u> |
| Contact No : <i>-</i> | Time in : <i>10:30 AM</i> | <input type="checkbox"/> Installation |
| E Mail ID : <i>-</i> | Time out : <i>12:00 PM.</i> | <input checked="" type="checkbox"/> Preventive Maintenance |
| | | <input checked="" type="checkbox"/> Applications |
| | | <input type="checkbox"/> Breakdown |
| | | <input type="checkbox"/> Emergency |
| | | <input type="checkbox"/> Others |

| | | |
|--|------------------|-------------|
| Warranty <input checked="" type="checkbox"/> | Service Contract | Charge Call |
|--|------------------|-------------|

| | | |
|----------------------------------|----------------------|-----------------------------|
| Particulars of Instrument | | Customer's Complaint |
| Instrument Name | <i>Medonic M32 B</i> | <i>Applications.</i> |
| Serial No. | <i>111 306</i> | |

| | | | |
|--|----------------------|-------------|---------------|
| Description of workdone/ Remarks <i>pm done and calibrated with Boule calibrator Lot.No. 2200434 and performed controls all the reports are satisfactory.</i> | Part Replaced | | |
| | Art No. | Description | Free / To pay |
| | | | |

Warranty : *-*

Customer Remarks :

| | | |
|--------------------|--------------------|-----------|
| Service Charge Rs. | Spare Replaced Rs. | Total Rs. |
|--------------------|--------------------|-----------|

This is to Certify that the above mentioned instrument has been installed/serviced/ Referred to service centers and working satisfactory

| | | |
|---|---------------------------|--|
| <i>M. Sreenivas</i> Name & Signature of Engineer | Date <i>20/06/2020</i> | <i>Aruna</i> Signature of Customer & Seal |
|---|---------------------------|--|

Franchisee:

H.No. 8-3-167/A/34, Vikaspuri Colony, Near A.G. Colony, SR Nagar, Hyderabad - 500 038.
INDIA Mobile: 9392995440, 9392995442, 9392995443, e-mail: primeagencies.bio@gmail.com

Boule Cal

CAL

IVD



Calibrator

| | |
|----------|------------|
| Lot | |
| 22004-34 | 2020-07-21 |

Open vial stability 5 days ¹

Medonic

| Parameter ² | RBC 10 ¹² /L | MCV fL | Plt 10 ⁹ /L | MPV fL | WBC 10 ⁹ /L | Hgb g/dL |
|------------------------|-------------------------|--------|------------------------|--------|------------------------|----------|
| CA620/530 | 4,18 | 77,1 | 260 | 9,8 | 8,9 | 12,0 |
| Medonic M-series / M32 | 4,18 | 80,0 | 240 | 9,8 | 8,6 | 12,1 |
| Range ± ³ | 0,10 | 2,0 | 15 | 1,5 | 0,3 | 0,2 |

Swelab

| Parameter ² | RBC 10 ¹² /L | MCV fL | Plt 10 ⁹ /L | MPV fL | WBC 10 ⁹ /L | Hgb g/dL |
|-------------------------|-------------------------|--------|------------------------|--------|------------------------|----------|
| AC9xxeo+ | 4,18 | 81,8 | 237 | 10,5 | 7,8 | 12,0 |
| Swelab Alfa / Alfa Plus | 4,18 | 81,1 | 245 | 10,3 | 8,6 | 12,1 |
| Range ± ³ | 0,10 | 2,0 | 15 | 1,5 | 0,3 | 0,2 |

1
bg Трайност на отворен флакон
hr Stabilnost otvorene bočice
cs Stabilita otevřené ampulky
da Åben r r stabilitet
et Stabiilsus avatud vialalis
fr Stabilit  en flacon ouvert
de Stabilit t ge ffneter Flaschen
el Διατήρηση δείγματος μετά την αποσφράγιση
it Stabilit  della fiala aperta
lv Atv rt flakonu stabilit ti
lt Stabilumas atidarius buteliuka
no Åpen r r/glass stabilitet
pl Trwa ość po otwarciu fioleki
pt Estabilidade ap s abertura do frasco
ro Termenul de valabilitate al fiolei desf cute
ru Стабильность вскрытого флакона
sr Proizvod stabilan nakon otvaranja pakovanja
sk Stabilita otvorenej injekcnej liekovky
sl Stabilnost odprte stekleni ke
es Estabilidad de la c psula abierta
sv H llbarhet f r  ppnad flaska
tr A ık ŐiŐe stabilitesi

2
Paramet r
Parametar
Parametr
Parameter
Parameeter
Param tre
Parameter
Par metro
Parametro
Parametrs
Parametras
Parameter
Kryterium
Par metro
Parametru
Parametr
Parametar
Parameter
Par metro
Parameter
Parametre

3
Obxvat
Raspon
Rozsah
Omr de
Vahemik
Intervalle
Bereich
Αναμενόμενο  ρος  
Intervallo
Diapazons
Intervalas
Omr de
Zakres
Intervalo
Interval
диапазон
Opseg
Rozsah
Obmo je
Intervalo
Intervall
Aralıđı



<http://www.medonic.se/>
www.swelab.com/

M. S. S. S.
20/06/2020

Boule
Boule Medical AB
Domnarvsgatan 4
SE-163 53 Sp nga, Sweden

25123-3

1504055



**maternityTM
matters**

Exclusive maternity wing from Pranaam Hospitals
(A unit of Pranaam Hospital Pvt. Ltd.)



Calibration Report

Printed: 20-06-2020 10:47 AM
Operator ID:
Instrument S/N: 111306
Instrument ID: 111306
Software Version: 2.3.2r
Authorization Level: Service

Calibration: Whole Blood
Calibrated by:
Calibration Time: 20-06-2020 10:47 AM

Calibration Result

| Param | Frev. | Current | Method | CV (%) |
|-------|-------|---------|--------|--------|
| RBC: | -3.7 | -4.0 | AT | 0.4 |
| MCV: | -0.5 | -9.3 | AT | 0.3 |
| PLT: | -11.2 | -19.1 | AT | 2.1 |
| MPV: | -20.3 | -22.5 | AT | 0.8 |
| HGB: | -17.8 | -17.6 | AT | 0.7 |
| WBC: | -6.7 | -16.4 | AT | 0.7 |
| RDW%: | | 6682 | | |
| RDW: | | 6800 | | |
| PDW: | | 7300 | | |
| PDW%: | | 3400 | | |

Sample Selection

| SEQ | ID |
|-------|----------|
| 12856 | 2200434+ |
| 12855 | 2200434+ |
| 12854 | 2200434+ |
| 12853 | 2200434+ |
| 12852 | 2200434+ |

Aruna
20/06/2020

RBC (10¹²/l)

| CV (%) | Mean | Target |
|--------|-------|--------|
| 0.4 | 4.19 | 4.18 |
| SEQ | Value | |
| 12856 | 4.19 | |
| 12855 | 4.19 | |
| 12854 | 4.21 | |
| 12853 | 4.20 | |
| 12852 | 4.17 | |

M. S. S.
20/06/2020

MCV (fl)

| CV (%) | Mean | Target |
|--------|-------|--------|
| 0.3 | 87.7 | 80.0 |
| SEQ | Value | |
| 12856 | 88.0 | |
| 12855 | 88.0 | |
| 12854 | 87.8 | |
| 12853 | 87.6 | |
| 12852 | 87.3 | |

Dr. madhavi