



Date: 06/1/2019

This is to certify for CLINITEK Status SL No: 286950 Urine chemistry analyzer performs system diagnostics check whenever instrument is switched 'ON'.


Additionally clinitek status also performs automatic self calibration check before and each test and this ensure satisfactory performance of the system for CLINITEK Status analyzer Siemens reading reagent strips.

The system doesn't generate any printed data for calibration & as it doesn't proceed with testing if calibration failed.

Hence for this instrument no need of calibration report

Thanking you,

For Siemens Healthcare Diagnostics Ltd.,

  
Raghavendra.D  
Key Account Manager – Service