

LAB MONTHLY SUMMARY

Lab Name SRI BABA DIAGNOSTIC CENTRE Lab No 13304
 Month November Year 2020
 Constituent Chemistry II
 Group

Date of Result Entered : 19/11/2020

Date of Report Published : 05/12/2020

Sl.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants		Your Value	SDI	U
						CV	SD			
1	GLUCOSE I	GOD-POD	Any Analyser	912	373.63	10.97	40.98	512 mg/dl	3.38	2.71
2	UREA I	UREASE BERTHELOT	Any Analyser	265	57.83	11.79	6.82	71 mg/dl	1.93	0.84
3	CREATININE I	JAFFE RATE-BLANKED & COMPENSATED/Alkaline Picrate	Any Analyser	448	7.76	13.08	1.01	4.8 mg/dl	-2.93	0.10
4	T.BILIRUBIN I	DIAZONIUM SALT (Colorimetric) / JENDRASSIK	Any Analyser	597	1.80	21.22	0.38	1.9 mg/dl	0.26	0.03
5	CHOLESTEROL I	CHOD-PAP	Any Analyser	767	113.22	11.86	13.43	212 mg/dl	7.36	0.97
6	TRIGLYCERIDE I	GPO-PAP / Enzymatic Colorimetric / End Point	Any Analyser	723	159.55	13.88	22.15	189 mg/dl	1.33	1.65

SDI Range	Interpretation
Within -1.0 to +1.0	Excellent.
Between ±1.0 to ±2.0	Good.
Between ±2.0 to ±3.0	Accept with caution. Warning Signal.
Beyond ±3.0	Unacceptable performance. Action Signal.

Page 1 of 1

Homogeneity and Stability of the sample is passed.

Data in CMC EQAS reports is confidential

Contact details:

Email: clinqc@cmcvellore.ac.in

Contact Number: 0416-2283102



Dr. Pamela Christudoss
 CMC EQAS Co-Ordinator
 Christian Medical College, Vellore

***** End of Report *****

	SRI BABA DIAGNOSTIC CENTRE	Form No.	SBDC / QAS-BC/230
	CORRECTIVE ACTION REPORT	Issue No.	01
		Rev No.	0
		Issue Date	12.1.2020

CAR No.: CAR/01 Date: 1.12.2020

Sample (Process (System (Cust. Complaint (
----------	-----------	----------	-------------------

Detail of Non-Conformance
We used Biosystems Analyser for Glucose & cholesterol. So the value was High.
N.S.I.

Signature: Date: 01.12.2020

Root Cause Analysis
Started using Anyo Analyser to get correct Report.
N.S.I.

Signature: Date: 01.12.2020

Corrective Action
We got calibration certificate & used biochemistry Control to get correct Result.
N.S.I.

Signature: Date: 01.12.2020

Verification & Effectiveness of Corrective Action
Analyser was serviced & got calibration certificate & counselling was given by company's service Manager.
N.S.I.

Signature: Date: 01.12.2020