M. Sadhasivam M. Phil. (Micro), PGDMLT

#40/25, Paloothu Street, Pookadai Corner, Tiruchengode. € 94874 93110, 94432 33110, skalpanamicrolab@gmail.com

08/12/2020

TO WHOMSOEVER IT MAY CONCERN

Dear Sir/Madam,

We are using MICRO SCOPE-BLISCO INDIA, S.No: 4786 in our laboratory. We have attached this equipment calibration certificate. Equipment photos we are attached. Kindly accept our request and make us favour.

Thank you.

Yours Sincerely

M.SADHASIVAM., MSc (Micro)., MPhil (Micro)., PGDMLT

Proprietor

KALPANA MICROBIOLOGY LAB











CALIBRATION CERTIFICATE

Name of the customer:

& Address

KALPANA MICROBIOLOGY LAB

No.40/25, Pookadai Corner, South

car Street, Thiruchengode -637211

Calibration Date

26-12-2020

Calibration validity

25-12-20201

Unit Descriptions:

Microscope

Blisco India Make:

S.No:

4786

Traceability details				
S.No	Test Equipment	Make	Valid Up To	
_1	Electrical safety analyzer	Rigelsafe test60	10-02-2021	
2	Digital Lux Meter	HTC /Lx -103	13-08-2021	
3	Thermo hygrometer	нтс	10-02-2021	

Room Temperature: 25°C ± 2 Room Humidity: 30 -50 %

(At the time of calibration)

ESA	Test	Result

S.No	Parameter	ESA Reading	Unit	Normal Range	Condition
1	Current	90	Ma	1 = 1	Normal
2	Frequency .	50.0	Volts	49-50	Normal
3	Mains voltage line – Neutral	230.3	Volts	220-240	Normal
4	Mains voltage line – Earth	230.3	Volts	220-240	Normal
5	Mains voltage Neutral – Earth	4.8	Volts	< 7	Normal

S.No	Knob Position	DUC Reading (LUX)	Distance (Cm)
1	Low	3270	30
2	High	19110	30

DATE: 26-12-2020

TEST BY: G. DINESH



NOTE: 1) Calibration Points & due date is given as per customer request 2) This certificate refer only to particular item submitted for certificate3) The calibration results reported are valid at the time of and under the stated conditions of the measurements4) Any error in the certificate should be brought to our lab within 45 days from the date of issue of certificate.