RESOLVING QC FAILURE

Steps taken to resolve out of range values for the month of January

- 1. Fresh reagents were prepared and samples were re-run.
- 2. Sample storage, handling and instrument operation was reviewed.
- 3. Maintenance of instruments were done prior to sample analysis.
- 4. Internal quality control was done.
- 5. Staff training was performed.





10

home.cmcvellore.ac. C





CHRISTIAN MEDICAL COLLEGE

In Service of the Nation Since 1900 ...



Lab Name SOBHI DIAGNOSTIC CENTRE

2021

Lab No 2739

Month January

Constituent Group Chemistry II

Click on the analyte to view Graphical Data

Details About Robust Analysis

Detail About Monthly Summary

Detail about SDI

All Analyser Result Print

Date of Result Entered: 17/01/2021

Date of Report Published: 05/02/2021

Si.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants		Your	SDI	U
						CV	SD	Value	301	
1	GLUCOSE II	GOD-POD II	Robonik (Semi automatic / Fully Automatic)	379	118.64	9.39	11.14	123.6 mg/dl	0.45	1.14
2	CREATININE II	Jaffes Kinetic- Alkaline Picrate	Robonik (Semi automatic / Fully Automatic)	263	1.30	13.52	0.18	1.4 mg/dl	0.57	0.02
3	T.BILIRUBIN II	Diazonium Salt (Colorimetric) / Jendrassik II	Robonik (Semi automatic / Fully Automatic)	233	2.01	23.12	0.46	0.9 mg/dl	-2.39	0.06
4	T-PROTEIN II	Biuret - Colorimetric II	Robonik (Semi automatic / Fully Automatic)	173	5.56	10.63	0.59	5 g/dl	-0.95	0.09
5	ALBUMIN II	BCG - Colorimetric II	Robonik (Semi automatic / Fully Automatic)	167	3.33	11.35	0.38	2.5 g/dl	-2.20	0.06
6	URIC ACID II	Enzymatic / Uricase Colorimetric II	Robonik (Semi automatic / Fully Automatic)	193	5.70	16.69	0.95	6.5 mg/dl	0.84	0.14
7	CHOLESTROL II	CHOD-PAP II	Robonik (Semi automatic / Fully Automatic)	353	114.48	12.54	14.36	160.2 mg/dl	3.18	1.53
8	TRIGLYCERIDE II	GPO-PAP / Enzymatic Colorimetric / End Point II	Robonik (Semi automatic / Fully Automatic)	198	142.26	14.28	20.32	151 mg/dl	0.43	2.89

SDI Range	Interpretation		
Within -1.0 to +1.0	Excellent.		
Between ±1.0 to ±2.0	Good.		
Between ±2.0 to ±3.0	Accept with caution. Warning Signal.		
Beyond ±3.0	Unacceptable performance. Action Signal.		

Links					
Modify Address	Modify Methods	Submit Lab Results	View Lab Monthly Summary		
Yearly Summary	Change Password	Home	Sign-out		



Best Viewed at pixels 1024 X 768







