

**CERTIFICATE OF CALIBRATION**

ADL/BO/SER/0100/20-21

DATE : 02/11/2017	SITE : Sparsha Diagnostics, Mangalore	CITY : Bangalore
STATE : KARNATAKA	COUNTRY : INDIA	S/N. : AAA15X2088
INSTRUMENT NAME	TBA25FR	

Section	Work	Clean / Check	Lubricate	Replace
Reaction Tray	Cleaning of Reaction try	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Sampling Pump	Cleaning	<input type="checkbox"/> <input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Reagent Pump	Cleaning	<input type="checkbox"/> <input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Probe Washing Pump	Cleaning	<input type="checkbox"/> <input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Reaction Cuvette washing	Cleaning	<input type="checkbox"/> <input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Sample Probe	Cleaning	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Reagent Probe	Cleaning	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Movement Check				
Reaction Tray	Checking of Cooling Temp.	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Sampling Pump	Movement Check	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Reagent pump	Movement Check	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Probe washing Check	Movement Check	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Reaction cuvette washing	Movement Check	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Sampling probe	Movement Check	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Reagent probe	Movement Check	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Rack Detector	Function Check	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Automated Sampler	Checking on reading of Sample Rack	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Stirrer Check	Up & Down Movement, Rotation	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Water supply & Drainage	Movement Check on Drainage pump	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Body	Checking of leakage	<input checked="" type="checkbox"/>		<input type="checkbox"/>
Total Inspection	Cell blank Test	<input checked="" type="checkbox"/>		<input type="checkbox"/>

Authorized Signatory:

AGAPPE DIAGNOSTICS LTD. | ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT
 Agappe Hills, Pattimattom (PO), Bhd., Ernakulam, Kerala - 683 562, India.
 Tel: +91 484 286 7000 | Email: agappe@agappe.in

EQUIPMENT PLANT
 X/549 CB, Block No. 32, KINFRA Small Industrial Park,
 Kollat, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

KOLKATA OFFICE
 406, Merin Mitra, Plot No 10, Block DN, Sector V,
 Salt Lake City, Kolkata - 700 091.
 Tel: +91 31 4003 0451 | Email: kolkataoffice@agappe.in

MUMBAI (REGISTERED OFFICE)
 401 & 402, 4th Floor, Jansangh Business Centre, 119,
 Sahar Road, Parsipawadi, Andheri (E-93), Mumbai - 400 099, India.
 Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

DELHI OFFICE
 D5M 54D, 5th Floor, DLF Tower, Shivaji Marg,
 New Delhi - 110 003, India.
 Tel: +91 11 4558 8416 | Email: delhioffice@agappe.in

BANGALORE OFFICE
 5-6, 1st Floor, Red Cross Bhavan, No.26 Race Course Road,
 Bangalore - 560 001. Tel: +91 80 2228 8288
 Email: bangaloreoffice@agappe.in




Temperature Check		Reaction Tray 37 Celsius	
Cuvette Check			
Wave Length	Gain	ADC	Range
340 nm	1.6	4116	8
380 nm	1.6	4153	8
404 nm	1.6	4199	6
416 nm	1.6	4246	6
450 nm	1.6	4197	7
476 nm	1.6	4161	7
500 nm	1.6	4206	7
524 nm	1.6	4218	6
548 nm	1.6	4188	8
572 nm	1.0	2618	5
604 nm	1.6	4192	6
628 nm	1.0	2605	6
660 nm	1.0	2585	5
700 nm	1.0	2588	6
748 nm	0.9	2364	6
804 nm	0.9	2432	5

This is to certify that the above mentioned product has been calibrated according to the standard operating procedure given by the Toshiba Medical Systems Corporation, Japan (Doc #2B586-297EN)

Calibrator used	Multicalibrator	Lot No.:	<u>30070039</u>	Exp.:	<u>12/21</u>
Control used	Norm	Lot No.:	<u>26471</u>	Exp.:	<u>28/2/23</u>
	Path	Lot No.:	<u>26472</u>	Exp.:	<u>28/2/23</u>
	Test	CV%			
	TP	0.8%		Date of Calibration:	<u>14/March/2021</u>
	GOT	0.6%		Calibration Validity:	<u>13/Sep/2021</u>

For AGAPPE DIAGNOSTICS LTD.

Authorized Signatory
Name: Sanoj T K
Designation: Product Specialist



AGAPPE DIAGNOSTICS LTD. ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT
Agappe Hills, Pattimattom IPO, Dist. Ernakulam, Kerala - 683 562, India.
Tel: +91 484 256 7000 | Email: agappe@agappe.in

EQUIPMENT PLANT
X/588 CB, Block No. 32, KINFRA Small Industrial Park, Nellikat, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

KOLKATA OFFICE
406, Merit Matrix, Plot No 10, Block DN, Sector V, Salt Lake City, Kolkata - 700 071.
Tel: +91 31 4001 0451 | Email: kolkataoffice@agappe.in

MUMBAI (REGISTERED OFFICE)
403 & 402, 4th Floor, Jansingh Business Centre, 119, Safar Road, Parsiwala, Andheri (East), Mumbai - 400 099, India.
Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

DELHI OFFICE
D5M 549, 5th Floor, DLF Tower, Shivaji Marg, New Delhi - 110 003, India.
Tel: +91 11 4558 8410 | Email: delhioffice@agappe.in

BANGALORE OFFICE
5-6, 11th Block, Red Cross Bhawan, No.26 Race Course Road, Bangalore - 560 001. Tel: +91 80 2228 8288
Email: bangaloreoffice@agappe.in