

18	ALP	PNP AMP kinetic	plus, C501, C502, C702 Roche C111, C311, Integra 400 plus, C501, C502, C702	246	82.32	7.87	6.48	56.5 U/L	-3.98	0.83
SDI Range			Interpretation							
Within -1.0 to +1.0			Excellent.							
Between ± 1.0 to ± 2.0			Good.							
Between ± 2.0 to ± 3.0			Accept with caution. Warning Signal.							
Beyond ± 3.0			Unacceptable performance. Action Signal.							

Homogeneity and Stability of the sample is passed.

Data In CMC EQAS reports is confidential

Contact details:

Email: clinqcc@cmcvellore.ac.in

Contact Number: 0416-2283102

*Checked by
V. Raju
02/01/21*

Paola Christudoss

Dr. Pamela Christudoss

CMC EQAS Co-Ordinator

Christian Medical College,
Vellore

***** End of Report *****