



CERTIFICATE OF CALIBRATION

ADL/BO/SER/0100/20-21

DATE : 14/03/2021	SITE: Sparsha Diagnostics, Mangalore	CITY: Bangalore
STATE :KARNATAKA	COUNTRY: INDIA	S/N. : AAA15X2088
INSTRUMENT NAME	TBA25FR	

Section	Work	Clean / Check	Lubricat e	Replac e
Reaction Tray	Cleaning of Reaction try			
Sampling Pump	Cleaning			
Reagent Pump	Cleaning			
Probe Washing Pump	Cleaning			
Reaction Cuvette washing	Cleaning			
Sample Probe	Cleaning			
Reagent Probe	Cleaning			
Movement Check				
Reaction Tray	Checking of Cooling Temp.			
Sampling Pump	Movement Check			
Reagent pump	Movement Check			
Probe washing Check	Movement Check			
Reaction cuvette washing	Movement Check			
Sampling probe	Movement Check			
Reagent probe	Movement Check			
Rack Detector	Function Check			
Automated Sampler	Checking on reading of Sample Rack			
Stirrer Check	Up & Down Movement, Rotation			
Water supply & Drainage	Movement Check on Drainage pump			
Body	Checking of leakage			
Total Inspection	Cell blank Test			

Authorized Signatory: AGAPPE DIAGNOSTICS LTD. 150 9001:2015 [EN 150 13485:2016 CERTIFIED COMPANY] CIN : UZ4239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT Agappe Hills, Pattimattom (PO), Dist, Ernakulam, Kerata - 683 562, India. Tel: + 91 484 286 7000 i Email: agappe@agappe.in

MUMBAI [REGISTERED OFFICE] 401 & 402, 4 Floor, Jaisingh Business Centre, II9, Sahar Road, Parsiwata, Andheri (East), Mumbai - 400 099, India. Tel: +91 22 4300 8000 | Email: mumbaiofficespagappe.in X/588-CB, Block No. 32, KINFRA Small Industrial Park, Nellad, Cochin, Kerala, India - 686 721, Tel: +91 484 276 7477.

DELHI OFFICE DSW 540, 5" Floor, DLF Tower, Shivaji Marg, New Delhi - 110 OIS, India. Tel: +91 II 4558 8416 i Email: delhioffice@agappe.in

KOLKATA OFFICE 496, Merlin Matrix, Plot No-10, Block-DR, Sector V, Salt Lake City, Kolkata - 700-091, Tet: +91-33 4003 0451 | Email: Rolkattaoffice@agappe.in

BANGALORE OFFICE
5-6. II^{ac} Floor, Red Cross Bhavan, No-26 Race Course Road,
Bangalore - 560 001. Tel; +91 80 2228 828
Email: bangaloreoffice@agappe.in

CELEBRATING TO



Temperature	Reaction Tray 37 Celsius			
urang mu Check				
Cuvette Check				
Wave Length	Gain	ADC	Range	
340 nm	1.6	4116	8	
380 nm	1.6	4153	8	
404 nm	1.6	4199	6	
416 nm	1.6	4246	6	
450 nm	1.6	4197	7	
476 nm	1.6	4161	7	
500 nm	1.6	4206	7	
524 nm	1.6	4218	6	
548 nm	1.6	4188	8	
572 nm	1.0	2618	5	
604 nm	1.2	4192	6	
628 nm	1.0	2605	6	
660 nm	1.0	2585	5	
700 nm	1.0	2588	6	
748 nm	0.9	2364	6	
804 nm	0.9	2432	5	

This is to certify that the above mentioned product has been calibrated according to the standard operating procedure given by the Toshiba Medical Systems Corporation, Japan(Doc #2B586-297EN)

Calibrator used	Multicalibrator	Lot No.: <u>3</u>	<u>30070039</u>	Exp. :12/21
Control used	Norm	Lot No.: 2	<u>6471</u>	Exp.:28/02/2023
	Path	Lot No.: 2	<u>6472</u>	Exp.:28/02/2023
	Test	CV%		
	TP	0.8%	Date of Calibration: 14/03/2021	
	GOT	0.6%	Calibration	Validity: 13/09/2021

For AGAPPE DIAGNOSTICS LTD.

Authorized Signatory Name: Sanoj T K

Designation: Product Specialist



AGAPPE DIAGNOSTICS LTD. 150 9001:2015 [EN 150 13485:2016 CERTIFIED COMPANY] CIN : UZ4239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT Rappe Hills, Pattinattom (PO), Dist. Ernakulam, Kerala - 683 562, India. Tel. • 91 484 286 7000 i Email: agappegagappe.m

MUMBAI [REGISTERED OFFICE] 401 & 402, 4" Floor, Jaisingh Business Centre, 119, Sahar Road, Parsiwada, Andheri (East), Mumbai - 400 099, India. Tel: +91.22 4300 8000 | Email: mumbaiofikespagappe.in EQUIPMENT PLANT X/588-CB, Block No. 32, KINFRA Small Industrial Park, Nellad, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

DEM 540, 5" Floor, DLF Tower, Shivaji Marg, New Delhi - 110 QE, India. Tel: •91 II 4558 84161 Email: delhioffice@agappe.in

406. Merlin Matrix, Plot No-10, Block-DN, Sector V, Salt Lake City, Kolkata - 700 091, Tel: +91 33 4003 0451 | Email: kolkattaofficedingappe.in

BANGALORE OFFICE S-6, II" Floor, Red Cross Bhavan, No:26 Race Course Road, Bangalore - 560 001, Tel: +91 80 2228 8288 Email: bangaloreoffice@agappe.in