

**MSACS**

GOVT. OF INDIA - NACO  
NOT FOR SALE

**WARNING**

1. READ OPERATING INSTRUCTIONS BEFORE USE
2. PROVIDE PROPER EARTHING.
3. FOR VIBRATION FREE PERFORMANCE
  - a. ENSURE LEVEL AND STABILITY.
  - b. BALANCE CENTRIFUGE TUBES.
  - c. ENSURE USE OF RUBBER CUSHIONS FOR GLASS TUBES.
4. BRING SPEED KNOB TO 'OFF' AND INCREASE SPEED GRADUALLY.
5. CHECK CARBON BRUSHES PERIODICALLY.
6. CHECK GLASS FUSE FOR NON-WORKING OF UNIT.

**R-4** LABORATORY CENTRIFUGE

**REMI**



MAINS



**BIOHAZARD**



Speed

