**CERTIFICATE OF CALIBRATION**

ADL/BO/SER/0020/20-21

Instrument : Mispa Nano Fully Automated Biochemistry Analyser

Customer Name : Rajagiri Victor Hospital,Goa

Date of installation : 13/MAY/2018

Serial No : WN-83107096

Alignment Check

Motion Check	Checked	Remark
Sample Probe	<input checked="" type="checkbox"/>	Ok
Reagent Probe	<input checked="" type="checkbox"/>	Ok
Mixing Bar	<input checked="" type="checkbox"/>	Ok
Reaction Disk	<input checked="" type="checkbox"/>	Ok
Reagent Disk	<input checked="" type="checkbox"/>	Ok
Sample Disk	<input checked="" type="checkbox"/>	Ok
Fluidic disk	<input checked="" type="checkbox"/>	Ok
Lamp	<input checked="" type="checkbox"/>	Ok

System	Checked	Remark
Main unit	<input checked="" type="checkbox"/>	Ok
Temp. unit	<input checked="" type="checkbox"/>	Ok
Mixing unit	<input checked="" type="checkbox"/>	Ok
Reaction unit	<input checked="" type="checkbox"/>	Ok
Reagent unit	<input checked="" type="checkbox"/>	Ok
Sample unit	<input checked="" type="checkbox"/>	Ok

AGAPPE DIAGNOSTICS LTD. ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT
 Agappe Hills, Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India.
 Tel: +91 484 286 7000 | Email: agappe@agappe.in

EQUIPMENT PLANT
 X/588-CB, Block No. 32, KINFRA Small Industrial Park,
 Nellad, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

KOLKATA OFFICE
 406, Merlin Matrix, Plot No-10, Block-DN, Sector V,
 Salt Lake City, Kolkata - 700 091.
 Tel: +91 33 4003 0451 | Email: kolkattoffice@agappe.in

MUMBAI [REGISTERED OFFICE]
 401 & 402, 4th Floor, Jaisingh Business Centre, 119,
 Sahar Road, Parsiwada, Andheri (East), Mumbai - 400 099, India.
 Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

DELHI OFFICE
 DSM 540, 5th Floor, DLF Tower, Shivaji Marg,
 New Delhi - 110 015, India.
 Tel: +91 11 4558 8416 | Email: delhioffice@agappe.in

BANGALORE OFFICE
 5-6, IInd Floor, Red Cross Bhavan, No:26 Race Course Road,
 Bangalore - 560 001. Tel: +91 80 2228 8288
 Email: bangaloreoffice@agappe.in

**OPTICAL SECTION**

Filter	Dark Current	Background
340	102	28144
405	98	32400
450	89	33316
510	87	35769
546	94	38529
578	86	38676
630	89	39232
670	89	40539.
Reference		

Power Supply : 5V 12V 24V

The equipment is calibrated for the different parameters by using Multi Calibrator

Lot.No: 30020180 Exp: 12/2021.

This is to certify that the above mentioned product has been calibrated according to the standard operating procedure given by the MINDRAY Medical Systems.

Date of Calibration : 25/SEP /2020

Next Due : 24/MAR/2021


For Agappe Diagnostics Ltd.,

Authorised signatory

Name: Antony Jose

Designation: Engineer - Customer Service

Antony



AGAPPE DIAGNOSTICS LTD. ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT
Agappe Hills, Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India.
Tel: + 91 484 286 7000 | Email: agappe@agappe.in

EQUIPMENT PLANT
X/588-CB, Block No. 32, KINFRA Small Industrial Park,
Nellad, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

KOLKATA OFFICE
406, Merlin Matrix, Plot No-10, Block-DN, Sector V,
Salt Lake City, Kolkata - 700 091.
Tel: +91 33 4003 0451 | Email: kolkataoffice@agappe.in

MUMBAI (REGISTERED OFFICE)
401 & 402, 4th Floor, Jaisingh Business Centre, 119,
Sahar Road, Parsiwada, Andheri (East), Mumbai - 400 099, India.
Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

DELHI OFFICE
DSM 540, 5th Floor, DLF Tower, Shivaji Marg,
New Delhi - 110 015, India.
Tel: +91 11 4558 8416 | Email: delhioffice@agappe.in

BANGALORE OFFICE
5-6, IInd Floor, Red Cross Bhavan, No.26 Race Course Road,
Bangalore - 560 001. Tel: +91 80 2228 8288
Email: bangaloreoffice@agappe.in