

op: Incubator O/R: Reagent Sup



DIAGNOSTICS - Calibrate Luminometer

DAC VALUES

High Voltage Setting: 3569
Discriminator DAC Setting: 1699
Analog Offset DAC Setting: 2026
Luminometer Reference
Intensity DAC Setting: 2606
Incubator Reference
Intensity DAC Setting: 2846



IRS
Calibration



Full
Calibration

Inop: Incubator O/R: Reagent Sup



DIAGNOSTICS - Calibrate Luminometer

DAC Values

Review the DAC values from before and after the calibration.

Setting	Before	After
High Voltage Setting:	3569	3569
Discriminator DAC Setting:	1699	1699
Analog Offset DAC Setting:	2026	2026
Luminometer Reference Intensity DAC Setting:	2606	2606
Incubator Reference Intensity DAC Setting:	2864	2846

OK

Help

... was successful.



IRS
Calibration



Full
Calibration



CALIBRATION CERTIFICATE

This is to certify that the system maintenance & calibration was performed on **07th APRIL 2021** for the VITROS ECIQ IMMUNOASSAY system installed at **CIHSR** laboratory, 4th mile, Dimapur bearing Serial No.30005029

The system's calibration includes checking the reproducibility performance of the instrument as per the guidelines provided by the manufacture.

Validity: Six months for Optics Calibration.

For *Ortho Clinical Diagnostics*.



Monuj Saikia
Sr. Service Engineer