

**CERTIFICATE OF CALIBRATION**Instrument : **Mispa Ace Fully Automated Clinical Chemistry Analyser**Customer Name : **Meditech Diagnostic Centre, West Mambalam, Chennai**Date of Installation : **09/July/2020**

S/N:WK-95106413

Alignment Check

Motion Check	Checked	Remarks
Sample Probe	<input checked="" type="checkbox"/>	OK
Reagent Probe	<input checked="" type="checkbox"/>	OK
Mixing Bar	<input checked="" type="checkbox"/>	OK
Reaction Disk	<input checked="" type="checkbox"/>	OK
Reagent Disk	<input checked="" type="checkbox"/>	OK
Sample Disk	<input checked="" type="checkbox"/>	OK
Lamp	<input checked="" type="checkbox"/>	OK

System	Checked	Remark
Main unit	<input checked="" type="checkbox"/>	OK
Temp. unit	<input checked="" type="checkbox"/>	OK
Mixing unit	<input checked="" type="checkbox"/>	OK
Reaction unit	<input checked="" type="checkbox"/>	OK
Reagent unit	<input checked="" type="checkbox"/>	OK
Sample unit	<input checked="" type="checkbox"/>	OK

AGAPPE DIAGNOSTICS LTD. ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT
 Agappe Hills, Palomattam (PO), Dist. Ernakulam, Kerala - 683 562, India.
 Tel: +91 484 286 7000 | Email: agappe@agappe.in

EQUIPMENT PLANT
 X/SR-CE, Block No. 32, KINFRA Small Industrial Park,
 Nellik, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

KOLKATA OFFICE
 404, Merlin Matrix, Plot No-10, Block-04, Sector V,
 Salt Lake City, Kolkata - 700 091.
 Tel: +91 33 4003 0451 | Email: kolkataoffice@agappe.in

MUMBAI (REGISTERED OFFICE)
 401 & 402, 4th Floor, Jai Singh Business Centre, 119,
 Sahar Road, Parkside, Andheri (East), Mumbai - 400 099, India.
 Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

DELHI OFFICE
 DSM 540, 5th Floor, DLF Tower, Shivaji Marg,
 New Delhi - 110 015, India.
 Tel: +91 11 4558 8416 | Email: delhioffice@agappe.in

BANGALORE OFFICE
 5-6, 11th Floor, Red Cross Bhavan, No.26 Race Course Road,
 Bangalore - 560 001. Tel: +91 80 2328 8288
 Email: bangaloreoffice@agappe.in

**OPTICAL SECTIONS**

Filter	Dark Current	Background
340	93	55420
405	97	51232
450	95	53412
510	98	54530
546	100	55315
578	99	54312
630	113	53201
670	110	55052
Reference		

Power Supply : 5V 12V 24V

The equipment is calibrated for the different parameters by using Multi

Calibrator Lot. No: 30020180 Exp: 07-2021

This is to certify that the above mentioned product has been calibrated according to the standard operating procedure given by the MINDRAY Medical Systems.

Date of Calibration : **08/Jan/2021**

Next Due on : **07/Jul/2021**

For Agappe Diagnostics Ltd,

M.T.S

Authorised signatory

Raja M

Engineer - Customer Service

AGAPPE DIAGNOSTICS LTD. ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT
Agappe Hills, Pattimattam (PO), Otl, Ernakulam, Kerala - 683 562, India.
Tel: +91 484 296 7000 | Email: agappind@agappe.in

MUMBAI (REGISTERED OFFICE)
403 & 402, 4th Floor, Jai Singh Business Centre, 119,
Sahar Road, Parkwest, Andheri (West), Mumbai - 400 099, India.
Tel: +91 22 4 100 9000 | Email: mumbai@agappe.in

EQUIPMENT PLANT
X/580 CB, Block No. 22, KINFDA Small Industrial Park,
Nellikud, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

DELHI OFFICE
OSM 540, 1st Floor, B1 F Tower, Shreej Marg,
New Delhi - 110 025, India.
Tel: +91 11 4558 8416 | Email: delhi@agappe.in

KOLKATA OFFICE
406, Merlin Matrix, Plot No-10, Block DM, Sector V,
Salt Lake City, Kolkata - 700 091.
Tel: +91 33 4001 0451 | Email: kolkata@agappe.in

BANGALORE OFFICE
5-6, 1st Floor, Red Cross Bhawan, No.26 Race Course Road,
Bangalore - 560 001. Tel: +91 80 2228 8260.
Email: bangalore@agappe.in