

#### Beckman Coulter India Pvt.

Regional Office
The DLF Galleria bearing unit no. 1201,
DLF Galleria, Plot No. BG-08,
Premises No. 02-0124, Action Area IB,
New Town, Kolkata – 700156.

<u>Tel: +91</u> 33-40050630/2 Fax: +91 33-40050631

## **CALIBRATION CERTIFICATE**

We here by certify that the fully automated chemistry analyzer Beckman Coulter AU480 has been installed by us and Calibrated. It has passed all the following quality checks.

Instrument serial number: 2018011701

Installed at : Sri Venkateshwara Institute of Medical Sciences

**Tirupati** 

Last Calibration Done on: 11.06.2021 Next Calibration Due on: 10.06.2022

### **Quality Checks - Passed**

- . Reagent and sample aspiration system.
- . Reagent Tray & STAT Table Refrigeration.
- . Cuvette Wheel Incubation.
- . Lamp & Optical Unit.
- . Water supply and Tubing system.
- . Calibration of Probes, Rack, Track & Carousel Movements.
- . Reagent and Sample Position Identification.
- . Water test and QC checks.
- . Calibration of all test parameters.
- . QC run for all test parameters.

We here by confirm that AU480 has passed all the above quality checks and is working in good condition.



#### Beckman Coulter India Pvt.

Regional Office The DLF Galleria bearing unit no. 1201, DLF Galleria, Plot No. BG-08, Premises No. 02-0124, Action Area IB,

<u>Tel: +91</u> 33-40050630/2 Fax: +91 33-40050631

New Town, Kolkata – 700156.

## **CALIBRATION CERTIFICATE**

Instrument serial number: 2018011701

Installed at : Sri Venkateshwara Institute of Medical Sciences

Tirupati

Last Calibration Done on: 11.06.2021 Next Calibration Due on: 10.06.2022

## **Photometric check**

Requirement : All O.Ds should be less than 1.5

Wavelength	Optical Density	Result
340	1.0943	OK
380	1.0121	OK
410	0.9651	OK
450	0.9545	OK
480	0.8656	OK
520	0.8865	OK
540	0.8659	OK
570	0.8459	OK
600	0.8395	OK
660	0.8278	OK
700	0.8261	OK
750	0.8136	OK
800	0.8120	OK



#### Beckman Coulter India Pvt.

Regional Office
The DLF Galleria bearing unit no. 1201,
DLF Galleria, Plot No. BG-08,
Premises No. 02-0124, Action Area IB,
New Town, Kolkata – 700156.

<u>Tel: +91</u> 33-40050630/2 Fax: +91 33-40050631

## **Incubation Temperture**

Range: 37C +0.5

On display the temperature is 37.0 stable

S.NO	TEST	MEASURED/CALIBRATED	REMARKS
	PARAMETERS		
01	INPUT VOLTAGE	230 V	ok
02	POWER	4 KVA	ok
03	INPUT FREQUENCY	50 Hz	ok
04	OPERATING TEMP.	20-25 °C	ok

Tool Used: Digital Multimeter

Make: Fluke Model: 17B+

Sl.no: 39591421WS

We hereby confirm that Beckman Coulter AU480 has been successfully calibrated and passed all the above quality checks and is working in good condition.

For Beckman Coulter India Pvt Ltd

# T V SAPTHAGIRIDHAR CUSTOMER SUPPORT ENGINEER

#### **Beckman Coulter India Pvt. Ltd**

Registered Office: Unit Nos. TF-B-07(A) to 15, Third Floor, B Wing, Art Guild House, Phoenix Market City, LBS Road, Kurla (West), Mumbai – 400 070. Tel: +91 22 3080 5000 Fax: +91 22 3080 5060 CIN: U33119MH2005PTC157177 www.beckmancoulter.com