

VIEW LAB MONTHLY SUMMARY

Lab Name	COSMOS CLINICS AND DIAGNOSTICS		Details About Robust Analysis							
Lab No	10077		Detail About Monthly Summary							
Month	July		Detail about SDI							
Year	2021									
Constituent Group	HbA1c									
Click on the analyte to view Graphical Data			All Analyser Result	Print						
Date of Result Entered : 10/07/2021										
Date of Report Published : 02/08/2021										
Sl.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants CV	Participants SD	Your Value	SDI	U
1	HbA1c	Turbidimetric Inhibition Immunoassay	Any Analyser	236	8.07	13.89	1.12	5.4 %	-2.38	0.15
SDI Range		Interpretation								
Within -1.0 to +1.0		Excellent.								
Between ±1.0 to ±2.0		Good.								
Between ±2.0 to ±3.0		Accept with caution. Warning Signal.								
Beyond ±3.0		Unacceptable performance. Action Signal.								

Links

Modify Address	Modify Methods	Submit Lab Results	View Lab Monthly Summary
Yearly Summary	Change Password	Home	Sign-out