

VIEW LAB MONTHLY SUMMARY

Lab Name	COSMOS CLINICS AND DIAGNOSTICS				Details About Robust Analysis						
Lab No	10077				Detail About Monthly Summary						
Month	September				Detail about SDI						
Year	2021										
Constituent Group	HbA1c										
				Click on the analyte to view Graphical Data				All Analyser Result		Print	
								Date of Result Entered : 16/09/2021			
								Date of Report Published : 05/10/2021			
SI.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants CV	Participants SD	Your Value	SDI	U	
1	HbA1c	Turbidimetric Inhibition Immunoassay	Any Analyser	241	6.83	11.97	0.82	6.6 %	-0.28	0.11	
SDI Range				Interpretation							
Within -1.0 to +1.0				Excellent.							
Between ±1.0 to ±2.0				Good.							
Between ±2.0 to ±3.0				Accept with caution. Warning Signal.							
Beyond ±3.0				Unacceptable performance. Action Signal.							

Links

Modify Address	Modify Methods	Submit Lab Results	View Lab Monthly Summary
Yearly Summary	Change Password	Home	Sign-out