

In Field Calibration Report

Standard F 200
Serial No: FA20EBITG3392


Operator
Name: SD Biosensor Healthcare Pvt.Ltd.

Instrument Owner
Name: Cosmos Clinics and Diagnostics, Bengaluru

Control Devices	Lot No.	Expiry Date
CAL 1	FC0120002	15-08-22
CAL 2	FC0120002	15-08-22
CAL 3	FC0120002	15-08-22
C - A1C - 1	CT01210111 - L1	26-07-22
C - A1C - 2	CT01210112 - L2	26-07-22

Operating Conditions	
Temperature	15- 32°C

Control Device	Result of Control Device Measurements		Accepted Range
	Measured Value		
	1	2	
CAL 1	Ok	Ok	N/A
CAL 2	Ok	Ok	N/A
CAL 3	Ok	Ok	N/A
CT01210111 - L1	5.7	5.7	4.3 % - 6.7%
CT01210112 - L2	9.3	9.4	7.5 %-11.5 %

Final Result		
The operator hereby confirm that the In Field Calibrations was performed in accordance to the instructions given in "Standard F 200 , In Field Calibration Procedure".		
Date:	Authorised Operator Signature	Instrument Responsible Signature
21.10.2021		

SD BIOSENSOR

SD BIOSENSOR HEALTHCARE PVT. LTD.

Registered & Corporate Office: Unit # 202 A-D, 2nd Floor, Tower-A, Unitech Signature Towers, South City-I, Gurugram-122001
Tel.: (+91) 124-454-0907, Toll Free No: (+91) 1800 10 23105 Website: www.sdbiosensor.com : Email: care@sdbiosensor.co.in
CIN - U33110HR2012PTC047034

Installation Certificate for Standard F200

This is to certify that Standard F200 instrument Serial No. FA20EBITG3392 is successfully installed and commissioned at Cosmos Clinics and Diagnostics, Bengaluru and the installation protocol/checklist has been successfully completed for the above instrument.

For **SD Biosensor Health Care Pvt. Ltd.**




Assistant Manager Application

Installation Report

Customer Details	Instrument Details
------------------	--------------------

Name: Cosmos Clinics and Diagnostics

Model: Standard F 200

Address: # 603, SMV Arcade, Andrahalli main road, Syndicate Bank layout, Bengaluru - 560091

Serial Number: FA20EBITG3392

Contact Person(s): Dr. Nirupama.S

Product of SD Biosensor

Contact No: 9945670868

Installation Date: 16-08-2021

Accessories and Consumable Supplied

1. Power adapter: Yes
2. Calibration set: Yes
3. Test kit, Buffer: Yes
4. User instruction guide: Yes

Installation Checks

1. Environment: Closed, AC
2. Power Adapter: Working
3. Temperature: 26 degrees Celsius (in limits between 15 degrees to 32 degrees Celsius)
4. UPS: Yes

Training Imparted	Assays Performed
-------------------	------------------

1. Operating Procedure
2. Maintenance
3. Dos and Don'ts
4. Calibration

Customer Remarks	Installer's Remarks
------------------	---------------------

Customer's Signature



Installer's Signature

Operations Qualification for Standard F 200

1. Verified Cal 1 Calibration (blank check) : Done
2. Verified Cal 2 Calibration (UV light source and detector check) : Done
3. Verified Cal 3 Calibration (RBG light source and three level band
Intensity (Strong positive, weak positive and negative) check) : Done

For **SD Biosensor Health Care Pvt. Ltd.**



Vikas Chawla
Assistant Manager Application



SD BIOSENSOR HEALTHCARE PVT. LTD.

Registered & Corporate Office: Unit # 202 A-D, 2nd Floor, Tower-A, Unitech Signature Towers, South City-I, Gurugram-122001
Tel.: (+91) 124-454-0907, Toll Free No: (+91) 1800 10 23105 Website: www.sdbiosensor.com : Email: care@sdbiosensor.co.in
CIN - U33110HR2012PTC047034

Performance Certificate – Standard F 200

Instrument Certification:

Study data has determined that the instrument described in this document either meets all criteria outlined in this Performance Qualification Protocol or exceptional conditions have been identified and documentation included.

It is recommended to do Instrument calibration by calibration set once in 30 days for checking optical parts of analyzer. There is no preventive maintenance required in this equipment. SD Biosensor will update software of system free of cost as when required to run machine at optimum levels.

This instrument is ready for specific usage.

Protocol Performed By:




Vikas Chawla
Assistant Manager Application