

kshat Biomedicals (p) '

Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph.: +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

Installation Certificate For AXIOM-21 Sr. No. 3020ET01360

This is to certify that instrument AXIOM-21 sr. no.:- 3020ET01360 is successfully installed and Commissioned at, JAIPUR and the installation protocol/ checklist have been successfully completed for the above instrument.

Date of Installation : 23/04/2018

AKSHAT BIOMEDICALS (P) LTD, Technical services Department

Name

: SAGAR SHARMA

Designation : Service Engineer

Customers Name

: RELIABLE DIAGNOSTIC CENTRE

Address

: KOTA, RAJASTHAN

Instrument Model

: AXIOM-21

Serial no.

: 3020ET01360

Date of Installation

: 23/04/2018

Initial Inspection of the unit carried out and the details are as follows:

System Condition Report: Found the system to have been delivered in satisfactory condition and no External physical damaged observed on the same, Package was kept in good condition as per the Directional Indicators like not tilt, indicating the system has not been subjected to mechanical shocks or stored in any manner, so as to cause any damage to the same.

Found all the required accessories as required.

External Requirement for installation:

Input voltage of 220v-240v/50Hz or 60 Hz

2. Recommended operating Temperature is 15-30 degree Celsius, with in Relative Humidity 30-85%And Atmospheric pressure 70-106kpa.



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotrnail.com

8.0.: 95/G1. Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

Carried out all the	installation proc	edures as per	the installation	n procedure and	checklists
Carried out all the	installation proc	euules as pei	the mistanation	, bi occurre	

Carried out all the necessary checks and alignments'.

Carried out all the necessary system checks and tests.

Handover the instrument for operators Training And Qualifications

AKSHAT BIOMEDICALS (P) LTD., Technical Services Department

Name : SAGAR SHARMA

Designation: Service Engineer



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph.: +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

Installation Certificate for AXIOM-21 Sr. No.:- 3020ET01360

This is to certify that the instrument AXIOM-21 Sr. No.:- 3020ET01360 is successfully installed and

Commissioned at RELIABLE DIAGNOSTIC CENTRE, KOTA, RAJASTHAN and the Installation Protocol/ checklist has been

Successfully completed for the above Instrument.

Date of installation

: 23/04/2018

AKSHAT BIOMEDICALS (P) LTD., Technical Services Department

Name

: Sagar Sharma

Designation : Service Engineer

Installation Qualification for AXIOM-21 Sr. No.:- 3020ET01360

Customer Name

: RELIABLE DIAGNOSTIC CENTRE

Address

: KOTA, RAJASTHAN

Instrument Name

: AXIOM-21

Serial No.

: 3020ET01360

Initial Inspection of the unit carried out and the details are as follows

System Condition Report: Found the system to have been delivered in satisfactory condition and no External physical damaged observed on the same, Package was kept in good condition as per the Directional Indicators like not tilt, indicating the system has not been subjected to mechanical shocks or stored in any manner, so as to cause any damage to the same.

Found all the accessories as required.



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotrnail.com

B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph. + 91.141.2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

1. Verified all the Mechanical Movements

: Done

2. Verified Hydraulics

: Done

3. Verified Electrical Systems

: Done

4. Verified the all reagents System

: Done

AKSHAT BIOMEDICALS (P) LTD., Technical Service Department

Name

: Sagar Sharma

Designation

: Service Engineer

Date

: 23/04/2018

Instrument Setup

1. Assembled the instrument accessories.

2. Removed the shipping Clamps.

3. Connected the REAGENTS.

4. Connected the REGENT TUBING.

5. Connected the power cord and connection cord.

6. Initialise the machine and follow the installation procedure.

Operational Instection

1. Checked and found mechanical movements OK.

2. Checked and found Hydraulics OK.

3. Checked and found electrical OK.

4. Checked with Control and Samples, Results Are found OK.

AKSHAT BIOMEDICALS (P) LTD., Technical Service Department

Name :Sagar Sharma
Designation : Service Engineer



kshat Biomedicals (p) I

Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph.: +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

Checked and found all controls to be within the acceptable SD.

System Certification:

Study data has determined that the system described in this document either meets the necessary criteria outlined in this performance qualification protocol, or exceptional conditions have been identified documentation include.

The system is ready for specific usage.

Protocol Performed By:

Name

: Sagar Sharma

Designation

: Service Engineer

PERFORMANCE QUALIFICATION DONE BY

Name

: LOKESH SHARMA

Designation: DIRECTOR

Operational Qualification:

System Certification:

Study data has determined that the system described in this document either meets all criteria outlined in this Operational Protocol, or exceptional have been identified and documentation include.

Exceptional Conditions, if any have been addressed.

The System is ready for specific usage.

Protocol Performed By

: Name

: Sagar Sharma

: Designation: Service Engineer



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com

8.0.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph.: +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

CALIBRATION RECORD

INSTEUMENT MODEL	AXIOM-21	Sr. No.	3020ET01360
TYPE	AXIOM-21 AUTOMATED HEMATOLOGY ANALYZER	INSTALLED AT	RELIABLE DIAGNOSTIC CENTRE
CALIBRATOR USED	HQ - CAL	Calibrator Lot No	KK102
O/LIBIOTI OIL OOL		Calibrator Expiry	05/09/2021

Enclosed:-

1. Assay Sheet of HQ-Cal AXIOM Hematology (Lot- KK102, Expiry- 2021-09-05)

Calibration done

05/07/2021

Next Calibration Due -

04/07/2022

FOR AKSHAT BIOMEDICALS (P) LTD.

(Authorized signatory)



Regd. Office: 352 Katewa Nagar New Sanganer Road: Soriala: Jaipur-302019 + akshatinda@hotmail.com

B.O.: 95/G1: Devi Nagar: New Sanganer Road: Sodala: Jaipur-302019 • Ph. + 91:141:2297075 | Factory: 121: Solitaira Industrial Park: RIICO: Bagru (Jaipur-302019)

Calibrator

Lot No. KK102

EXPIRY

05-09-2021

Parameter	WBC	RBC	НСВ	MCV	PLT
Target	9.23	4.57	13.6	87.5	272
1	9.02	4.61	13.7	87.8	265
√2	9.14	4.65	14	87.6	266
√3	9.09	4.58	13.8	87.6	275
√4	8.93	4.68	13.9	87.9	266
v /5	9.14	4.54	13.8	87.8	271
٧6	9.23	4.63	13.9	87.6	264
√7	9.12	4.66	13.9	87.5	273
٧8	9.18	4.70	14.1	87.7	274
v9	9.19	4.56	13.8	88.0	264
V10	9.03	4.57	13.7	87.3	267
V11	9.12	4.54	13.8	87.5	269
Mean	9.11	4.61	13.84	87.66	269.25
CV(%)	0.95	1.11	0.8	0.2	1.8
Calibration Factor	99.00	101	101	100	98
Calibration Factor:	100.00	100.00	100.00	100.00	100.00



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com

8.0.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph : +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

External Requirement for Installation:

- 3. Input voltage of 220V-240V/50Hz or 60Hz.
- 4. Recommended operating Temperature is 15-30 degree Celsius with in relative Humility 30-85% and Atmospheric pressure 70-106kpa.

Installation Qualification for AXIOM-21 Sr. No. :- 3020ET01360

Carried out all the installation procedures as per the installation procedure and checklists.

Carried out all the necessary checks and alignments.

Carried out all the necessary system checks and Qualifications.

Handover the instrument for Operation Training and Qualification.

AKSHAT BIOMEDICALS (P) LTD., Technical Service Departments

Name : Sagar Sharma

Designation: Service Engineer



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com B.O.: 95/G1, Devi Nagar, NN Sola A Cold Addil (2) N 30200 LA La F4/22 Ad 5 (2) Algaire Industrial Park, RIICO, Bagru (Jaipur)

INSTRUMENT NAME:

AXIOM-21 AUTOMATED HEMATOLOGY ANALYZER

Model No. :-

AXIOM-21

SR. NO:-

3020ET01360

Prepared By:-

Sagar Sharma_

Date:- 23/04/2018

Customer Name

: RELIABLE DIAGNOSTIC CENTRE

Address

: KOTA, RAJASTHAN

Instrument Model

: **AXIOM-21**

Serial No.

: 3020ET01360

Date Of Installation

: 23/04/2018

The instrument was installed and was found to be working satisfactory. Preliminary Customer Training was provided and standardization of the parameters were done. The results were found to be within the expected range and system found to be working satisfactorily.

AkshatBiomedicals (P) Ltd

Technical Services Department

Name

: Sagar Sharma Customer Name: RELIABLE DIAGNOSTIC CENTRE

Designation: Service Engineer

Contact person: VIJAY NAGAR

Signature

Designation

: Lab Technician



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com

B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph.: +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

OPERATIONAL QUALIFICATION (OQ)

INSTRUMENT NAME: AXIOM-21 AUTOMATED HEMATOLOGY ANALYZER

Model No. :-

AXIOM-21

SR. NO :-

3020ET01360

Prepared By:- Sagar Sharma Date:- 23/04/2018

System Certification:

Study data has determined that the system described in this documents either meets all criteria outlined in this operational protocol, or exceptional conditions have been identified and documentation included.

Exceptional conditions, if any have been addressed.

The System is ready for specific usage.

Protocol Performance

: SagarShrama

Designation

: Service Engineer

Customer Authorization

: RELIABLE DIAGNOSTIC CENTRE KOTA, RAJASTHAN

Name: Sagar Sharma

Designation: Service Engineer

Company Representative Name And Sign

Customer Name And Sign

Date:- 23/04/2018 Date:-23/04/2018

Operations Qualification for AXIOM-21 Sr. No.:- 3020ET01360



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com

B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph.: +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

Performance Qualification for AXIOM-21 Sr. No.: 3020ET01360

Calibration Parameters

Checked and found all the control [parameters to be within the acceptable CV limits and in range.

Checked and found all the control to be within the acceptable SD.

System Certification:

Study data determined that the system described in this documents either meets the necessary criteria outlined in this performance Qualification Protocol, or exceptional conditions have been identified and documents included.

The system is ready for specific usage.

Protocol Performance By:

Name

: SAGAR SHARMA

Designation

: SERVICE ENGINEER

PERFORMANCE QUALIFICATION DONE BY

Name

: LOKESH SHARMA

Designation

: DIRECTOR



Regd. Office: 352, Katewa Nagar, New Sanganer Road. Sodala, Jaipur-302019 • akshatindia@hotmail.com B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur-302019 • Ph. +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

PROCEDURAL QUALIFICATION (PQ)

INSTRUMENT NAME:- AXIOM-21 FULLY AUTOMATED HEMATOLOGY ANALYZER

Model No. :-

AXIOM-21

SR. NO :-

3020ET01360

Prepared By:- Sagar Sharma Date:- 23/04/2018

Installation Qualification for AXIOM-21 Sr. No. : - 3020ET01360

Customer Name

: RELIABLE DIAGNOSTIC CENTRE

Address

: KOTA, RAJASTHAN

Instrument Model

: **AXIOM-21**

Serial No.

: 3020ET01360

Date of installation

: 23/04/2018

The instrument was installed and was found to be working satisfactory. Preliminary Customer Training was provided and standardization of the parameters were done. The results were found to be within the expected range and system found to be working satisfactorily.

AKSHAT BIOMEDICALS (P) LTD., Technical Services Departments

Name

: Sagar SharmaCustomer Name: RELIABLE DIAGNOSTIC CENTRE

Designation : Service Engineer

Contact Person: VIJAY NAGAR

Signature

Designation : Lab Technician



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur-302019 • akshatindia@hotmail.com B.O.: 95/G1, Devi Nagar, New Sanganer Road, Sodala, Jalpur-302019 • Ph.: +91 141 2297075 | Factory: 121, Solitaire Industrial Park, RIICO, Bagru (Jaipur)

Test Name: Precision Testing

Purpose: To comply with specification and find deviation If any, if required to apply correction.

Method: To run normal sample 11 times consecutively and identify SD, CV and establish that it lies as per the manufacture's specs. To omit first run and to consider rest 10 runs and calculate SD & CV.

Sample No.	WBC	RBC	HGB	НСТ	MCV	PLT
t-1	7.70	4.80	13.9	46.6	95.4	139
t-2	7.84	4.50	13.7	48.5	95.6	140
t-3	7.65	4.70	14.7	47.0	95.4	136
t-4	8.14	3.90	14.9	47.0	95.5	147
t-5	7.45	4.87	12.5	46.6	95.7	133
t-6	8.15	4.83	14.7	45.9	95.2	140
t-7	7.95	4.87	14.9	46.5	95.5	133
t-8	7.91	5.21	15.5	49.7	95.5	140
t-9	7.65	4.90	14.7	46.7	95.4	143
t-10	8.21	4.91	15.0	46.8	95.5	134
Mean	7.86	4.74	14.45	47.16	95.47	138.5
SD	0.6453	0.2327	0.2712	1.2670	1.2670	0.1337
CV%	1.538	2.451	1.823	2.686	0.140	3.286

Authorization:

Name

:SAGAR SHARMA

Designation: SERVICE ENGINEER

Date

23-04-2018



Regd. Office: 352, Katewa Nagar, New Sanganer Road, Sodala, Jaipur 302019 • akshatindia@hotmail.com B.O.: 95.G1, Devi Nagar, New Sanganer Road, Sodala, Jaipur 302019 • Ph. + 91.141.2297075 | Factory. 121. Solitaire Industrial Park, RIICO. Bagru (Jaipur)

NAME OF THE INSTITUTION:- RELIABLE DIAGNOSTIC CENTRE

CARRYOVER WORKSHEET.

Instrument Id: AXIOM-21, SN = 3020ET01360

Name of technical person: VIJAY NAGAR

Supervisor review: VIJAY NAGAR Date performed: 05/07/2021

Carryover calculation:

Carryover % = (L1-L3)/(H3-L3)*100

Interpretation of results.:

1) If (L1-L3) is equal to a Negative number or zero Then stop as there is no Evidence of carryover.

The percentage
 Carryover effect should be
 Less than the imprecision
 Of the assay listed by
 Manufacture.

				Analytes		
Instrument	WBC	RBC	Plt	Hb	Hct	
H1	14.5	4.62	244	11.6	35.5	
H2	14.7	4.73	226	11.8	36.1	
Н3	14.7	4.75	242	11.8	36.3	
L1	0	0	0	0	0	
L2	0	0	0	0	0	
L3	0	0	0	0	0	
% age Carryover	0.00	0.00	0.00	0.00	0.00	