

MANJERI CLINICAL LABORATORY & IDEAL X-RAYS, MANJERI
PPM CHECKLIST FOR ELECTROLYTE ANALYZER

CBS 300 ELECTROLITE ANALYZER		SERIAL NO: BE010802	
I. SPECIAL PRECAUTIONS			
<p>If there is evidence of blood or serum contamination, submit the device for cleaning and decontamination before inspecting it.</p> <p>Wear appropriate Personnel Protection Equipment (PPE) during work.</p> <p>Refer to the Safety procedure for additional precautions and guidance.</p>			
II. TEST APPARATUS REQUIRED			
1.ELECTROLYTE ANALYZER		Model: CBS 300	Serial:BE010314
Description	status	Description	status
1. AC Plug	Passed	8. Rotary Valve	Passed
2. Line Cord	Passed	9. Fluid Flow	Passed
3. Fuse	Passed	10. Indicator	Passed
4. Tubes	Passed	11. Display	Passed
5. Connectors	Passed	12. User Calibration	Passed
6. Electrode Value	Passed	13. Pump Tube	Passed
7. Controls	Passed	14. Pinch Valve	Passed
		15. Air Sensor	Passed
III. PREVENTIVE MAINTENANCE TASKS			
Description	status	Description	status
1. Clean the whole tubings	Done	3. Replace Deteriorated components	Done
2. Clean the rotary valave	Done	4. Lubricate	Done