

Calibration Certificate, Semi Automatic BioChemistry Analyzer

Name	STAR - 20 (SEMIAUTO BIOCHEMISTRY ANALYZER)	
Hospital	YASH DIAGNOSTIC CENTER , VISHESHARPUR PACHAHATIYA	U.P. – (222001)
Date	20 July 2021	

This is certify that STAR-20 (SEMIAUTO BIOCHEMISTRY ANALYZER) have calibrated & produced following data..

Action Details:-

AD VALUE - (49200 To 51550) Limits - 35000 to 55000
 OFFSET VALUE- (34 To 126) Limits - 0.0 to 4000
 Water Blank Abs - 0.005 Limits - < 0.01

• Temp Value - < 5% Variation Limits - < 5%

Instrument Identification	on			QC Pas	sed
Serial no.		Date of Calibration	20.07.2021	Yes	Yes
Model	STAR-20	Due Date of Calibration	19.07.2022	Yes	Yes

After calibration Q.C. Test for all level passed successfully.

Test results is coming ok	✓	Remarks	Done
Muss = 0			
Service Engineer			
20 00			



Sales & Service

Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc. H-36, Muradganj, Naiganj, Jaunpur - 222002

Email.: pragysenterprises@gmail.com

Ref. :

Date: 01-10-208

Performance Qualification for Semy Automatic Bio Chemstry Analyser Star 20 ES1001905P9004

Calibration Parameters

Checked and found all the control [parameters to be with in the acceptable CV limits and in range.

Checked and found all controls to be within The acceptable SD.

System Certification:

Study data has determined that the system described in this document either meets thenecessary criteria outlined inthis Performance Qualification Protocol, or exceptional conditions have been identified and documentation included.

The system is ready for Specific Usage.

Protocol Performed By:

Name

: SAURAJ SINGH

Designation :Application

PERFORMANCE QUALIFICATION DONE BY

Name

: SAURAJ SINGH

Designation

: Application



Sales & Service

Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc. H-36, Muradganj, Naiganj, Jaunpur - 222002

Email.: pragysenterprises@gmail.com

Ref. :

Date: 01- 10 2081

External Requirement for Installation:

- Input voltage of 220V-240V/50Hz or 60 Hz.
- 2. Recommended operating Temperature is 15-30 degree Celsius, with in Relative Humidity 30-85% and Atmospheric pressure 70-106kPa.

Installation Qualification for Semy Automatic Biochemstry Analyser STAR 20 Sr. ES1001905P9004

Carried out all the installation procedures as per the installation procedure and checklists.

Carried out all the necessary checks and alignments.

Carried out all the necessary system checks and tests.

Handover the Instrument for operatos Training And Qualifications

For Pragya Enterprises

Technical Services Department



Name

: Sauraj Singh

Designation: Area Manager Sales & Service



Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc. H-36, Muradganj, Naiganj, Jaunpur - 222002

Email. : pragyaenterprises@gmail.com

Ref. :

Operational Qualification:

System Certfication:

Study data has determeined that the system described in this document either meets all criteria outlined in this Operational Protocol, or exceptional conditions have been identified and documentation included.

Exceptional Conditions, If any Have been Addressed.

The Syatem is ready for specific usage.

Protocol Performed by

: For Pragya Enterprises

Name: Sauraj Singh

Designation: Area Manager Sales & Service

Customer Authorization

:: YASH DIAGNOSTIC CENTER

VISHESHARPUR PACHAHATIYA JAUNPUR 222001

Name:

Designation:

Company Jaura / Singh

Date:

Customer Name And Sign

Date:



Sales & Service

Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc. H-36, Muradgani, Naigani, Jaunpur - 222002

Email.: pragyaenterprises@gmail.com

Ref. :

Date: 01-10-2021

Operation Qualification for Semy Automatic Bio Chemstry Analyzer Star 20

1. Verified all the Mechanical Movements : Done

2. Verifeied Hydraulics : Done

3. Verified Electrical Systems : Done

4. Verified the all reagents Systems : Done

For Pragya Enterprises

Technical Services Department

Name

: SAURAJ SINGH

Designation: Area Manager Sales & Service

Date

: 01-10-2021



Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc.

H-36, Muradganj, Naiganj, Jaunpur - 222002

Email.: pragyaenterprises@gmail.com

Ref. :

Date: 01-10-2081

Installation Qualification for Semy Automatic Bio Chemstry Analyzer

This is to certify that the instrument Star 20 Sr. No:- ES 100195P9004 is successfully installed the Installation Protocol/ checklist has been Successfully completed for the above Instrument.

Date of Installation : October 2018

For Pragya Enterprises

Technical Services Department

Name

: Sauraj Singh

Designation : Area Manager Sales & Service



Sales & Service

Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc. H-36, Muradganj, Naiganj, Jaunpur - 222002

Email.: pragyaenterprises@gmail.com

Ref. :

Date :

Operational Inspection:

- 1. Checked and found Mechanical Movements OK.
- 2. Checked and found Hydraulics OK.
- 3. Checked and found Electricals OK.
- 4. Checked with Controls And Samples, Results Are found OK.

For Pragya Enterprises

Technical Services Department



Name

: SAURAJ SINGH

Designation: Area Manager Sales & Service



Sales & Service

Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc.

H-36, Muradganj, Naiganj, Jaunpur - 222002

Email.: pragyaenterprises@gmail.com

Ref. :

Date: 01-10-2021

Installation Qualification for Semy Automatic BioChemstry Analyzer

Customer Name

: YASH DIAGNOSTIC CENTER

Address

: VISHESHARPUR PACHAHATIYA JAUNPUR - U.P.

Contact Person

: Mr. DHARMENDRA JI

Instrument Model

: STAR 20

Serial No.

: ES 1001905P9004

Date Of installation

: **OCTOBER 2018**

The instrument was installed And was found to be working satisfactory. Preliminary Customer Training was provided And standardization of the parameters Were done. The results were found to be within the expected range and system found to be working satisfactorily.

For Pragya Enterprises

Customer: Yash Diagnostic Center

. Technical Services Department

Name

: Sauraj Singh

Contact Person

: Mr. Anil Singh

Designation

: Area Manager Sales & Service

Designation

: Lab Incharge

Signature

Signature

Sand Sive

GSTIN: 09CKLPS0218M1ZO

Pragya Enterprises

Sales & Service

Medical Equipments Hospital Furniture & O.T. Equipments, Chemicals and Laboratory wares etc. H-36, Muradganj, Naiganj, Jaunpur - 222002

Email.: pragyaenterprises@gmail.com

Ref.:

Date : 01-10-200

Installation Qualification for Semy Automatic Biochemstry Analyser Star 20 No:-

ES1001905P9004

Customer Name

: YASH DIAGNOSTIC

Address

: VISHESHARPUR JAUNPUR 222001 – U.P.

Instrument Name

: STAR 20

Serial No

: ES1001905P9004

Initial Inspection of the unit carried out and the details are as follows:

System Condition Report:

Found the system to have been delivered in satisfactory condition and no external physical damaged observed on the same, Package was kept in good Condition as per the directional indicators like not tilt, indicating the system has not been subjected to mechanical shocks or stored in any manner, so as to cause any damage to the same.

Found all the required accessories as required.

Intallation Procedure And Checklist Attached for the records.

