



QC Data After PM.

Result Reprint

Report Type : Controls

Sr #	Lot #	Consumable	Test	Result Unit	Flag	Curve #	Result Date	Mean	SD	Interval (3SD)
1	B082182	Erba Norm	SGPTD	36.5 U/L		36391	24-May-2022 16:35:12	38.390	2.560	30.71 - 46.07
2	B082182	Erba Norm	SGOTD	51.3 U/L		36392	24-May-2022 16:35:30	51.504	3.230	41.814 - 61.194
3	B082182	Erba Norm	ALPU	131 U/L		36393	24-May-2022 16:35:49	125.920	8.520	100.38 - 151.48
4	B082182	Erba Norm	GGT	44.3 U/L		36394	24-May-2022 16:36:07	45.100	2.800	36.7 - 53.5
5	B082182	Erba Norm	BID	0.59 mg/dl		36395	24-May-2022 16:36:24	0.606	0.050	0.456 - 0.756
6	B082182	Erba Norm	BIT	1.47 mg/dl	-1SD	36396	24-May-2022 16:36:42	1.642	0.130	1.252 - 2.032
7	B082182	Erba Norm	UREA	37.4 mg/dl	+1SD	36397	24-May-2022 16:37:01	35.322	1.700	30.222 - 40.422
8	B082182	Erba Norm	CRENZ	1.09 mg/dl		36398	24-May-2022 16:37:18	1.140	0.070	0.93 - 1.35
9	B082182	Erba Norm	GLU	98.8 mg/dl		36399	24-May-2022 16:37:37	101.004	4.780	86.664 - 115.344
10	B082182	Erba Norm	TRIG	117.8 mg/dl		36400	24-May-2022 16:37:55	121.430	5.670	104.42 - 138.44
11	B082182	Erba Norm	CHOL	145 mg/dl		36401	24-May-2022 16:38:13	140.130	6.880	119.49 - 160.77
12	B082182	Erba Norm	HDLc	40.7 mg/dl		36402	24-May-2022 16:38:31	43.146	2.660	34.566 - 51.728
13	B082182	Erba Norm	UA	6.2 mg/dl		36403	24-May-2022 16:38:49	6.426	0.300	5.526 - 7.326
14	B082182	Erba Norm	CA	8.9 mg/dl	+1SD	36404	24-May-2022 16:39:07	8.212	0.400	7.012 - 9.412
15	B082182	Erba Norm	PRO	5.50 g/dl	-1SD	36405	24-May-2022 16:39:25	6.074	0.310	5.144 - 7.004
16	B082182	Erba Norm	ALBD	3.40 g/dl		36406	24-May-2022 16:39:43	3.556	0.180	3.016 - 4.096

QC Data After PM.

Result Reprint

Report Type : Controls

Sr #	Lot #	Consumable	Test	Result Unit	Flag	Curve #	Result Date	Mean	SD	Interval (3SD)
1	B082182	Erba Norm	SGPTD	36.5 U/L		36391	24-May-2022 16:35:12	38.390	2.560	30.71 - 46.07
2	B082182	Erba Norm	SGOTD	51.3 U/L		36392	24-May-2022 16:35:30	51.504	3.230	41.814 - 61.194
3	B082182	Erba Norm	ALPU	131 U/L		36393	24-May-2022 16:35:49	125.920	8.520	100.36 - 151.48
4	B082182	Erba Norm	GGT	44.3 U/L		36394	24-May-2022 16:36:07	45.100	2.800	36.7 - 53.5
5	B082182	Erba Norm	BID	0.59 mg/dl		36395	24-May-2022 16:36:24	0.606	0.050	0.456 - 0.756
6	B082182	Erba Norm	BIT	1.47 mg/dl	-1SD	36396	24-May-2022 16:36:42	1.642	0.130	1.252 - 2.032
7	B082182	Erba Norm	UREA	37.4 mg/dl	+1SD	36397	24-May-2022 16:37:01	35.322	1.700	30.222 - 40.422
8	B082182	Erba Norm	CRENZ	1.09 mg/dl		36398	24-May-2022 16:37:18	1.140	0.070	0.93 - 1.35
9	B082182	Erba Norm	GLU	98.8 mg/dl		36399	24-May-2022 16:37:37	101.004	4.780	86.664 - 115.344
10	B082182	Erba Norm	TRIG	117.8 mg/dl		36400	24-May-2022 16:37:55	121.430	5.870	104.42 - 138.44
11	B082182	Erba Norm	CHOL	145 mg/dl		36401	24-May-2022 16:38:13	140.130	6.880	119.49 - 160.77
12	B082182	Erba Norm	HDLc	40.7 mg/dl		36402	24-May-2022 16:38:31	43.148	2.860	34.568 - 51.728
13	B082182	Erba Norm	UA	6.2 mg/dl		36403	24-May-2022 16:38:49	6.426	0.300	5.526 - 7.326
14	B082182	Erba Norm	CA	8.9 mg/dl	+1SD	36404	24-May-2022 16:39:07	8.212	0.400	7.012 - 9.412
15	B082182	Erba Norm	PRO	5.50 g/dl	-1SD	36405	24-May-2022 16:39:25	6.074	0.310	5.144 - 7.004
16	B082182	Erba Norm	ALBD	3.40 g/dl		36406	24-May-2022 16:39:43	3.556	0.180	3.016 - 4.096

TRANSASIA

TECHNICAL SERVICE REPORT

No.

1411438

DATE: 24/6/22

CUSTOMER DETAILS		INSTRUMENT DETAILS		SERVICE STATUS	
NAME: Dr. Vohra's MediLab	MODEL: EM-200	<input type="checkbox"/> WARRANTY	<input type="checkbox"/> R&R		
ADDRESS: Panchkula	SR. NO.: 2002214	<input type="checkbox"/> AMC	<input type="checkbox"/> CMC		
		<input type="checkbox"/> CHARGED CALL			
		TYPE OF CALL			
		<input type="checkbox"/> INSTALLATION			
		<input type="checkbox"/> P.M. VISIT <input type="checkbox"/> I <input type="checkbox"/> II <input type="checkbox"/> III <input type="checkbox"/> IV			
		<input type="checkbox"/> APPLICATION SUPPORT			
		<input type="checkbox"/> BREAKDOWN			
TEL NO:	JOB COMPLETED	DOWN TIME:			
NAME OF THE OPERATOR:	RESPONSE TIME	COUNTER READING:			
PROBLEM REPORTED:	TRAVEL TIME				
OBSERVATIONS:					
ACTION TAKEN: Done preventive Maintenance of the Instrument, wash cassette clean probe, try tubing, biopsy section, replace lamp, Run Quality control, sample, checked Instrument functioning satisfactory					
SITE CONDITION : LINE-NEUTRAL VOLT. : _____ NEUTRAL-EARTH VOLT. : _____ LINE-EARTH VOLT. : _____					
BRAND OF REAGENT USED :		TO BE FILLED IN BY CUSTOMER			
<input type="checkbox"/> FOLLOWING PARTS HAVE BEEN REPLACED		<input type="checkbox"/> PREVENTIVE MAINTENANCE CARRIED OUT SATISFACTORILY.			
<input type="checkbox"/> FOLLOWING PARTS NEED TO BE REPLACED. PLEASE APPROVE		<input type="checkbox"/> FAULT RECTIFIED & INSTRUMENT IS WORKING SATISFACTORILY.			
NO.	DESCRIPTION	QTY.	COST	TOTAL	<input type="checkbox"/> WE HEREBY APPROVE RS. _____ FOR PARTS
					<input type="checkbox"/> COMMENTS (IF ANY):
TOTAL Rs.					
INVOICE NO. :		DATE :		BRANCH H. O.	
FOLLOW-UP ACTION (Required if any) :		RECEIVED ON :		CHECKED BY :	
ENGINEER'S/PRODUCT SPECIALIST'S SIGNATURE :		JOB CARD NO. :			
TIME :	NAME: Nikhil Singh				

NOTE: Parts replaced are chargeable except during warranty Period. Consumables like printer head, lamp, tubing, paper rolls etc. & breakable parts are not covered by warranty and hence are chargeable. Parts replaced due to negligence in operation will also be charged in every case.

AT TRANSASIA, CUSTOMER SATISFACTION IS OUR PRIME CONCERN. IN CASE YOU HAVE ANY SUGGESTIONS PLEASE CONTACT : GENERAL MANAGER (TECHNICAL SERVICE), MUMBAI. TEL. : 4030 9000 FAX : (022) 4030 9090

TRANSASIA BIO-MEDICALS LTD.

Head Office : Transasia House, 8, Chandivali Studio Road, Andheri (E), Mumbai - 400 072. Tel. : (022) 4030 9000 Fax : (022) 2857 3030
Mumbai : Toll Free 1800 103 8226

Indore : Tel. : (0731) 403 3948

Lucknow : Tel. : (0522) 443 0900

Chandigarh : Tel. : (0172) 464 3837

Jaipur : Tel. : 93092 63147

Patna : Tel. : 77669 00553

Kochi : Tel. : (0484) 402 0511

Chennai : Tel. : (044) 2822 7149

Pune : Tel. : (020) 2616 2658

Hyderabad : Tel. : (040) 4430 9000

Guwahati : Tel. : 91369 50480

Chuttack : Tel. : (0671) 232 4555

Hubli : Tel. : (0836) 485 0900

Kolkata : Tel. : (033) 2215 7839 Fax : 2215 8035

Dehl : Tel. : (011) 2578 5451/2573 2223

Ahmedabad : Tel./Fax : (079) 02640 7030/4032 1903

Bangalore : Tel./Fax : (080) 2556 8044

DOC NO : SC00-406 / ISS - 5

www.transasia.co.in

ISO 13485 & ISO 9001 Certified Co.

CUSTOMER COPY

Erba
Mannheim