

Transasia Bio-Medicals Ltd., Transasia House, 8 Chandivali Studio Road, Andheri (East), Mumbai - 400 072 Tel: +91 22 4030 9000 Fax: +91 22 2857 3030 Email: transasia@transasia.co.in CIN: U33110MH1985PLC036198



Custon	ner Name:			CALIE Lahoti Medi	RATION		ATE	Data	12 84 2000		
		mated Bioc	-			iu.		Date	13 May 2022	4	
		matea bioc	neillisti y A	naiyzer One	<u>:111 7</u>			-			
SL. NO.	N110714										
•	1				Functiona	I Check:					
1	Pump Calibr	Pump Calibration						2062			
	(Calibrate th	e Peristaltic	oump.								
2	Key Functions						Pass				
	(Confirm all	keys function	& keys not h	ard to press)							
3	Printer Test Pass										
	(Confirm that all characters printed w/o gap & are clearly redeable										
4	Buzzer						Pass				
	(Confirm buzzer beeps when press any key)										
5	VOLTAGE C	CHECK	Required				Actual Voltage				
	LAMP VOLT TEST +12V		12V (+/- 0.5V)				11.68V				
	LINE VOLT		230V (+/- 10)V)			228v				
	SMPS		17V/ DC to 19V /DC 5.01V (+/- 0.5V)				17.7V 5.1V				
	TEST +5V										
	TEST +12 V 12.00V (+/- 0.5V)						12.1V				
6	6 PRE & LOG CHANNEL READINGS FOR THE FILTER										
Filter		340	405	450	505	546	578	630	670	Range	
Gains		169k	169k	170k	170k	170k	172k	171k	169k	30k ~ 200k	
Offset		Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	NA	

All above test has been performed and it passed .

This is to certify that the above mentioned Automated Biochemistry Analyzer IS Technically Calibrated and it is working OK.

Technical Sevice Deptt. -Transasia Bio-Medicals Ltd.

