19-02-12



#### Installation Qualification

Instrument Type: Blood Cell Counter	Manufacture Name: ERMA Japan INC
	Supplier Name: AGD Biomedicals (P)
	LTD

### **Operational Verification Report**

Instrument Type: Blood Cell Counter	Name of the Customer: Care Hospital, Bareilly
Serial Number: 29210	Model Number: ERMA PCE 210
Date of Installation: 11.03.2011	Date of receipt: 10.03.2011

Power Check:	
Power Supply:	220 V AC
Earth Voltage	OK

Physical Check of the equipment

Checked	No damage found
Checked	OK
Checked	No damage found
	Checked

System configuration equipment

Positioning of the equipment – flat stable surface	Checked	Found OK
Adequate space around the equipment	Checked	Found OK
Interference to other equipment	Checked	No Interference



Result:

Installed the ERMA PCE 210 satisfactorily and is ready for operational verification

Installed by: Mr Harivansh Yadav	Checked by:	box Marika hod.
Signature: Dalimet	Signature:	
Date: 11.03.2011	Date:	11/02/2114



### Performance Qualification

Instrument Type: Blood Cell	Manufacture Name: ERMA Japan
Counter	INC
Model Name: ERMA PCE 210	Supplier Name: AGD Biomedicals (P)
	LTD
Installed On: 11.03.2011	Sales & Service Engineer: Mr
	Harivansh Yadav

# Performance Verification Report

Name of the Hospital: Care Hospital,	Instrument Type: Blood Cell Counter
Bareilly Model Number: ERMA PCE 210	Serial Number: 29210
Date of receipt: 10.03.2011	Date of Installation: 11.03.2011

## Performance Check:

Background Check	OK
Pressure Check	OK
HGB Unit Check	OK
Syringe Movement Check	OK
Printer Check	OK



### Result:

The performance verification of ERMA PCE 210 has been completed satisfactorily and Instrument is ready to be placed in service.

HOD:
Signature:
Date:
Checked by:
Date:

Installation, Operational and Performance Qualification Certificate
The equipment ERMA PCE 210 Serial No 29210 is installed and
commissioned. The operational and performance verification is done and
observations are:

- All system components are working in normal condition
- All tests were carried out in accordance with qualification protocol and standard operating procedure.
- The performance of the ERMA PCE 210 is as per the acceptance criteria.