

CERTIFICATE OF CALIBRATION

ADL/KO/SER/0015/22-23

Instrument: Mispa CXL PRO Fully Automated Biochemistry Analyser

Customer Name: Sevana Hi-Tech Lab, Parappur, Thrissur

Date of Installation : 09/11/2021

S/N: 210T180CS0345K

Alignment Check

Motion Check	Checked	Remark
Sample Probe	<input checked="" type="checkbox"/>	OK
Reagent Probe	<input checked="" type="checkbox"/>	OK
Mixing Bar	<input checked="" type="checkbox"/>	OK
Reaction Disk	<input checked="" type="checkbox"/>	OK
Reagent Disk	<input checked="" type="checkbox"/>	OK
Sample Disk	<input checked="" type="checkbox"/>	OK
Cuvette Check	<input checked="" type="checkbox"/>	OK
Lamp	<input checked="" type="checkbox"/>	OK

System	Checked	Remark
Main unit	<input checked="" type="checkbox"/>	OK
Temp. unit	<input checked="" type="checkbox"/>	OK
Mixing unit	<input checked="" type="checkbox"/>	OK
Reaction unit	<input checked="" type="checkbox"/>	OK
Reagent unit	<input checked="" type="checkbox"/>	OK
Sample unit	<input checked="" type="checkbox"/>	OK

AGAPPE DIAGNOSTICS LTD. ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC1154

CORPORATE OFFICE / REAGENT PLANT
Agappe Hills, Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India.
Tel: +91 484 286 7000 | Email: agappe@agappe.in

EQUIPMENT PLANT
K/588-CB, Block No. 12, KINFRA Small Industrial Park,
Nellad, Cochin, Kerala, India - 686 723. Tel: +91 454 276 7477.

KOLKATA OFFICE
406, Merlin Matrix, Plot No-10, Block-DN, Sector
Salt Lake City, Kolkata - 700 091.
Tel: +91 33 4003 0451 | Email: kolkataoffice@agappe.in

MUMBAI (REGISTERED OFFICE)
403 & 402, 4th Floor, Jaisingh Business Centre, 119,
Sahar Road, Parsiwada, Andheri (East), Mumbai - 400 099, India.
Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

DELHI OFFICE
D5M 540, 5th Floor, DLF Tower, Shivaji Marg,
New Delhi - 110 015, India.
Tel: +91 11 4558 8416 | Email: delhioffice@agappe.in

BANGALORE OFFICE
5-6, 11th Floor, Red Cross Bhavan, No.26 Race Co
Bangalore - 560 001. Tel: +91 80 2228 8288
Email: bangaloreoffice@agappe.in

OPTICAL SECTION

Filter	Light Intensity
340	9500
380	9750
405	9850
450	10,100
480	10,250
505	11,000
546	11,500
570	11,750
600	12,010
660	12,500
700	13,000
800	13,200

Power Supply :

12V 24V

The equipment is calibrated for the different parameters by using Multi Calibrator

Lot.No: 32060079 Exp: 06/2024

This is to certify that the above mentioned product has been calibrated according to the standard operating procedure given by the Agappe Diagnostics.

Date of Calibration : **09/05/2022**Next Due : **08/11/2022**

For, Agappe Diagnostics Ltd,

Authorised signatory

Syam Mohan

Engineer - Customer Service

AGAPPE DIAGNOSTICS LTD.

ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115

CORPORATE OFFICE / REAGENT PLANT
 Agappe Mills, Pattinattom IPOI, Dist. Ernakulam, Kerala - 683 562, India.
 Tel: +91 484 286 7000 | Email: agappe@agappe.in

MUMBAI (REGISTERED OFFICE)
 401 & 402, 4th Floor, Jaisingh Business Centre, 119,
 Sahar Road, Parsiwada, Andheri (East), Mumbai - 400 099, India.
 Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

EQUIPMENT PLANT
 K/388-CB, Block No. 32, KINFRA Small Industrial Park,
 Nelladi, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

DELHI OFFICE
 DSN 540, 5th Floor, DLF Tower, Shivaji Marg,
 New Delhi - 110 015, India.
 Tel: +91 11 4558 8416 | Email: delhioffice@agappe.in

KOLKATA OFFICE
 406, Merlin Matrix, Plot No-30, Block-DN, Sector
 Salt Lake City, Kolkata - 700 093.
 Tel: +91 33 4003 0451 | Email: kolkataoffice@agappe.in

BANGALORE OFFICE
 S-6, 11th Floor, Red Cross Bhavan, No.26 Race Coe,
 Bangalore - 560 001. Tel: +91 80 2228 8288
 Email: bangaloreoffice@agappe.in