



Case No.	CAS-0014673900	Instrument Model	ELECTROLYTE ANALYZER W/O STARTERKIT 9180
Order No.	ORD-0018944526	Instrument Serial No.	12730
Contract Type	PREVENTATIVE MAINTENANCE & PARTS	Finance Status	CASH
Lab/Inst./hosp.Name	SRL Diagnostic-Dr. Avinash Phadke		
Customer No.	0052086132		
Contact Name :	Thanuja		
Contact Number :	+912225623742		
Address :	Sumati, Ground floor,		
City :	Mulund West		

Call Received Date/Time:	26.12.2021 17:30	Call Attended Date/Time:	09.03.2022 10:00
--------------------------	------------------	--------------------------	------------------

Job Type	PM Visit
----------	----------

Job Description	
-----------------	--

Action Summary	Cause: Preventive Maintenance Workdone:1) checked the system & checked all modules, 2) Carried out preventive maintenance satisfactorily as per checklist along with p.m. kit, 3) Performed the calibration & done samples Spare consumable used: Gmmi no. 03074064001 Item: Harness Main Tubing set Qty: 01 Batch: 354121 Note: Replaced under comfort AMC Verification:Machine working ok Customer Satisfaction Rating (1-5):5
----------------	--

Spare Part Replaced

Part No	Parts Description	Batch No	Batch Expiry Date	Qty	Invoice Type
03074064001	HARNESS, MAIN TUBING, 918X			1	Contract

Time Report

Effective Visit Date : 09.03.2022	Complete Date : 09.03.2022
-----------------------------------	----------------------------

Date	Type	Time
2022.03.09	Travel Time--Standard	1
2022.03.09	Working Time--Standard	1
Total:		2


Customer's Signature Name : Thanuja	Service Engineer/Application Specialist Name : Machindar Parappa Bilkar
--	--

Thanuja

Machindar Parappa Bilkar

Date: 09.03.2022	Date: 09.03.2022
------------------	------------------

29/7/22





control ①

	Na ⁺	K ⁺	Cl ⁻
low ⁻	135	3.7	94.2
High ⁻	149	4.20	110
Target ⁻	142	3.95	102

[Handwritten signature]

Done By - Vaishali K.
Checked by - Priyanka



Analyzer


Shot on OnePlus
By vaishali

SRL DIAGNOSTIC
MUMBAI

Control (2)

	Na^+	K^+	Cl^-
low -	120	5.72	78.9
High -	133	6.44	92.3
Target -	127	6.08	85.6

M




Control-①

Low - Na⁺ K⁺ Cl⁻
 High - 135 3.71 94.2
 Target - 142 3.95 102

Na	K	Cl
145	4.1	98

80 Electrolyte Analyzer



- Done by Vaishali Kalambe,
 Checked by Priyanka

Dial RPS On-line
 Customer Support Centre
 MON to SAT 8am to 8pm

SRL DIAGNOS

31/7/22


② control ②

Na K Cl

Target - 120 5.72 78.9

High - 133 6.44 92.3


Low - 127 5.08 85.6



Electrolyte Analyzer

✓
[Signature]

- Done by Vaishali Kalambe.




SRL DIAL
MUL
HE
[Signature]
MALHODI

SERIAL NO. NAME OF TECHNICIAN
 01 MS. TEJASWI KADAM
 02 MS. SHRADDHA VALI NDE
 03 MS. KHYATI MUNDHE

Na K Cl
 143 3.9 97

Electrolyte Analyzer



SRL DIA
 MUL
 HE

Control ①

	Na ⁺	K ⁺	Cl ⁻
low	135	3.71	94.2
high	149	4.20	110
target	142	3.95	102

Done By - Vaishali K.
 Checked by - Priyanka

SEMI-AUTOMATIC SERUM ELECTROLYTE ANALYZER

SR NO. NAME OF TECHNICIAN
MS. JESSICA KADAM

↓ Na ↑ K ↓ Cl
128 6.2 85

Electrolyte Analyzer

SRL

Dialysis Online
Customer Support Centre

Control ②

	Na ⁺	K ⁺	Cl ⁻
Low -	120	5.72	78.9
High -	133	6.44	92.3
target -	127	6.08	85.6

Checked by Poojanka

CONTROLLED COPY
LIN: 29
Issue Date: 20.06.2019
Auth.: Supt 5
QA Function

DAILY QC LOG FOR ELECTROLYTE ANALYZER

MONTH / YEAR :- July-2022
NAME OF THE INSTRUMENT :- 9180 Electrolyte

EXPIRY DATE :- 31/07/2023
INSTRUMENT SERIAL NO. :- 12730

	LEVEL - 1			LEVEL - 2			LEVEL - 3			INITIAL	SUPERVISOR'S SIGN
	Na	K	Cl	Na	K	Cl	Na	K	Cl		
Units	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L	mmol/L		
Low	135	3.7	94.2	120	5.7	78				<u>Supt</u>	<u>[Signature]</u>
High	149	4.2	110	133	6.4	92					
I.V	142	3.9	102	127	6.0	85					
1											
2											
3											
4											
5											
6											
7											
8											
9											
10											
11											
12											
13											
14											
15											
16											
17											
18											
19											
20											
21	147	4.1	102	129	6.2	86				<u>Supt</u>	<u>[Signature]</u>
22	147	4.1	100	128	6.2	85					
23	141	3.8	95	123	5.9	82					
24	140	3.9	95	125	5.9	82					
25	141	3.9	95	124	5.9	82					
26	141	3.9	95	126	6.0	83					
27	142	4.0	96	124	6.0	82					
28	142	4.0	97	125	5.9	84					
29	142	3.8	97	125	5.8	83					
30	143	3.9	97	128	6.2	85					
31	145	4.1	98	122	5.9	81					

SECTION HEAD SIGN: [Signature]

Doc No.: SRLD/BIO/F/10	Daily QC Log for Electrolyte		Page 1 of 1
Issue No.: 01	Issue Date: 20.06.2019	Amend No.: 00	<u>Supt 5</u>
Issued By: Mr. Surendra Gupta (DGM-Quality Assurance)			